

Instructions

Lateral Trunk Support Double Flap

Important Information

The Double Flap Lateral Trunk Support brings the trunk to midline position.

- The Inner flaps mobilize the trunk.
- The Outer flaps centralize the trunk and keep the arms in the front of the chair.
- The Scolio strap aligns the spine.




STEP 1

Attach the Lateral Trunk Support to the chair by wrapping the two rear straps around the back of the chair and securing with Velcro.



STEP 2

To attach the double flap, wrap the strap on the left-hand outer flap around the side tubing of the frame and secure with Velcro.



STEP 3

Thread the strap on the left-hand inner flap into the plastic loop on the outer flap working from front to back.



STEP 4

Adjust and secure strap. (Repeat steps 3 and 4 on right-hand flap.)




STEP 5

Bring the long scoliosis strap through the plastic loop on the inner flap on the right-hand side and across the chest.



STEP 6

Thread the strap through the plastic loop on the inner flap on the left-hand side.



STEP 7

Fasten the strap with Velcro.



STEP 8

If additional tightness is desired, pull the scoliosis strap further and Velcro it beyond the inner flap.

Note:

The Lateral Trunk Support – Double Flap can be used to support either the left or right side of the trunk.

- **To support the right side** of the trunk, attach the Support to the chair with the double-sided Velcro strap on the right.
- **To support the left side** of the trunk, attach the Support to the chair with the double-sided Velcro strap on the left.