

### Folding Instructions



**Caution: Keep fingers free of folding mechanism.**



**TIP:** The Safari Tilt **MUST** be fully upright and in a locked position before folding.



**SPECIAL INSTRUCTIONS FOR 3" LEGREST EXTENSION:**

Legs must be slightly elevated to allow for proper folding. When unfolding, raise the ELR's\* to a 180° angle (parallel to the push handles).

\*ELR = Elevating Legrest



**STEP 1**

Place foot on rear lock brace and bring the chair to the most upright position (until you hear a "click").



**STEP 2**

Flip up armrests.



**STEP 3**

Remove seat cushion.



**STEP 4**

Swing footplates to the side.

**Note:** One-piece footplate must be removed from chair before folding.



**STEP 5**

Unlock rear lock brace by pushing up on center of brace with foot.




**STEP 6**

Lean chair back on push handles and lower chair to ground. Walk to front of chair and push down on tubing above front tires until chair is folded.



**STEP 7**


For a more compact fold, pull up on slide lock and fold the ELR's toward the front of chair.



**STEP 8**

Fasten the closure strap. **Do not use the closure strap to carry.**

## Unfolding Instructions



**STEP 1**

Unfasten closure strap.




**STEP 2**

Pull up on the slide locks and raise the ELR's until they are at a 90 – 180 degree angle to the folded chair.



**STEP 3**

Stand the folded chair vertically on rear wheels. Grasp the push handle closest to you and at the same time push down on the ELR closest to you using a forward motion.



**STEP 4**

Push down with your foot on the rear lock brace until it is in the locked position.



**STEP 5**

Attach the seat cushion by pushing firmly into the back of seat, then pressing down onto Velcro ®.



**STEP 6**

Flip down armrests.



**STEP 7**

Swing footplates forward.