

Wheelchair Features / Model	EZ Rider	EZ Rider Planar	Convertible	Convertible Planar	Life Rider	Vivo	Cruiser	Cruiser Classic	Cruiser Planar	Cruiser Scout	Cruiser Profiler	CruisAire	Metro	CuddleBug	Rodeo	Safari	Mountee
Made in the USA with Lifetime USA Warranty	✓	✓	✓	✓	5-year	5-year	✓	✓	✓	✓	✓	5-year	✓	✓	✓	✓	✓
Stable on uneven terrain for increased safety	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Replacement parts available - USA warehouse	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
10-degree fixed tilt (upright) for safety and comfort	✓	✓	✓	✓	✓	✓											
30-degree fixed tilt for safety and comfort							✓	✓	✓	✓	✓	✓	✓				
Adjustable tilt wheelchair for comfort & repositioning														✓	✓	✓	✓
Adjustable back angle for increased comfort						✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	
Adjustable recline for pressure relief						✓								Option		Option	
One-piece height adjustable push handle							✓		Option					✓	✓		
Attendant hand brakes for increased control & safety							Option			✓						Option	
High / Low height adjustment for interaction with peers														✓			
Removable large rear wheels for self propelling			✓	✓													
Adjustable strap back for increased comfort						✓	✓				✓				✓		
Comprehensive strap back system-custom positioning											✓						
Solid seat and back for added support		✓		✓					✓					✓			✓
Seat cushion for comfort, stability and support	Option		Option				Option	Option		Option	✓					✓	✓
Elevating legrests for comfort and positioning														✓	✓	✓	✓
2-piece height adjustable, swing-away footplates	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			Option	Option
One-piece height adjustable footplate														✓	✓	✓	✓
Toddler size wheelchair - 6 months to 4 years														✓			
Pediatric size wheelchair	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Adult size wheelchair	✓	✓	✓	✓			✓	✓	✓	✓			✓			✓	
Compact fold to fit in small spaces	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Lightweight, easy to fold up, store & travel	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Comfort and positioning for all-day seating	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
HCPCS coded as a wheelchair for funding	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓	✓	✓
Adjustable seat depth for growth	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓
Multiple in-depth positioning options	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓
Large variety of fabric colors	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓
Meets ANSI/RESNA safety standards	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
WC19 crash tested transit wheelchair	Option	Option	Option	Option			Option	Option	Option	Option	Option		Option	Option	Option	Option	Option
Wheelchair with all-terrain tires							Option		Option	✓						✓	✓
Under seat storage basket for convenience							✓	✓	✓	✓			Option	✓	✓		
Washable fabric for EZ cleaning	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Quick release wheels with hub brakes	✓	✓					✓		✓				✓	✓		✓	