

# Convoid Vivo

## Retail Price / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_

Contact: \_\_\_\_\_

Bill To:		Ship To:	
Address:		Address:	
City/State/Zip:			
Phone:	Fax:	City/State/Zip:	
Email:		Phone:	
<b>Vivo Standard Features:</b>			
10-Degree Fixed Tilt	2-Point Positioning Belt	Closure Strap	
95 - 120 Degrees of Adjustable Recline	6" Front, 8" Rear Solid Tire	Adjustable Back Upholstery	
Removable Adj. Swing-Away Footplates	Toggle Wheel Locks	Made in the USA	
<b>Standard Models:</b>		<b>Price</b>	
<input type="checkbox"/> VIVO VV12 - 903426		<b>\$900</b>	
<input type="checkbox"/> VIVO VV14 - 903556		<b>\$920</b>	

**Upholstery Colors - Pac Cloth**

- |  |            |
|--|------------|
| <input type="checkbox"/> Electric Blue | <b>Std</b> |
| <input type="checkbox"/> Plum          | <b>Std</b> |
| <input type="checkbox"/> Turquoise     | <b>Std</b> |
| <input type="checkbox"/> Red           | <b>Std</b> |
| <input type="checkbox"/> Purple        | <b>Std</b> |
| <input type="checkbox"/> Orange        | <b>Std</b> |
| <input type="checkbox"/> Black         | <b>Std</b> |

VV12	VV14
903487	903487
903488	903488
903489	903489
903490	903490
903491	903491
903492	903492
904175	904175

**Frame Option**

- |  |              |              |
|--|--------------|--------------|
| <input type="checkbox"/> Height Adjustable, Swing-Away Armrests (pair) | <b>Price</b> | <b>\$203</b> |
|--|--------------|--------------|

VV12	VV14
903438	903557

**Pelvis Positioning**

- |   |              |
|---|--------------|
| <input type="checkbox"/> 3-Pt Positioning Belt with Depth Adjustable Crotch Strap | <b>\$49</b>  |
| <input type="checkbox"/> Short Seat Depth - 8.25" to 10.25"                       | <b>N/C</b>   |
| <input type="checkbox"/> Foam Wedge   | <b>\$152</b> |

VV12	VV14
904017	904017
904121	904122
903427	903559

**Foot & Leg Positioning**

- |  |              |              |
|--|--------------|--------------|
| <input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)  | <b>Price</b> | <b>\$156</b> |
| <input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair) | <b>\$156</b> |              |
| <input type="checkbox"/> Footplate Securement Strap                | <b>\$30</b>  |              |
| <input type="checkbox"/> Medial Thigh Support (Abductor)           | <b>\$162</b> |              |
| <input type="checkbox"/> Lateral Thigh Support (Adductor)          | <b>\$162</b> |              |

VV12	VV14
903415	903415
903416	903416
903584	903584
903429	903561
903430	903562

**Trunk Positioning**

- |  |              |              |
|--|--------------|--------------|
| <input type="checkbox"/> H-Harness with Padded Covers (3-Pt Positioning Belt Req.) | <b>Price</b> | <b>\$143</b> |
| <input type="checkbox"/> Anterior Trunk Support                                    | <b>\$191</b> |              |

VV12	VV14
903988	903988
903455	903564

**Head Supports**

- |   |              |              |
|---|--------------|--------------|
| <input type="checkbox"/> Occi Headwings     | <b>Price</b> | <b>\$228</b> |
| <input type="checkbox"/> Headrest Extension | <b>\$130</b> |              |

VV12	VV14
903433	903567
903460	903566

**Accessories**

- |   |              |              |
|---|--------------|--------------|
| <input type="checkbox"/> Headrest Cover (Canopy)                              | <b>Price</b> | <b>\$228</b> |
| <input type="checkbox"/> Grab Bar ^   | <b>\$151</b> |              |
| <input type="checkbox"/> Travel Bag (For Wheelchair Only)                     | <b>\$206</b> |              |
| <input type="checkbox"/> Upper Extremity Support Surface (Padded Tray) ^      | <b>\$101</b> |              |
| <input type="checkbox"/> Adjustable Removable Hardware (Req. for Padded Tray) | <b>N/C</b>   |              |

VV12	VV14
903436	903569
903437	903570
902022	902022
903435	903568
904119	904120

^ Armrests required.

PO Box 4209, Palos Verdes, CA 90274 ~ 1-888-CONVAID (266-8243), 310-618-0111 ~ Fax: (310) 618-2166