

## Convertible Features

HCPCS: E1236












Convertible

- 10° Fixed Tilt  
*Promotes an upright position for individuals with good trunk strength and upper body stability*
- Quick Release Large Rear Wheels  
*Lets you choose between independent or dependent mobility*
- Use Large Wheels (Independent Mobility)  
*To self propel for exercise and build upper body strength*
- Remove Large Wheels (Dependent Mobility)  
*To maneuver chair in and out of tight spaces (restaurant, theater, etc.)*
- Lightweight & Compact Folding Frame  
*Is easy to lift, carry and transport – fold and go!*
- Comprehensive Positioning Options  
*Address all types of positioning needs*
- Custom Modifications  
*Achieve exact fit for any size or shape*
- Multiple Frame Adjustments  
*Extend the life of the chair for a growing child*
- Crash Tested Transit Models  
*Adhere to national safety standards, and provide safe transport in bus or van*



compact folded  
**CONVERTIBLE**

# Convertible Specifications

Measurement	Chair Model	 Seat Width*	 Seat Depth**	 Seat Back Height	 Seat to Footplate*	 Seat to Floor*	 Seat-to-Back Angle	 Seat Angle	 Headrest Extension	 Shoulder Strap Height	 Weight Capacity Regular/Transit	 Weight of Chair <sup>Δ</sup>
in.	<b>CV 12</b>	12	10 - 13	17.5	2.5 - 19	17	100°	10°	5	11 - 16	75 lbs / 66 lbs	31 lbs
cm	<b>CV 12</b>	30.5	25 - 33	44	6 - 48	43			13	28 - 40.5	34 kg / 30 kg	14 kg
in.	<b>CV 14</b>	14	12 - 15	20	2.5 - 21	19.5	100°	10°	8	10 - 16	100 lbs / 100 lbs	36 lbs
cm	<b>CV 14</b>	35.5	30.5 - 38	51	6 - 53	49.5			20	25.5 - 40.5	45.5 kg / 45.5 kg	16 kg
in.	<b>CV 16</b>	16	13 - 16	18.75	2.5 - 21	21	100°	10°	8	10 - 16	170 lbs / 170 lbs	39 lbs
cm	<b>CV 16</b>	40.5	33 - 40.6	47.5	6 - 53	53			20	25.5 - 40.5	77 kg / 77 kg	18 kg
in.	<b>CV 18</b>	18	14 - 17	19.5	2 - 22	22	100°	10°	10.5	13 - 17	250 lbs / 170 lbs	36.5 lbs
cm	<b>CV 18</b>	45.5	35.6 - 43.2	49.5	5 - 56	56			27	33 - 43	114 kg / 77 kg	16.5 kg


in. - Inches, cm - Centimeters \* Plus or minus .5 inch (1.27 cm) due to fabric. \*\* 2 inch cut-out for minimum seat depth. Available as an option. Δ - Weight of chair is without front rigging.

## Upholstery Colors


Removable, Washable Fabric




Ocean Blue




Purple










Navy Blue



Teal



Raspberry

Measurement	Chair Model	 Overall Width (with large rear wheels)	 Overall Height (with large rear wheels)	 Overall Length (with large rear wheels)	 Height of Rear Wheels	 Folded Length	 Folded Width	 Folded Height
in.	<b>CV 12</b>	25.5	36.5	43	20	43	12.5	14.5
cm	<b>CV 12</b>	65	93	109	51	109	32	37
in.	<b>CV 14</b>	26.5	36	45	22	44	12	14.25
cm	<b>CV 14</b>	67	92	114.3	56	112	30.5	36.2
in.	<b>CV 16</b>	28.5	37	47.5	24	45.5	15	16
cm	<b>CV 16</b>	72	94	120.7	61	115.5	38	40.6
in.	<b>CV 18</b>	30	37	48.5	24	45	14.5	15.25
cm	<b>CV 18</b>	76	94	122.2	61	114	36.8	38.7

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.

Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

## Limited Lifetime Warranty

Frame and X-Braces.....Lifetime of original retail buyer  
Fabric and webbing.....One year  
Other components.....One year

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.

