

## CruisAire Features

HCPCS:E1236

CRUISAIRE

- **30° Fixed Tilt**  
*Promotes upper body positioning for individuals with limited upper body strength or trunk stability*
- **Ultra-Lightweight Frame**  
*Is easy to lift and carry*
- **Compact Fold**  
*Offers unbeatable portability – easy to stow and transport*
- **Comprehensive Positioning Options**  
*Provide comfort and optimal positioning*
- **Multiple Frame Adjustments**  
*Extend the life of the chair for a growing child*
- **Rigorous Product Testing**  
*Ensures maximum safety, strength, stability and endurance*



# CruisAire Specifications

Measurement	Chair Model	Seat Width*	Seat Depth**	Seat Back Height	Seat to Footplate*	Seat to Floor*	Seat-to-Back Angle	Seat Angle	Headrest Extension	Shoulder Strap Height	Weight Capacity	Weight of Chair <sup>Δ</sup>
in.	<b>CA 11</b>	11	8 - 11	19	4 - 15	22	85° - 90° - 95° - 100°	30°	N/A	9 - 15	66 lbs	14 lbs
cm	<b>CA 11</b>	28	20.3 - 28	48	10 - 38	56				23 - 38	30 kg	6.3 kg
in.	<b>CA 13</b>	13	10 - 13	21	6 - 16.5	24	85° - 90° - 95° - 100°	30°	5.5	12.5 - 18.5	75 lbs	15 lbs
cm	<b>CA 13</b>	33	25 - 33	53	15 - 42	61			14	32 - 47	34 kg	6.8 kg

in. - Inches, cm - Centimeters \* Plus or minus .5 inch (1.27 cm) due to fabric. \*\* 2 inch cut-out for minimum seat depth. Available as an option. <sup>Δ</sup> - Weight of chair is without front rigging.

Measurement	Chair Model	Overall Width	Overall Height	Overall Length***	Folded Length	Folded Width	Folded Height
in.	<b>CA 11</b>	20	37	33	46	12.5	11
cm	<b>CA 11</b>	51	94	84	117	32	28
in.	<b>CA 13</b>	21.5	34	34	46.5	14	13
cm	<b>CA 13</b>	55	94	86.4	118	36	33

\*\*\* without footrest

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.


Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

## Limited Lifetime Warranty

Frame and X-Braces.....Five years  
 Fabric and webbing.....One year  
 Other components.....One year

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.

## Upholstery Colors



## CruisAire Anatomic

The CruisAire Anatomic Back enhances posture and improves line of vision. Recommended by therapists, the Anatomic Back provides lumbar support, which stabilizes the neck and head to assist the development of righting reflexes. It comes with ABS stiffeners, which can be custom contoured.

