

# Convaid Cruiser (Classic) - HCPCS E1236

## Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_

Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone: Fax:	City/State/Zip:
Email:	Phone:

### Cruiser Classic Standard Features:

Two-Piece Push Handle	6"Front, 8"Rear Solid Tires (CX10/12)	30-Degree Fixed Tilt
Silver Powder Coated Frame	8"Front, 8"Rear Solid Tires (CX14/16)	2-Point Positioning Belt
Removable Adj. Swing-Away Footplates	8"Front, 11"Rear Solid Tires (CX18)	Toggle Wheel Locks

#### Standard Models:

<input type="checkbox"/> Cruiser CX10 - 903314	<b>Price</b>	<b>\$1,876</b>
<input type="checkbox"/> Cruiser CX12 - 902845	<b>Price</b>	<b>\$1,876</b>
<input type="checkbox"/> Cruiser CX14 - 900490	<b>Price</b>	<b>\$1,876</b>
<input type="checkbox"/> Cruiser CX16 - 900145	<b>Price</b>	<b>\$2,036</b>
<input type="checkbox"/> Cruiser CX18 - 902594	<b>Price</b>	<b>\$2,036</b>

#### Transit Models\* (WC19 Compliant):

<input type="checkbox"/> Cruiser Transit CX10T - 903316	<b>Price</b>	<b>\$1,876</b>
<input type="checkbox"/> Cruiser Transit CX12T - 902847	<b>Price</b>	<b>\$1,876</b>
<input type="checkbox"/> Cruiser Transit CX14T - 900491	<b>Price</b>	<b>\$1,876</b>
<input type="checkbox"/> Cruiser Transit CX16T - 900173	<b>Price</b>	<b>\$2,036</b>
<input type="checkbox"/> Cruiser Transit CX18T - 902596	<b>Price</b>	<b>\$2,036</b>

#### TEXTILENE Upholstery Package

- Ocean Blue **Std**
- Navy Blue **Std**
- Purple **Std**
- Teal **Std**
- Raspberry **Std**

CX10	CX12	CX14	CX16	CX18
903463	903463	903463	903463	903463
903464	903464	903464	903464	903464
903465	903465	903465	903465	903465
903466	903466	903466	903466	903466
903467	903467	903467	903467	903467

#### Cruiser Classic Required Item

- 2-Piece Push Handle Package **Std**

CX10	CX12	CX14	CX16	CX18
904124	904125	904126	904127	<b>Std</b>

#### Transit Required Items\*

- \*Transit Option (Manufacturer Installation Required) **Price \$200**
- \*Head Support - Select Style Below:
  - Headrest Extension **\$130**
  - Contoured Firm Headrest (Planar Option Only/Hrdw Req) **\$244**
  - Adj. Removable Hardware (Req for Contoured Headrest) **\$247**
- \*5-Pt Harness **OR** \*H-Harness / 3-Pt Belt - Select Style Below:
  - 5-Point Harness with Padded Covers - CX10/12\*\* **\$143**
  - H-Harness w/ Padded Covers -CX14/16/18 (3-Pt Belt Req) **\$143**
  - 3-Pt Pos. Belt w/ Depth Adj. Crotch Strap -CX14/16/18 **\$49**
- \*Foot Positioners - Select Size Below:
  - Foot Positioners - Small 7" - 10" (pair) **\$156**
  - Foot Positioners - Medium 8" - 11" (pair) **\$156**
  - Foot Positioners - Large 9" - 12" (pair) **\$156**

#### HCPCS

CX10	CX12	CX14	CX16	CX18
904016	904016	904016	904016	904016
600041	600005	600006	600007	600008
904032	904032	904033	904034	904034
904035	904035	904035	904035	904035
903682	903683	N/A	N/A	N/A
N/A	N/A	903988	903989	903989
904017	904017	904017	904018	904018
903415	903415	903415	903415	903415
903416	903416	903416	903416	903416
903417	903417	903417	903417	903417

\* The "Transit Required Items" must be ordered with transit model as were used in WC19 Safety Standard Testing.

\*\* CX10T/12T ordered with reinforced upholstery, anatomic back, cushion, or planar option must be ordered with H-harness & 3-pt positioning belt, instead of 5-pt harne.

#### Transit Optional Items

- Q'Straint Transit Lap Belt **Price \$59**
- Sure-Lok Transit Lap Belt **Price \$59**

#### HCPCS

CX10	CX12	CX14	CX16	CX18
903094	903094	903094	903094	903094
903232	903232	903232	903232	903232

#### Frame Option

- Anatomic Back\* (Heavy Duty Reinforced Upholstery Required) **Price \$295**
- Telescoping Handles Straight - 4" (pair) **\$204**
- Telescoping Handles for Anatomic - 4" (pair) **\$204**

#### HCPCS

CX10	CX12	CX14	CX16	CX18
903320	902851	902319	902320	903096
903422	903251	903253	903255	903257
903423	903252	903254	903256	903258

\*Not Available with Planar Option

#### Planar Option Required Items (Seat & Back) \*

- Planar Seat Modification **Price \$344**
- Planar Back with Attachment (Incl 2 backrest ext. w/ hrdw) **\$319**
- Cushion - select one cushion below:
  - Support - General Use Cushion **\$122**
  - Position - Medial Thigh Support/Anti Thrust **\$347**
  - Align - Lateral Pelvis/Lateral Thigh Support **\$347**
- Solid Seat Insert (Required with Cushions) **\$119**

#### HCPCS

CX10	CX12	CX14	CX16	CX18
904019	904020	904021	904022	904023
904024	904025	904026	904027	904028
903639	903642	903645	903648	903651
903641	903644	903647	903650	903653
903640	903643	903646	903649	903652
904029	904029	904029	904029	904029

\* Planar option must be ordered with planar seat modification, planar back, cushion and solid seat insert.

**Accessories for Planar Option \***

- Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req) **Price \$137**
- Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req) **Price \$137**
- Right Support:**
  - Small = 3"x 5"
  - Medium = 3.75"x 5.5"
  - Large = 4.75"x 6"
- Left Support:**
  - Small = 3"x 5"
  - Medium = 3.75"x 5.5"
  - Large = 4.75"x 6"
- Right Adj. Removable Hardware (Req for Lateral Support) **\$178**
- Left Adj. Removable Hardware (Req for Lateral Support) **\$178**
- Contoured Firm Headrest (Hardware Required) **\$244**
- Adj. Removable Hardware (Req for Contoured Firm Headrest) **\$247**
- Flat Firm Backrest Extension (Hardware Req) **\$228**
- Adj. Removable Hardware (Req for Flat Firm Backrest Ext) **\$223**

\* Most accessories listed below are also available with the planar option.

**Pelvis Positioning**

- 3-Pt Positioning Belt with Depth Adjustable Crotch Strap **Price \$49**
- Heavy Duty Reinforced Upholstery ^ **\$268**
- Custom Seat - Cut Back 1" ^ **\$108**
- Custom Seat - Cut Back 2" ^ **\$108**
- Custom Seat Depth Extension 2" ^ **\$94**
- Custom Seat Depth Extension 3" ^ **\$94**
- Foam Wedge **\$152**
- Support - General Use Cushion \* **\$122**
- Position Cushion - Medial Thigh Support/Anti Thrust \* **\$347**
- Align Cushion - Lateral Pelvis/Lateral Thigh Support \* **\$347**
- Solid Seat Insert (Required with Cushions) **\$119**

^ Item N/A with Planar Option. \* Solid seat insert required with solid cushion. Seat cut back, seat depth extension & tray are N/A with solid cushion.

**Foot & Leg Positioning**

- Foot Positioners - Small 7" - 10" (pair) **Price \$156**
- Foot Positioners - Medium 8" - 11" (pair) **\$156**
- Foot Positioners - Large 9" - 12" (pair) **\$156**
- Angle Adjustable Footplates (pair) **\$179**
- Footplate Securement Strap **\$30**
- Medial Thigh Support-Abductor (N/A Position/Align Cushion) **\$162**
- Lateral Thigh Support-Adductor (N/A Position/Align Cushion) **\$162**
- Calf Panel\* **\$61**

\* Calf pannel cannot be used when footplate extension tube is rotated, and in the highest position.

**Trunk Positioning**

- H-Harness with Padded Covers (3-Pt Positioning Belt Req) **Price \$143**
- SA Adj. Lateral Support Single Flap w/ Scoli Strap ^ **\$137**
- SA Adj. Lateral Support Double Flap w/ Scoli Strap ^ **\$137**
- Full Torso Swing-Away Vest Support ^ **\$263**
- Reducer Seat Insert (Reduces Seat Width 2" & Depth 1") ^ \* **\$353**

^ Item not available with Planar Option. \* Item must be removed to fold the chair.

**Head Supports**

- Padded Headwings (N/A with Planar Option) **Price \$244**
- Occi Headwings (N/A with Planar Option) **\$228**
- Headrest Extension (N/A with Planar Option) **\$130**

**Tire Options**

- 11" Rear Solid Tire **Price \$155**
- 6"x2" Front Ribbed, 12.5" Rear Solid Knobby Tire **\$299**
- 8"x2" Front Ribbed, 12.5" Rear Solid Knobby Tire **\$299**
- 8"x2" Front Ribbed, 12.5" Rear Pneumatic Knobby Tire **\$299**

**Accessories**

- Headrest Cover (Canopy) **Price \$228**
- Upper Extremity Support Surface (Clear Tray) ^ **\$111**
- Adj. Removable Hardware (Required for Clear Tray) **\$227**
- Upper Extremity Spprt Surface (Height Adj Black Plastic Tray)\* **\$135**
- Adj. Removable Hardware (Req for Height Adj. Black Tray) **\$227**
- Under Seat Storage Basket **\$68**
- Travel Bag (For Wheelchair Only) **\$206**
- Utility Bag - Small **\$71**
- Utility Bag - Large **\$71**
- Saddle Bags (pair) **\$119**
- Curb Tipper (elevates front of chair to go up curbs) **\$98**
- Caster Lock for 2" Wide Front Wheels **\$94**

^ Hardware required. Item not available with Planar Option or Solid Cushion. \* Hardware required.

**HCPCS**

**E0956**  
**E0956**

**E1028**

**E1028**

**E0955**

**E1028**

**E0955**

**E1028**

**HCPCS**

**E0978**

**E1297**

**E1297**

**E1297**

**E1297**

**E2601**

**E2607**

**E2607**

**E0992**

**HCPCS**

**K0040**

**E0957**

**E0956**

**K0038**

**HCPCS**

**E0960**

**E0956**

**E0956**

**HCPCS**

**E0955**

**E0955**

**E0966**

**HCPCS**

**HCPCS**

**E0950**

**E1028**

**E0950**

**E1028**

**K0073**

CX10	CX12	CX14	CX16	CX18
Select Size for Right and Left Lateral Supports Below:				
902965	902965	902965	902965	902965
902969	902969	902969	902969	902969
902973	902973	902973	902973	902973
904030	904030	904030	904031	904031
904030	904030	904030	904031	904031
904032	904032	904033	904034	904034
904035	904035	904035	904035	904035
904089	904090	904091	904092	904093
904057	904057	904057	904057	904057

CX10	CX12	CX14	CX16	CX18
904017	904017	904017	904018	904018
903322	902853	900443	900139	902601
903525	903527	903529	903531	903533
903323	902854	900445	900198	902602
903994	903996	903998	904000	N/A
903995	903997	903999	904001	N/A
903333	902864	900654	900146	902604
903639	903642	903645	903648	903651
903641	903644	903647	903650	903653
903640	903643	903646	903649	903652
904029	904029	904029	904029	904029

CX10	CX12	CX14	CX16	CX18
903415	903415	903415	903415	903415
903416	903416	903416	903416	903416
903417	903417	903417	903417	903417
903445	903268	903268	903269	903269
903584	903584	903584	903584	903584
903326	902857	900480	900137	902607
903327	902858	900482	900188	902608
904111	904112	904113	904114	904115

CX10	CX12	CX14	CX16	CX18
903988	903988	903988	903989	903989
903328	902859	900678	900140	902614
903329	902860	900679	900204	902615
903330	902861	900680	900339	902616
904211	904211	904110	904110	N/A

CX10	CX12	CX14	CX16	CX18
903331	902862	900681	900142	902618
904076	904076	904316	904315	904315
600041	600005	600006	600007	600008

CX10	CX12	CX14	CX16	CX18
900407	900407	900407	900407	N/A
N/A	903082	N/A	N/A	N/A
N/A	N/A	903084	903084	903084
N/A	N/A	903083	903083	903083

CX10	CX12	CX14	CX16	CX18
903336	902867	900162	900184	902623
903078	903385	903386	903387	903388
903352	903353	903354	903355	903356
N/A	N/A	903696	N/A	N/A
N/A	N/A	903697	N/A	N/A
904046	904046	904047	904047	904048
902022	902022	902022	902022	902022
903227	903227	903227	903227	903227
903228	903228	903228	903228	903228
904094	904094	904094	904094	904094
903916	903916	903916	903916	903916
N/A	902920	902920	902920	902920

HCPCS codes are provided for your convenience and NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.