

# Cruiser Classic Planar Features

HCPCS:E1236

PLANAR

- 30° Fixed Tilt  
*Promotes upper body positioning for individuals with limited upper body strength or trunk stability*
- Solid Seat & Back  
*Provide added support for more in-depth seating and positioning*
- Promotes Good Trunk Posture  
*For better respiration and improved digestion*
- Comprehensive Positioning Options  
*Provide maximum comfort and optimal positioning*

- Solid Seat and Back Remove  
*For compact fold and easy transport*
- Frame Adjusts for Growth  
*To provide optimal positioning, and extend the life of the chair for a growing child*
- Rigorous Product Testing  
*Ensures maximum safety, strength, stability, and endurance*
- Crash Tested Transit Models  
*Adhere to national safety standards, and provide safe transport in bus or van*

PLANAR  
CX














compact folded  
**CRUISER  
PLANAR**

Positioning options and accessories sold separately.

Adjustable Swing-Away Footplates

# Cruiser Classic Planar Specifications



Measurement	Chair Model											
		Seat Width*	Seat Depth*	Seat Back Height	Seat to Footplate*	Seat to Floor*	Seat-to-Back Angle	Seat Angle	Headrest Extension	Shoulder Strap Height	Weight Capacity Regular/Transit	Weight of Chair <sup>Δ</sup>
in.	<b>CX 10</b>	10	8 - 10	18	7 - 22	25	85° - 90° - 95°	30°	5	12 - 18	75 lbs / 66 lbs	19 lbs
cm	<b>CX 10</b>	25.5	20 - 25.5	45.5	18 - 56	63.5			13	30.5 - 45.5	34 kg / 30 kg	8.6 kg
in.	<b>CX 12</b>	12	8 - 10	18	7 - 22	26	85° - 90° - 95°	30°	4.5	11 - 18	75 lbs / 66 lbs	19 lbs
cm	<b>CX 12</b>	30.5	20 - 25.5	45.5	18 - 56	66			11.4	28 - 45.5	34 kg / 30 kg	8.6 kg
in.	<b>CX 14</b>	14	11 - 14	29.5	6.5 - 26	26	85° - 90° - 95°	30°	5	15 - 24	100 lbs / 100 lbs	21 lbs
cm	<b>CX 14</b>	35.5	28 - 36	75	16.5 - 66	66			13	38 - 61	45.5 kg / 45.5 kg	9.5 kg
in.	<b>CX 16</b>	16	11 - 14	26	12 - 24	24	85° - 90° - 95°	30°	8	15.5 - 22.5	170 lbs / 170 lbs	24.5 lbs
cm	<b>CX 16</b>	40.5	28 - 35.5	66	30.5 - 61	61			20	39 - 57	77 kg / 77 kg	11 kg
in.	<b>CX 18</b>	18	14 - 17	26	13 - 25	24.5	85° - 90° - 95°	30°	8	16 - 22	250 lbs / 170 lbs	30 lbs
cm	<b>CX 18</b>	45.5	35.5 - 43	66	33 - 63.5	62			20	40.5 - 56	114 kg / 77 kg	13.6 kg







in. - Inches, cm - Centimeters \* Plus or minus .5 inch (1.27 cm) due to fabric. Δ - Weight of chair is without seating system and front rigging.

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.

Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

## Upholstery Colors



Measurement	Chair Model						
		Overall Width	Overall Height	Overall Length	Folded Length	Folded Width	Folded Height
in.	<b>CX 10</b>	19	36	43	46	13	14
cm	<b>CX 10</b>	48	91	109	117	33	36
in.	<b>CX 12</b>	21	41	45	46	13	14
cm	<b>CX 12</b>	53	104	114	117	33	35.5
in.	<b>CX 14</b>	23	45	47	46	14	14
cm	<b>CX 14</b>	58	114	119.4	117	36	36
in.	<b>CX 16</b>	26	37.5	39.5	48	14	15
cm	<b>CX 16</b>	66	95	100	122	36	38
in.	<b>CX 18</b>	28	38	40	48	15	16.5
cm	<b>CX 18</b>	71	96.5	102	122	38	42

## Limited Lifetime Warranty

Frame and X-Braces.....Lifetime of original retail buyer  
 Fabric and webbing.....One year  
 Other components.....One year

Planar Seat & Back.....Five years  
 Other components.....18 months

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.

