

Cruiser Planar Features

HCPCS:E1236

PLANAR

- 30° Fixed Tilt
Promotes upper body positioning for individuals with limited upper body strength or trunk stability
- Solid Seat & Back
Provide added support for more in-depth seating and positioning
- Promotes Good Trunk Posture
For better respiration and improved digestion
- Comprehensive Positioning Options
Provide maximum comfort and optimal positioning

- Solid Seat and Back Remove
For compact fold and easy transport
- Frame Adjusts for Growth
To provide optimal positioning, and extend the life of the chair for a growing child
- Rigorous Product Testing
Ensures maximum safety, strength, stability, and endurance
- Crash Tested Transit Models
Adhere to national safety standards, and provide safe transport in bus or van

PLANAR

CX



Standard Solid Tires
CX10/12/14/16: 7.5" x 2" Front, 11.5" x 2.5" Quick Release Rear Wheels
CX18: 8" Front, 11" Rear

Positioning options and accessories sold separately.

Cruiser Planar Specifications

| Measurement | Chair Model | Seat Width* | Seat Depth* | Seat Back Height | Seat to Footplate* | Seat to Floor* | Seat-to-Back Angle | Seat Angle | Headrest Extension | Shoulder Strap Height | Weight Capacity Regular/Transit | Weight of Chair ^Δ |
|-------------|--------------|-------------|-------------|------------------|--------------------|----------------|--------------------|------------|--------------------|-----------------------|---------------------------------|------------------------------|
| in. | CX 10 | 10 | 9 - 11 | 18 | 5 - 28 | 27.5 | 85° - 90° - 95° | 30° | 5 | 12 - 18 | 75 lbs / 66 lbs | 36 lbs |
| cm | CX 10 | 25.5 | 23 - 28 | 45.5 | 12.5 - 71 | 70 | | | 12.5 | 30.5 - 45.5 | 34 kg / 30 kg | 16 kg |
| in. | CX 12 | 12 | 8 - 10 | 18 | 7 - 22 | 26 | 85° - 90° - 95° | 30° | 4.5 | 11 - 18 | 75 lbs / 66 lbs | 36 lbs |
| cm | CX 12 | 30.5 | 20 - 25.4 | 45.5 | 18 - 56 | 66 | | | 11.4 | 28 - 45.5 | 34 kg / 30 kg | 16 kg |
| in. | CX 14 | 14 | 11 - 14 | 29.5 | 6.5 - 26 | 26 | 85° - 90° - 95° | 30° | 5 | 13 - 24 | 100 lbs / 100 lbs | 38 lbs |
| cm | CX 14 | 35.5 | 28 - 36 | 75 | 16.5 - 66 | 66 | | | 13 | 33 - 61 | 45.5 kg / 45.5 kg | 17 kg |
| in. | CX 16 | 16 | 13 - 16 | 27 | 7 - 25 | 26 | 85° - 90° - 95° | 30° | 8 | 15.5 - 22.5 | 170 lbs / 170 lbs | 45 lbs |
| cm | CX 16 | 40.5 | 33 - 40.5 | 68.5 | 17.5 - 63.5 | 66 | | | 20 | 39 - 57 | 77 kg / 77 kg | 20.5 kg |
| in. | CX 18 | 18 | 14 - 17 | 26 | 13 - 25 | 24.5 | 85° - 90° - 95° | 30° | 8 | 16 - 22 | 250 lbs / 170 lbs | 30 lbs |
| cm | CX 18 | 45.5 | 35.5 - 43 | 66 | 33 - 63.5 | 62 | | | 20 | 40.5 - 56 | 114 kg / 77 kg | 13.6 kg |

in. - Inches, cm - Centimeters * Plus or minus .5 inch (1.27 cm) due to fabric. Δ - Weight of chair is without seating system and front rigging.

| Measurement | Chair Model | Overall Width | Overall Height | Overall Length | Push Handle Adjustments [°] | Folded Length | Folded Width | Folded Height |
|-------------|--------------|---------------|----------------|----------------|--------------------------------------|---------------|--------------|---------------|
| in. | CX 10 | 21.25 | 37.5 | 53.5 | 28 - 43 | 47.5 | 17 | 16.5 |
| cm | CX 10 | 54 | 95 | 136 | 71 - 109 | 120.5 | 43 | 42 |
| in. | CX 12 | 23.25 | 38 | 44 | 29 - 46 | 48 | 15 | 17 |
| cm | CX 12 | 59 | 96.5 | 112 | 74 - 117 | 122 | 38 | 43 |
| in. | CX 14 | 25.5 | 40 | 49 | 28 - 45 | 47 | 16 | 17 |
| cm | CX 14 | 65 | 102 | 124.5 | 71 - 114 | 119.4 | 41 | 43 |
| in. | CX 16 | 27.5 | 42.5 | 60 | 31 - 48 | 48.5 | 19 | 16 |
| cm | CX 16 | 70 | 108 | 152.5 | 78 - 122 | 123 | 48 | 40.5 |
| in. | CX 18 | 28 | 38 | 40 | N/A | 48 | 15 | 16.5 |
| cm | CX 18 | 71 | 96.5 | 02 | N/A | 122 | 38 | 42 |

[°] Adjustable Push Handle is not available on 18 inch models

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.


Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

Limited Lifetime Warranty

Frame and X-Braces.....Lifetime of original retail buyer
Fabric and webbing.....One year
Other components.....One year

Planar Seat & Back.....Five years
Other components.....18 months

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.



Upholstery Colors

