

CuddleBug Tilt Features

HCPCS:E1232

C
U
D
D
L
E
B
U
G

- 5° - 30° of Adjustable Tilt
Provides pressure relief and upper body support for individuals with insufficient trunk stability
- 90° - 170° of Adjustable Recline
Opens the seat-to-back angle for pressure relief, better respiration, improved digestion, and sleeping comfort
- 7" of Seat-to-Floor Adjustment (High/Low)
Changes height of chair for interaction with peers and table access
- Early Intervention Wheelchair
Addresses the positioning needs of children from 6 months - 4 years
- Tilt-in-Space Base with Unique Compact Fold
Folds compactly for easy transport and storage - fold and go!
- In-Depth Positioning
Provide maximum comfort and optimal positioning
- New Wider Wheels
Provide added stability and a comfortable ride on uneven terrain
- Quick Release Rear Wheels
Remove for transport and storage
- WC19 Crash Tested Transit Models
Adhere to national safety standards, and provide safe transport in bus or van



Standard Solid Tires

CB 1/2: 7.5" x 2" Front, 11.5" x 2.5" Quick Release Rear Wheels

Positioning options and accessories sold separately.




compact folded
CuddleBug

CuddleBug Tilt Specifications

CuddleBug Standard Features

- 5°-30° of adjustable tilt
- Depth adjustable seat with foam hip inserts
- Removable / washable seat cover
- Breathable padded seat insert
- One-piece angle adjustable footplate
- Elevating legrest
- Headrest cover (canopy)
- High /Low settings – 7” of seat-to-floor adjustment
- Silver powder coated frame
- One-piece height adjustable push handle
- Under-seat storage tray
- Swivel front wheels
- 7.5” x 2” front, 11.5” x 2.5” rear solid wheels
- Quick release rear wheels
- Wheel locks
- Calf panel
- 2-pt positioning belt

Upholstery Colors



Due to manufacturer's processes, we cannot guarantee the color of the upholstery will match from lot to lot.

CuddleBug Specifications

Small Module CB 1	Large Module CB 2
Seat width: 6" – 9"	Seat width: 9" – 12"
Headrest adjustment: 10.5" - 23"	Headrest adjustment: 11" - 25"
Maximum seat width (without hip inserts): 9"	Maximum seat width (without hip inserts): 12"
Width between laterals (min/max): 4" - 8"	Width between laterals (min/max): 6.5" - 10.5"
Seat depth: 5" – 8"	Seat depth: 7" – 10"
Seat back height: 17"	Seat back height: 17"
Seat to footplate: 4" – 12"	Seat to footplate: 4" – 12"
Seat to floor: 16" – 23"	Seat to floor: 16" – 23"
Tilt adjustment: 5° - 30°	Tilt adjustment: 5° - 30°
Recline adjustment: 90° - 170°	Recline adjustment: 90° - 170°
Weight capacity standard: 35 lbs	Weight capacity standard: 66 lbs
Weight capacity transit: 35 lbs	Weight capacity transit: 66 lbs
Weight of chair: 41 lbs	Weight of chair: 44 lbs
Folded Dimension	Folded Dimension
Folded length: 32 1/2"	Folded length: 33 1/2"
Folded width: 22"	Folded width: 24"
Folded height: 17 3/4" - 19 3/4" (± 1/2")	Folded height: 17 3/4" - 19 3/4" (± 1/2")

Limited Lifetime Warranty

Frame and X-Braces.....Lifetime of original retail buyer
 Fabric and webbing.....One year
 Other components.....One year

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.

