

# EZ Rider Planar Features

HCPCS:E1236

EZ RIDER PLANAR

- 10° Fixed Tilt  
*Promotes an upright position for individuals with good trunk strength and upper body stability*
- Solid Seat & Back  
*Provide added support for more in-depth seating and positioning*
- Planar Seating Promotes Good Trunk Posture  
*For better respiration and improved digestion*
- Solid Seat and Back Remove  
*For compact fold and easy transport*
- Comprehensive Positioning Options  
*Provide maximum comfort and optimal positioning*
- Frame Adjusts for Growth  
*To provide optimal positioning, and extend the life of the chair for a growing child*
- "Quick Release" Rear Wheels  
*Easily remove for storage or transport*
- Crash Tested Transit Models  
*Adhere to national safety standards, and provide safe transport in bus or van*














Standard Solid Tires

**EZ12/14/16:** 7.5" x 2" Front, 11.5" x 2.5" Quick Release Rear Wheels  
**EZ18:** 8" Front, 11" Rear

Positioning options and accessories sold separately.

compact folded  
**EZ RIDER  
PLANAR**

# EZ Rider Planar Specifications


Measurement	Chair Model											
		Seat Width*	Seat Depth**	Seat Back Height	Seat to Footplate*	Seat to Floor*	Back Angle	Seat Angle	Headrest Extension	Shoulder Strap Height	Weight Capacity Regular/Transit	Weight of Chair $\Delta$
in.	<b>EZ 12</b>	12	11-14	16.5	5- 20.5	21.5	100°	10°	5	12 - 21	75 lbs / 66 lbs	32 lbs
cm	<b>EZ 12</b>	30.5	28 - 35.5	42	12.5 - 52	54.5			12.5	30.5 - 53.3	34 kg / 30 kg	14.5 kg
in.	<b>EZ 14</b>	14	11-14	19	6.5-23	23	100°	10°	5	13 - 24	100 lbs / 100 lbs	35 lbs
cm	<b>EZ 14</b>	35.5	28 - 35.5	48	16.5 - 58.5	58.5			12.5	33 - 61	45.5 kg / 45.5 kg	16 kg
in.	<b>EZ 16</b>	16	11-14	20	6.5-25	25	100°	10°	6	14.5 - 25	170 lbs / 170 lbs	37 lbs
cm	<b>EZ 16</b>	40.5	28 - 35.5	51	16.5 - 63.5	63.5			15	36.8 - 63.5	77 kg / 77 kg	17 kg
in.	<b>EZ 18</b>	18	16 - 19	28	4.5 - 24.5	24.5	100°	10°	6	17 - 28	250 lbs / 170 lbs	49 lbs
cm	<b>EZ 18</b>	45.5	40.6 - 48.3	71	11.4 - 62.2	62.2			15	43.2 - 71	114 kg / 77 kg	22.2kg







in. - Inches, cm - Centimeters \* Plus or minus .5 inch (1.27 cm) due to fabric. \*\* 2 inch cut-out for minimum seat depth. Available as an option.  $\Delta$  - Weight of chair is without front rigging.


Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.

Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

## Upholstery Colors



Measurement	Chair Model						
		Overall Width	Overall Height	Overall Length	Folded Length	Folded Width	Folded Height
in.	<b>EZ 12</b>	21.5	38	43	45.5	14	18
cm	<b>EZ 12</b>	54.5	96.5	109	115.5	35.5	45.5
in.	<b>EZ 14</b>	23	38	46	45.5	13.5	17
cm	<b>EZ 14</b>	58.5	96.5	117	115.5	34	43
in.	<b>EZ 16</b>	26	36	49	45.5	14	16
cm	<b>EZ 16</b>	66	91.5	124.5	115.5	35.5	40.5
in.	<b>EZ 18</b>	25	48	47.5	44.5	14.5	15.25
cm	<b>EZ 18</b>	63.5	122	120.7	113	36.8	38.7



## Limited Lifetime Warranty

Frame and X-Braces.....Lifetime of original retail buyer  
Fabric and webbing.....One year  
Other components.....One year

Planar Seat & Back.....Five years  
Other components.....18 months

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.

