

Lite Rider Features












LITE RIDER

- **10° Fixed Tilt**
Provides an upright position for individuals with good trunk strength and upper body stability
- **Ultra-Lightweight Wheelchair**
Is easy to lift and carry
- **Compact Folding Frame**
Easily stows in a car, bus or airplane
- **Comprehensive Positioning Options**
Address all positioning needs
- **Size Adjustable Frame**
Extends the life of the chair for a growing child
- **Rigorous Product Testing**
Ensures maximum safety, strength, stability, and endurance
- **Removable, Washable Fabric**
Is Anti-Bacterial and Fire Retardant









Positioning options and accessories sold separately.

Lite Rider Specifications

Measurement	Chair Model	 Seat Width*	 Seat Depth**	 Seat Back Height	 Seat to Footplate*	 Seat to Floor*	 Seat-to-Back Angle	 Seat Angle	 Headrest Extension	 Shoulder Strap Height	 Weight Capacity	 Weight of Chair Δ
in.	LR 12	12	10 - 13	17.5	2.5 - 19	17	100°	10°	5	11 - 16	75 lbs	16 lbs
cm	LR 12	30.5	25 - 33	44	6 - 48	43			13	28 - 41	34 kg	7 kg
in.	LR 14	14	11.5 - 14.5	20	1.75 - 19	19	100°	10°	8	11 - 18	100 lbs	18 lbs
cm	LR 14	35.5	29 - 37	51	4 - 48	48			20	28 - 46	45.5 kg	8 kg

in. - Inches, cm - Centimeters * Plus or minus .5 inch (1.27 cm) due to fabric. ** 2 inch cut-out for minimum seat depth. Available as an option. Δ - Weight of chair is without front rigging.

Measurement	Chair Model	 Overall Width	 Overall Height	 Overall Length	 Folded Length	 Folded Width	 Folded Height
in.	LR 12	19.5	36.5	43	43	13.5	14.5
cm	LR 12	50	93	109	109	34	37
in.	LR 14	20.5	37	44	44	12	14.5
cm	LR 14	52	94	112	112	30	37

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.


Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

Limited Lifetime Warranty

Frame and X-Braces.....Lifetime of original retail buyer
 Fabric and webbing.....One year
 Other components.....One year

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.

Upholstery Colors



U.S. Patent 4,118,065
 This product may be covered by the
 following U.S. patent - 6,099,075
 Additional patents pending.

visit us at www.convoid.com