

Mountee Features

HCPCS Codes
 Base: E1234
 Planar Seat: E2292
 Planar Back: E2291

E
E
T
N
D
O
M

- 5° - 45° of Adjustable Tilt
Provides pressure relief and repositions individuals with insufficient upper body strength or trunk stability
- Seating System (Flat or Contoured)
Provides comfort and postural support during extended use
- Elevating Legrests (option)
Allow for changes in position and promote better circulation
- Seating System Removes from Base
For compact fold - easy to transport and store

- Comprehensive Positioning Options
Provide maximum comfort and optimal positioning
- Frame Adjusts for Growth
To provide optimal positioning, and extend the life of the chair for a growing child
- New Wider Wheels
Provide added stability and a comfortable ride on uneven terrain
- Crash Tested Transit Model (shown)
Adhere to national safety standards, and provide safe transport in bus or van



One-Piece Footplate (standard)



compact folded
MOUNTEE

MT12/14/: 7.5" x 2" Front, 11.5" x 2.5" Quick Release Rear Wheels

Positioning options and accessories sold separately.

Mountee Specifications

Measurement	Chair Model	Frame Size	Seat Width*	Lateral Support Width Adjustment	Lateral Support Height Adjustment	Seat Depth	Seat Back Height	Seat to Footplate*	Seat to Floor*	Armrest to Seat	Seat Angle (from horizontal)	Seat-to-Back Angle	Tilt Adjustment (from vertical)
in.	MT 12	14	8 - 12	6 - 12	9 - 11.75	11 - 15	17	10 - 16	26	7.75 - 10.75	5°	90°	5° - 45°
cm	MT 12	35.5	20 - 30.5	15 - 30.5	23 - 30	28 - 38	43	25.5 - 40.5	66	19.7 - 27.3	-----	-----	-----
in.	MT 14	16	8 - 14	7 - 14	10 - 12.5	13 - 17	18.5	10 - 17	26	7.75 - 10.25	5°	90°	5° - 45°
cm	MT 14	40.6	20 - 35.5	18 - 35.5	25.5 - 32	33 - 43	47	25 - 43	66	19.7 - 26	-----	-----	-----

in. - Inches, cm - Centimeters * Plus or minus .5 inch (1.27 cm) due to fabric.

Measurement	Chair Model	Firm Backrest Extension	Shoulder Strap Height	Pommel Height	Overall Width	Overall Height	Overall Length	Folded Length	Folded Width	Folded Height	Weight Capacity Standard/Transit	Weight of Seating System	Weight of Base ^Δ
in.	MT 12	5	17	2.5	24.25	45.5	39	50	14.75	14	75/66 lbs	17 lbs	32 lbs
cm	MT 12	12.5	43	6	61.5	115.5	99	122	37.5	35.5	34/30 kg	7.7 kg	14.5 lbs
in.	MT 14	5	18.5	2.5	25.5	45.5	42.5	51	16.25	14	100/100 lbs	18.5 lbs	36.5 lbs
cm	MT 14	12.5	47	6	65	115.5	108	121.5	41	35.5	45/45 kg	8.4 kg	16.6 lbs

^Δ - Weight of Base is without front rigging.

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.

Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

Limited Lifetime Warranty

Frame and X-Braces.....Lifetime of original retail buyer
 Fabric and webbing.....One year
 Other components.....One year

Seating System.....Five years hardware
 Other components.....18 months

Occi-Flex Headrest.....Five years
 Other components.....One year

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.

Upholstery Colors

