

# Cruiser Profiler Features

HCPCS: E1238

# PROFILER

- 30° Fixed Tilt  
*Promotes upper body positioning for individuals with limited upper body strength or trunk stability*
- Adjustable Strap Back System  
*Provides customized positioning – straps conform to the shape of the back and address rotation issues*
- Strap System Easily Readjusts  
*For growth, postural correction and to alleviate painful pressure points*
- Frame Folds Compactly  
*For easy storage and transport*
- In-Depth Positioning Options  
*Provide optimal positioning for increased therapeutic benefits – all day seating*
- Multiple Frame Adjustments  
*Extend the life of the chair for a growing child*
- Crash Tested Transit Models  
*Adhere to national safety standards, and allow safe transport in bus or van*

Occi-Flex Cushioned Headrest

Anterior Trunk Support


3-Point Positioning Belt with Depth Adjustable Crotch Strap

Height Adjustable Flip-Up Armrests

Anti-Thrust Solid Seat Cushion

Toggle Wheel Locks

Adjustable Swing-Away Footplates



compact folded  
**PROFILER**

# Profiler Specifications

Measurement	Chair Model	Frame Size	Seat Width*	Seat Width with Hip Guide	Strap Contour Depth	Seat Depth	Seat Back Height	Armrest Height	Seat to Footplate*	Seat to Floor*	Seat Angle	Shoulder Strap Height
in.	<b>PR 14</b>	16	14	12	7	10.5 - 19.5	21	5 - 8	26	25	30°	21.5 - 20.5
cm	<b>PR 14</b>	40.5	35.5	30.5	17.8	27 - 40	53	13 - 20	66	63.5		31.75 - 52.07
in.	<b>PR 16</b>	18	16	14	7	12 - 19	23	6 - 9	26	26.5	30°	11.75 - 21.75
cm	<b>PR 16</b>	45.5	40.5	35.5	17.8	30.5 - 48	58.5	15 - 23	66	67.3		29.84 - 55.3

in. - Inches, cm - Centimeters \* Plus or minus .5 inch (1.27 cm) due to fabric. \*\* 2 inch cut-out for minimum seat depth. Available as an option. Δ - Weight of chair is without front riggings.



## Adjustable Strap Back System

The Profiler's back support straps individually adjust to provide a custom fit. Because each strap is adjustable, the resulting support will accommodate almost any conceivable contour requirement. Straps are covered with a padded seat cover for the user's comfort.

Measurement	Chair Model	Overall Width	Overall Height	Overall Length	Folded Length	Folded Width	Folded Height	Weight Capacity Regular/Transit	Weight of Chair
in.	<b>PR 14</b>	25	38	47	48	17	17	100 / 100 lbs	51 lbs
cm	<b>PR 14</b>	63.5	96.5	119.5	122	43	43	45.5 / 45.5 kg	23 kg
in.	<b>PR 16</b>	29	41.5	51	53	20	16	170 / 170 lbs	55 lbs
cm	<b>PR 16</b>	73.5	105.5	129.5	135	51	40.5	77 / 77 kg	25 kg

in. - Inches, cm - Centimeters

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.

Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

## Limited Lifetime Warranty

Frame and X-Braces.....Lifetime of original retail buyer  
Fabric and webbing.....One year  
Other components.....One year

Occi-Flex Headrest.....Five years  
Other components.....One year

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.

## Upholstery Colors

Removable, Washable Fabric

Panther Black	True Red	Purple	Nile Green
Forest Green	Mediterranean Blue	Deep Burgundy	Pink

