

R82 Pony - E8002

USD Retail Price List/ Order Form


Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email	

US Configuration

Effective: 05.01.21 - Rev 3

- Pony
- Head Support
- Abdominal Pad
- Hip Support
- Arm Support



A (Total Height) _____
B (Shoulder Height) _____
C (Chest Height) _____
D (Buttock Height) _____
E (Chest Width) _____
F (Shoulder Width) _____
 Age _____
 Weight _____

Measurements



	Size 0 Inches	Size 1 Inches	Size 2 Inches
Seat height (L)	9¼ - 14¾"	14¾ - 19½"	18¾ - 25"
Width (G)	23"	27"	30"
Length (H)	26¼"	28"	40½"
Height (I)	20¼ - 27"	23½ - 32"	26½ - 40½"
Angle range (A)	32°	32°	32°
Weight	15½ lbs.	22 lbs.	30 lbs.
Max. load	66 lbs.	88 lbs.	110 lbs.


Pony Standard Features



- Rear Drag Wheels with Directional Stabilizer
- Leg & Foot Separator
- Height Adjustable Seat with Spring
- Height Adjustable Removable Back Support
- Height Adjustable Chest Support
- 5 Wheels (1 Rear & 2 Front with Locks)
- Height & Depth Adjustable Handle Bar

NOTES:

Pony	Description	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Pony Walking Aid, Yellow	0	86400	\$2,383	E8002
	<input type="checkbox"/> Pony Walking Aid, Purple	1	86401	\$2,383	E8002
	<input type="checkbox"/> Pony Walking Aid, Lilac	2	86402	\$2,383	E8002

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Hip & Thigh Supports (Pr.)	6" x 3"	0	86408	\$338	
	<input type="checkbox"/> Hip & Thigh Supports (Pr.)	6" x 4"	1	86410	\$338	
	<input type="checkbox"/> Hip & Thigh Supports (Pr.) ••	6" x 4"	2	86412	\$338	

••Custom Order; Approx. lead time: 6-8 weeks

Minimum/Maximum Width	Size 0	Size 1	Size 3
	6" x 8¾"	6" x 8¾"	6" x 10¾"

Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Abdominal Cushion (Ea.) ••		1	86420	\$200	
<input type="checkbox"/> Abdominal Cushion (Ea.) ••		2	86422	\$200	

••Custom Order; Approx. lead time: 6-8 weeks




<input type="checkbox"/> Armrest/Tray (Ea.)	0	86427	\$497
<input type="checkbox"/> Armrest/Tray (Ea.)	1	86428	\$497
<input type="checkbox"/> Armrest/Tray (Ea.) ••	2	86429	\$497

••Custom Order; Approx. lead time: 6-8 weeks



<input type="checkbox"/> Vertical Handle for Mounting in the Middle (Pr.)	One Size	86432	\$118
---	----------	-------	-------

Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Head Support No. 1, Black PU		One Size	6171550	\$210	
<input type="checkbox"/> Head Support No. 2, Black PU		One Size	6172550	\$210	
<input type="checkbox"/> Head Support No. 3, Black PU		One Size	6173550	\$210	



<input type="checkbox"/> Anatomic Head Support ♦	A:5½ x B:5½ x C:9½"	Small	5099001	\$261	E0955
<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	N/C	
<input type="checkbox"/> Anatomic Head Support ♦	A:6 x B:5½ x C:10"	Medium	5099002	\$261	E0955
<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	N/C	
<input type="checkbox"/> Anatomic Head Support ♦	A:6½ x B:5½ x C:12½"	Large	5099003	\$261	E0955
<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	N/C	

♦ Turtle Bar System Short or Long Required; 98619-25 or 98619-26

<input type="checkbox"/> Headband System for Anatomic Head Support Only	One Size	99500	\$236	E0955
---	----------	-------	-------	-------



Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------








Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable

<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$236	E1028
<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$236	E1028

♦ Required with Anatomic Head Support; 5099001, 5099002 or 5099003

<input type="checkbox"/> Connecting Rod for Head Support 7¾"	1	861424-3	\$264
<input type="checkbox"/> Connecting Rod for Head Support 13"	2	861424-4	\$264

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Leg Protection		One Size	86450	\$226	
	<input type="checkbox"/> Push Handle		One Size	86416	\$261	

Spare Parts	Description	Measurements	Size	Part No.	USD Price	Qty
	<input type="checkbox"/> Chest Support, Complete		0	86310	\$523	_____
	<input type="checkbox"/> Chest Support, Complete		1	86311	\$523	_____
	<input type="checkbox"/> Chest Support, Complete		2	86312	\$523	_____
	<input type="checkbox"/> Straps for Chest Support		0	5540250	\$37	_____
	<input type="checkbox"/> Straps for Chest Support		1	5540300	\$47	_____
	<input type="checkbox"/> Straps for Chest Support		2	5540350	\$57	_____
	<input type="checkbox"/> Handle		0	86320	\$118	_____
	<input type="checkbox"/> Handle		1	86321	\$118	_____
	<input type="checkbox"/> Handle		2	86322	\$118	_____
	<input type="checkbox"/> Seat		0	86330	\$256	_____
	<input type="checkbox"/> Seat		1	86331	\$256	_____
	<input type="checkbox"/> Seat		2	86332	\$256	_____
	<input type="checkbox"/> Back Support		0	86434	\$210	_____
	<input type="checkbox"/> Back Support		1, 2	86435	\$210	_____
	<input type="checkbox"/> Directional Stabilizer for Rear Wheel		0	86438	\$89	_____
	<input type="checkbox"/> Directional Stabilizer for Rear Wheel		1	86440	\$89	_____
	<input type="checkbox"/> Directional Stabilizer for Rear Wheel		2	86442	\$89	_____
	<input type="checkbox"/> Directional Stabilizer for Front Wheel		0, 1, 2	86448	\$379	_____
	<input type="checkbox"/> Drag Wheel		One Size	86418	\$84	_____

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Prices in USD & is subject to change without notice
For Quotes and Orders please send completed Order Form to: us.info@r82.com

Etac
13137 Bleinheim Lane
Matthews, North Carolina 28105
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF5900-ENG-REV.03

Creating Possibilities