



UPRIGHT EZ RIDER

Standard Specifications



















MEASUREMENTS:

INCHES

HCPCS: E1236

- * Plus or minus .5 inches due to fabric
- △ Weight of chair is without front rigging
- ^ Seat back height
 - (C) Cordura upholstery
 - (T) Textilene upholstery
- Adjustable push handle is not available on 18 inch models

Note: Convaid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

	EZ 12	EZ 14	EZ 16	EZ 18
Seat Width* 	12	14	16	18
Seat Depth* 	9 - 14	11 - 16	12 - 17	14.5 - 19.5
Seat Back Height^ 	21C/19T	24C/22T	26C/24T	27C/25T
Seat to Floor 	19	20	20	21.5
Seat-to-Back Angle 	100°	100°	100°	100°
Seat Angle 	10°	10°	10°	10°
Overall Height 	38	38.5	40	37
Overall Length 	36	37.5	43.5	48
Overall Width 	21.5	23	26	26
Folded Length 	44.5	46	48	45.5
Folded Width 	14	13.5	14	14.5
Folded Height 	17.5	18	17	15.5
Headrest Extension 	8	8	8	10
Shoulder Strap Height 	11 - 17	14 - 20	15 - 21	16 - 22
Seat to Footplate 	3 - 20	3 - 21	3 - 22	3 - 22
Weight Capacity Regular/Transit 	75/66 lbs	100/100 lbs	170/170 lbs	250/170 lbs
Weight of Chair △ 	25 lbs	30 lbs	30 lbs	29 lbs
Push Handle Adjustment * 	29.5 - 46.5	29 - 47	28 - 45	N/A



UPRIGHT EZ RIDER

Standard Specifications

MEASUREMENTS:

CENTIMETERS

HCPCS: E1236

- * Plus or minus .5 inches due to fabric
- △ Weight of chair is without front rigging
- ^ Seat back height
- (C) Cordura upholstery
- (T) Textilene upholstery
- Adjustable push handle is not available on 18 inch models

Note: Convaid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

	EZ 12	EZ 14	EZ 16	EZ 18
Seat Width* 	30.5	35.5	40.5	45.5
Seat Depth* 	23 - 36	28 - 41	30.5 - 43	37 - 49.5
Seat Back Height^ 	53C/48T	61C/56T	66C/61T	68C/63.5T
Seat to Floor 	48	51	51	55
Seat-to-Back Angle 	100°	100°	100°	100°
Seat Angle 	10°	10°	10°	10°
Overall Height 	96.5	98	101.5	94
Overall Length 	91.5	95	110.5	122
Overall Width 	54.6	54.8	66	66
Folded Length 	113	117	122	116
Folded Width 	36	34	36	37
Folded Height 	44.5	46	43	39
Headrest Extension 	20	20	20	25
Shoulder Strap Height 	28 - 43	36 - 51	38 - 53	41 - 56
Seat to Footplate 	8 - 51	8 - 53	8 - 56	8 - 56
Weight Capacity Regular/Transit 	34/30 kg	45.5/45.5 kg	77/77 kg	114/77 kg
Weight of Chair △ 	11 kg	14 kg	14 kg	13 kg
Push Handle Adjustment * 	75 - 118	73.5 - 119	71 - 114	N/A

