



# TILT-IN-SPACE RODEO

## Standard Specifications

MEASUREMENTS:

INCHES

### HCPCS: E1234

\* Plus or minus .5 inches due to fabric

\*\* Reduced height when footplates are flipped

Δ Weight of chair is without front rigging

^ One Piece Angle Adjustable Footplate

^^ Individual Angle Adjustable Footplates

Note: Convaid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

	RD 10	RD 12	RD 14	RD 16
Seat Width*	10	12	14	16
Seat Depth*	9.5 - 13.5	11 - 15	13 - 17	15 - 19
Seat Back Height	21	22	24	25
Seat to floor	20	21	22	22
Recline Adjustment	90° or 110°	90° or 110°	90° or 110°	90° or 110°
Tilt Adjustment	5° - 45°	5° - 45°	5° - 45°	5° - 45°
Overall Height	43	43	44	45
Overall Length	30.5	35	35	36
Overall Width	21	23	25	27
Folded Length	33.5	34	34	40
Folded Width	21	23	25	27
Folded Height	17.5	17	18	18
Headrest Extension	6	8	8	8
Shoulder Strap Height	13 - 19	15 - 21	16 - 22	18 - 24
Footplate Adjustment Angle ^	±20°	±20°	±20°	±20°
Seat to Footplate^	9 - 14	9 - 14	11 - 16	13 - 18
	5.25 - 11.25**	5.25 - 11.25**	7.25 - 13.25**	9.25 - 15.25**
Seat to Footplate^^	7 - 13	7 - 13	9.5 - 15.5	12 - 18
	4.5 - 9.25**	4.5 - 9.25**	4.25 - 11**	7.25 - 13.25**
Weight Capacity Regular/Transit	75/66 lbs	75/66 lbs	100/100 lbs	170/130 lbs
Weight of Chair Δ	29 lbs	33 lbs	37 lbs	39 lbs
Push Handle Adjustment	32 - 48.5	34.5 - 51	35 - 52	36 - 53.5



# TILT-IN-SPACE RODEO

## Standard Specifications



















MEASUREMENTS:

CENTIMETERS

### HCPCS: E1234

- \* Plus or minus .5 inches due to fabric
- \*\* Reduced height when footplates are flipped
- △ Weight of chair is without front rigging
- ^ One Piece Angle Adjustable Footplate
- ^^ Individual Angle Adjustable Footplates

Note: Convaid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

	RD 10	RD 12	RD 14	RD 16
Seat Width* 	25.5	30.5	35.5	40.5
Seat Depth* 	24 - 34	28 - 38	33 - 43	38 - 48
Seat Back Height 	53	56	61	63.5
Seat to floor 	51	53	56	56
Recline Adjustment 	90° or 110°	90° or 110°	90° or 110°	90° or 110°
Tilt Adjustment 	5° - 45°	5° - 45°	5° - 45°	5° - 45°
Overall Height 	109	109	112	114
Overall Length 	77.5	89	89	91.5
Overall Width 	53	58.5	63.5	68.5
Folded Length 	85	86	86	101.5
Folded Width 	53	58.5	63.5	68.5
Folded Height 	44.5	43	46	46
Headrest Extension 	15	20	20	20
Shoulder Strap Height 	33 - 48	38 - 53	41 - 56	46 - 61
Footplate Adjustment Angle ^ 	±20°	±20°	±20°	±20°
Seat to Footplate^ 	23 - 36 13 - 28.5**	23 - 36 13 - 28.5**	28 - 41 18 - 34**	30.5 - 46 23.5 - 39**
Seat to Footplate^^ 	18 - 33 11.5 - 23.5**	18 - 33 11.5 - 23.5**	24 - 39.5 11 - 28**	30.5 - 46 18.5 - 34**
Weight Capacity Regular/Transit 	34/30 kg	34/30 kg	45.5/45.5 kg	77/58.9 kg
Weight of Chair △ 	13 kg	15 kg	17 kg	18 kg
Push Handle Adjustment 	81 - 123	88 - 129	89 - 132	93 - 136

