

# Convoid Carrot 3 Booster Seat - Child Restraint System

## Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_

Contact: \_\_\_\_\_ ATP/RTS: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

### Carrot 3 Booster Seat Standard Features:

- Carrot 3 Backrest
- Recline Adjustment (90° - 180°)
- Head Support: 8.5" x 10.5" (H x L)
- 4 Point Positioning Harness
- H-Harness Attaching Bracket
- Booster Seat Base
- Lower Anchors Connector System
- 2" Backrest Extension
- Shoulder Protector Wings

### Standard Models:

	Price
<input type="checkbox"/> Carrot 3 Booster Seat CRT3BLACK - CCRT3BS-10	<b>\$2,868</b>
<input type="checkbox"/> Carrot 3 Booster Seat CRT3BLUE - CCRT3BS-20	<b>\$2,868</b>
<input type="checkbox"/> Carrot 3 Booster Seat CRT3RED - CCRT3BS-30	<b>\$2,868</b>

### Carrot Booster Seat Recommended Item

	Price
<input type="checkbox"/> Tether strap (Black)	<b>\$93</b>
<input type="checkbox"/> Top Securement Strap*	<b>\$55</b>
<input type="checkbox"/> Lower Securement Strap*	<b>\$68</b>

*\*For vehicles not equipped with anchor points.*

Black	Blue	Red
CCT3TS0-10	CCT3TS0-10	CCT3TS0-10
CCT3HSF-10	CCT3HSF-10	CCT3HSF-10
CCT3LFS-10	CCT3LFS-10	CCT3LFS-10

### Accessories

	Price
<input type="checkbox"/> Booster Support Tray (Crash Tested)	<b>\$302</b>
<input type="checkbox"/> 3D Body Balance Pads L/R Set	<b>\$157</b>
<input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°)	<b>\$118</b>
<input type="checkbox"/> Backrest Extension - 4" (Extends back height to 32.8")	<b>\$95</b>
<input type="checkbox"/> Large Crotch Belt	<b>\$248</b>

Black	Blue	Red
CCRT3BT-19	CCRT3BT-29	CCRT3BT-39
CCT3BBP-00	CCT3BBP-00	CCT3BBP-00
CCT3SA0-10	CCT3SA0-20	CCT3SA0-30
CCT3BE4-10	CCT3BE4-20	CCT3BE4-30
CCT3CBX-00	CCT3CBX-00	CCT3CBX-00

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or Email [convoidsales.us@etac.com](mailto:convoidsales.us@etac.com)**

Prices in USD & subject to change, ex works Torrance, California.


Effective: 6/1/2022

Convoid  
2830 California Street  
Torrance, CA 90503  
1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF2100-ENG-REV.06



**CLIENT MEASUREMENTS:**



**Carrot 3 Base Model Specifications:**

Seat Width:	16.5"
Seat Depth:	15" - 16.5" *
Backrest Height Adjustment:	28" - 33"
Head Support Height:	8.5"
Backrest Height:	29" - 34"
Backrest Extension Height (Standard):	2"
Backrest Extension Height (Option):	4"
Shoulder Wing Height (Standard):	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"
Product Weight:	21 lbs.
Recline Adjustment:	90° - 180°
Product Overall Width:	20" - 22"

**User Measurement Guide:**

Minimum User Weight:	79 lbs.
Maximum User Weight:	165 lbs.
Minimum User Height:	55"
Maximum User Height:	70"

\*Factory Setting

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166