

Convaid Rodeo Tilt

Canadian Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

Rodeo Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- One Piece Height Adjustable Push Handle
- One Piece Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire
- Quick Release Rolko Rear Tires
- Single Foot Activated Wheel Lock Bar
- Zippered Pockets in Seat Bottom Upholstery
- 2-Point Positioning Belt
- Seat Depth Growth of 5"
- Elevating Legrest
- HCPCS code E1161 is only valid on size 16 (RD16)
- Tool-less Upholstery Removal
- Made in the USA

Standard Models:

	Price
<input type="checkbox"/> Rodeo RD10 - C904340	\$3,323
<input type="checkbox"/> Rodeo RD12 - C902689	\$3,323
<input type="checkbox"/> Rodeo RD14 - C902691	\$3,323
<input type="checkbox"/> Rodeo RD16 - C903134	\$3,323

All prices in CAD

Upholstery Colors

	Price
<input type="checkbox"/> Panther Black	\$0
<input type="checkbox"/> Royal Blue	\$0
<input type="checkbox"/> Sassy Purple	\$0
<input type="checkbox"/> Princess Pink	\$0
<input type="checkbox"/> Forest Green	\$0
<input type="checkbox"/> Cherry Red	\$0
<input type="checkbox"/> Chocolate Brown	\$0
<input type="checkbox"/> Aqua Blue	\$0
<input type="checkbox"/> Apple Green	\$0

Frame Option

	Price
<input type="checkbox"/> Recline Option (Requires Manufacturer Installation)	\$612
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	\$303
<input type="checkbox"/> Attendant Hand Brakes (pair) ¹	\$314
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires ¹	\$128
<input type="checkbox"/> Curb Tipper	\$133
<input type="checkbox"/> Rear Anti-Tip Tubes (Pair)	\$154

¹Manufacturer Installation Required.

Transit Required Items

	Price
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	\$275
<input type="checkbox"/> Headrest Extension	\$179
<input type="checkbox"/> H-Harness with Padded Covers	\$198
<input type="checkbox"/> 3-Pt Positioning Belt	\$68
Foot Positioners - Select Size Below:	
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$215
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$215
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$215

Tire Options

	Price
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	\$272
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$409

Head Supports

	Price
<input type="checkbox"/> Padded Headwings	\$325
<input type="checkbox"/> Occi Headwings	\$310
<input type="checkbox"/> Headrest Extension	\$179

Trunk Positioning

	Price
<input type="checkbox"/> H-Harness with Padded Covers (3-Pt Positioning Belt Req)	\$198
Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)	
<input type="checkbox"/> X Small= 2.75"x 3.75"	\$157
<input type="checkbox"/> Small = 3"x 5"	\$157
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$157
<input type="checkbox"/> Large = 4.75"x 6"	\$157
Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)	
<input type="checkbox"/> X Small= 2.75"x 3.75"	\$157
<input type="checkbox"/> Small = 3"x 5"	\$157
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$157
<input type="checkbox"/> Large = 4.75"x 6"	\$157
<input type="checkbox"/> Right S/A Adj. Hardware for Lateral Trunk Support	\$316
<input type="checkbox"/> Left S/A Adj. Hardware for Lateral Trunk Support	\$316
<input type="checkbox"/> Full Torso Support Vest	\$348

RD10	RD12	RD14	RD16
C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855
C904189	C904189	C904189	C904189
C904753	C904753	C904753	C904753
C904754	C904754	C904754	C904754

RD10	RD12	RD14	RD16
C904342	C904343	C904344	C904345
C904346	C904346	C904346	C904346
		C904370	C904371
C904650	C904650	C904650	C904650
C905310	C905310	C905310	C905310
C904649	C904649	C904649	C904649

RD10	RD12	RD14	RD16
C904341	C903167	C903168	C903169
C903332	C902863	C900484	C900138
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

RD10	RD12	RD14	RD16
C904117	C904117	C904117	C904117
C904118	C904118	C904118	C904118

RD10	RD12	RD14	RD16
C904365	C902673	C902653	C903147
C905270	C905270	C905271	C905271
C903332	C902863	C900484	C900138

RD10	RD12	RD14	RD16
C903988	C903988	C903988	C903988
C904385	C904385	C904385	C904385
	C904071	C904071	C904071
	C904072	C904072	C904072
	C904073	C904073	C904073
C904465	C904465	C904465	C904465
	C904466	C904466	C904466
	C904467	C904467	C904467
	C904468	C904468	C904468
C904362	C904362	C904362	C904362
C904363	C904363	C904363	C904363
C904364	C903261	C903260	C903259

- Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1") **\$482**
**RD10 requires X Small laterals.*

C904211	C904211	C904110	C904110
---------	---------	---------	---------

Pelvic Positioning

- 3-Pt Positioning Belt **\$68**
- Support - General Use Cushion*1 **\$330**
- Position Cushion - Medial Thigh Support/Anti Thrust*1 **\$636**
- Align Cushion - Lateral Pelvis/Lateral Thigh Support*1 **\$636**
- Incontinence Cushion Liner (Cushion req)¹ **\$99**

RD10	RD12	RD14	RD16
C904017	C904017	C904017	C904017
C904347	C904350	C904353	C904356
C904348	C904351	C904354	C904357
C904349	C904352	C904355	C904358
C902324	C902324	C902324	C902324

**Manufacturer Installation Required.*

**Solid Seat Insert included with all Cushions.*

Foot & Leg Positioning

- Medial Thigh Support (Abductor) **\$222**
- Lateral Thigh Support (Adductor) **\$222**
- Padded Elevating Legrest Cover **\$21**
- Legrest Extension - 2"*1 **\$75**
- Calf Panel **\$84**
- Individual Angle Adjustable Footplates¹ **\$244**
- Full Padded Footbox (One-Piece Footplate Required) **\$444**
- Foot Positioners - Small 7" - 10" (pair) **\$215**
- Foot Positioners - Medium 8" - 11" (pair) **\$215**
- Foot Positioners - Large 9" - 12" (pair) **\$215**

RD10	RD12	RD14	RD16
C904359	C902667	C902647	C903140
C904360	C902668	C902648	C903141
C904619	C904619	C904619	C904619
C904095	C904095	C904096	C904097
C904361	C904361	C904361	C904361
	C904957	C904958	C904959
C904953	C904954	C904955	C904956
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

**Manufacturer installation required if model was made prior to June 2019.*

**2" Legrest Extension interferes with caster swivel when seat is in upright position.*

Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

Transit Options

- Q'Straint Transit Lap Belt **\$123**
- Sure-Lok Transit Lap Belt **\$123**

RD10	RD12	RD14	RD16
C903094	C903094	C903094	C903094
C903232	C903232	C903232	C903232

Accessories

- Headrest Cover (Canopy) **\$303**
- Extended Headrest Cover - No Windows (Canopy) **\$342**
- Extended Headrest Cover - with Windows (Canopy) **\$404**
- Upper Extremity Support Surface (Clear Tray)* **\$152**
- Adj. Removable Hardware (Required for Clear Tray) **\$309**
- Upper Extremity Support Surface (Foam Tray) Small[^] **\$209**
- Upper Extremity Support Surface (Foam Tray) Medium[^] **\$209**
- Solid Side Panels **\$280**
- Mesh Side Panels **\$196**
- Medical Necessity Storage Basket **\$405**
- Under Seat Storage Basket - Tilt-In-Space **\$103**
- IV Pole (height adjustable and collapsible)² **\$205**
- Oxygen Tank Bag **\$238**
- Utility Bag **\$51**
- Rain Cover Size 1 (Canopy Req) **\$205**
- Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req) **\$205**
- Rain Cover Size 3 (Extended Canopy Req) **\$205**
- Mosquito Net Size 1 (Canopy req) **\$205**
- Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req) **\$205**
- Mosquito Net Size 3 (Extended Canopy req) **\$205**

RD10	RD12	RD14	RD16
C904366	C902679	C902659	C903152
C904728	C904729	C904730	C904731
C904744	C904745	C904746	C904747
C904367	C903405	C903401	C903402
C903369	C903369	C903369	C903369
C903669	C903669		
		C903670	C903670
C904388	C902672	C902652	C903145
C904389	C902680	C902681	C903146
C904476	C904476	C904477	C904477
C904529	C904529	C904530	C904530
C904566	C904566	C904566	C904566
C902870	C902870	C902870	C902870
C904531	C904531	C904531	C904531
C904939	C904939	C904939	
		C904940	C904940
			C905283
C904655	C904655	C904655	
		C904656	C904656
			C905319

**Flip-Up Armrests and Hardware required.*

²Anti-tip Tubes Required

[^]Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo mobility base was tested with Rodeo seating module equipped with head support, H-Harness with padded Covers, 3-point restraint system, foot positioners, and 4 transit anchor brackets.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specially designed for the wheelchair to be secured in a forward facing configuration and not in any other direction.
For more information on RESNA WC-19 and WC-20 standards please visit <http://www.recrwts.pitt.edu/WC19.html>.

Standard Shipping Rates: Ground \$69 • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 6/1/2022

OF0102-ENG-REV.07

Convaid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

