

# Convaid Rodeo Tilt - HCPCS E1232/E1161

## Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_

Contact: \_\_\_\_\_ ATP/RTS: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

### Rodeo Tilt Standard Features:

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>• 5° to 45° of Adjustable Tilt</li> <li>• Silver Powder Coated Frame</li> <li>• One Piece Height Adjustable Push Handle</li> <li>• One Piece Angle Adjustable Footplate</li> <li>• 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire</li> <li>• Quick Release Rolko Rear Tires</li> <li>• Single Foot Activated Wheel Lock Bar</li> </ul> | <ul style="list-style-type: none"> <li>• Zippered Pockets in Seat Bottom Upholstery</li> <li>• 2-Point Positioning Belt</li> <li>• Seat Depth Growth of 5"</li> <li>• Elevating Legrest</li> <li>• HCPCS code E1161 is only valid on size 16 (RD16)</li> <li>• Tool-less Upholstery Removal</li> <li>• Made in the USA</li> </ul> |
|--|---|

#### Standard Models:

	Price
<input type="checkbox"/> Rodeo RD10 - C904340	\$2,574
<input type="checkbox"/> Rodeo RD12 - C902689	\$2,574
<input type="checkbox"/> Rodeo RD14 - C902691	\$2,574
<input type="checkbox"/> Rodeo RD16 - C903134	\$2,574

<u>Upholstery Colors</u>	<u>Price</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Panther Black	\$0	C903854	C903854	C903854	C903854
<input type="checkbox"/> Royal Blue	\$0	C903850	C903850	C903850	C903850
<input type="checkbox"/> Sassy Purple	\$0	C903851	C903851	C903851	C903851
<input type="checkbox"/> Princess Pink	\$0	C903759	C903759	C903759	C903759
<input type="checkbox"/> Forest Green	\$0	C903857	C903857	C903857	C903857
<input type="checkbox"/> Cherry Red	\$0	C903855	C903855	C903855	C903855
<input type="checkbox"/> Chocolate Brown	\$0	C904189	C904189	C904189	C904189
<input type="checkbox"/> Aqua Blue	\$0	C904753	C904753	C904753	C904753
<input type="checkbox"/> Apple Green	\$0	C904754	C904754	C904754	C904754

<u>Frame Option</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Recline Option (Requires Manufacturer Installation)	\$481	E1225	C904342	C904343	C904344	C904345
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	\$238	E0973	C904346	C904346	C904346	C904346
<input type="checkbox"/> Attendant Hand Brakes (pair) <sup>1</sup>	\$246				C904370	C904371
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires <sup>1</sup>	\$98	K0073	C904650	C904650	C904650	C904650
<input type="checkbox"/> Curb Tipper	\$102		C905310	C905310	C905310	C905310
<input type="checkbox"/> Rear Anti-Tip Tubes (Pair)	\$120	E0971	C904649	C904649	C904649	C904649

<sup>1</sup>Manufacturer Installation Required.

<u>Transit Required Items</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	\$209		C904341	C903167	C903168	C903169
<input type="checkbox"/> Headrest Extension	\$135	E0966	C903332	C902863	C900484	C900138
<input type="checkbox"/> H-Harness with Padded Covers	\$150	E0960	C903988	C903988	C903988	C903988
<input type="checkbox"/> 3-Pt Positioning Belt	\$51	E0978	C904017	C904017	C904017	C904017
Foot Positioners - Select Size Below:						
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$163	K0108	C903415	C903415	C903415	C903415
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$163	K0108	C903416	C903416	C903416	C903416
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$163	K0108	C903417	C903417	C903417	C903417

<u>Tire Options</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	\$208		C904117	C904117	C904117	C904117
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$311		C904118	C904118	C904118	C904118

<u>Head Supports</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> Padded Headwings	\$254	E0955	C904365	C902673	C902653	C903147
<input type="checkbox"/> Occi Headwings	\$238	E0955	C905270	C905270	C905271	C905271
<input type="checkbox"/> Headrest Extension	\$135	E0966	C903332	C902863	C900484	C900138

<u>Trunk Positioning</u>	<u>Price</u>	<u>HCPCS</u>	<u>RD10</u>	<u>RD12</u>	<u>RD14</u>	<u>RD16</u>
<input type="checkbox"/> H-Harness with Padded Covers (3-Pt Positioning Belt Req)	\$150	E0960	C903988	C903988	C903988	C903988
Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)						
<input type="checkbox"/> X Small= 2.75"x 3.75"	\$123	E0956	C904385	C904385	C904385	C904385
<input type="checkbox"/> Small = 3"x 5"	\$123	E0956		C904071	C904071	C904071
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$123	E0956		C904072	C904072	C904072
<input type="checkbox"/> Large = 4.75"x 6"	\$123	E0956		C904073	C904073	C904073
Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)						
<input type="checkbox"/> X Small= 2.75"x 3.75"	\$123	E0956	C904465	C904465	C904465	C904465
<input type="checkbox"/> Small = 3"x 5"	\$123	E0956		C904466	C904466	C904466
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$123	E0956		C904467	C904467	C904467
<input type="checkbox"/> Large = 4.75"x 6"	\$123	E0956		C904468	C904468	C904468
<input type="checkbox"/> Right S/A Adj. Hardware for Lateral Trunk Support	\$248	E1028	C904362	C904362	C904362	C904362
<input type="checkbox"/> Left S/A Adj. Hardware for Lateral Trunk Support	\$248	E1028	C904363	C904363	C904363	C904363
<input type="checkbox"/> Full Torso Support Vest	\$248		C904364	C903261	C903260	C903259
<input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")	\$368		C904211	C904211	C904110	C904110

\*RD10 requires X Small laterals.

<u>Pelvic Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 3-Pt Positioning Belt	\$51	E0978
<input type="checkbox"/> Support - General Use Cushion* <sup>1</sup>	\$251	E2601
<input type="checkbox"/> Position Cushion - Medial Thigh Support/Anti Thrust* <sup>1</sup>	\$486	E2607
<input type="checkbox"/> Align Cushion - Lateral Pelvis/Lateral Thigh Support* <sup>1</sup>	\$486	E2607
<input type="checkbox"/> Incontinence Cushion Liner (Cushion req) <sup>1</sup>	\$75	E2619

RD10	RD12	RD14	RD16
C904017	C904017	C904017	C904017
C904347	C904350	C904353	C904356
C904348	C904351	C904354	C904357
C904349	C904352	C904355	C904358
C902324	C902324	C902324	C902324

<sup>1</sup>Manufacturer Installation Required.

\*Solid Seat Insert included with all Cushions.

<u>Foot &amp; Leg Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Medial Thigh Support (Abductor)	\$169	E0957
<input type="checkbox"/> Lateral Thigh Support (Adductor)	\$169	E0956
<input type="checkbox"/> Padded Elevating Legrest Cover	\$15	
<input type="checkbox"/> Legrest Extension - 2"* <sup>1</sup>	\$57	
<input type="checkbox"/> Calf Panel	\$64	E0995
<input type="checkbox"/> Individual Angle Adjustable Footplates <sup>1</sup>	\$186	K0040
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	\$339	K0108
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$163	K0108
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$163	K0108
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$163	K0108

RD10	RD12	RD14	RD16
C904359	C902667	C902647	C903140
C904360	C902668	C902648	C903141
C904619	C904619	C904619	C904619
C904095	C904095	C904096	C904097
C904361	C904361	C904361	C904361
	C904957	C904958	C904959
C904953	C904954	C904955	C904956
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

<sup>1</sup>Manufacturer installation required if model was made prior to June 2019.

\*2" Legrest Extension interferes with caster swivel when seat is in upright position.

Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

<u>Transit Options</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$94	E0978
<input type="checkbox"/> Sure-Lok Transit Lap Belt	\$94	E0978

RD10	RD12	RD14	RD16
C903094	C903094	C903094	C903094
C903232	C903232	C903232	C903232

<u>Accessories</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$238	
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$260	
<input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy)	\$307	
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray)*	\$116	E0950
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$237	E1028
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small <sup>^</sup>	\$158	E0950
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium <sup>^</sup>	\$158	E0950
<input type="checkbox"/> Solid Side Panels	\$220	
<input type="checkbox"/> Mesh Side Panels	\$155	
<input type="checkbox"/> Medical Necessity Storage Basket	\$307	
<input type="checkbox"/> Under Seat Storage Basket - Tilt-In-Space	\$78	
<input type="checkbox"/> IV Pole (height adjustable and collapsible) <sup>2</sup>	\$157	K0105
<input type="checkbox"/> Oxygen Tank Bag	\$186	E2208
<input type="checkbox"/> Utility Bag	\$40	
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$157	
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req)	\$157	
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$157	
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$157	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req)	\$157	
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$157	

RD10	RD12	RD14	RD16
C904366	C902679	C902659	C903152
C904728	C904729	C904730	C904731
C904744	C904745	C904746	C904747
C904367	C903405	C903401	C903402
C903369	C903369	C903369	C903369
C903669	C903669		
		C903670	C903670
C904388	C902672	C902652	C903145
C904389	C902680	C902681	C903146
C904476	C904476	C904477	C904477
C904529	C904529	C904530	C904530
C904566	C904566	C904566	C904566
C902870	C902870	C902870	C902870
C904531	C904531	C904531	C904531
C904939	C904939	C904939	
		C904940	C904940
			C905283
C904655	C904655	C904655	
		C904656	C904656
			C905319

\*Flip-Up Armrests and Hardware required.

<sup>2</sup>Anti-tip Tubes Required

<sup>^</sup>Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo mobility base was tested with Rodeo seating module equipped with head support, H-Harness with padded Covers, 3-point restraint system, foot positioners, and 4 transit anchor brackets.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specially designed for the wheelchair to be secured in a forward facing configuration and not in any other direction.

For more information on RESNA WC-19 and WC-20 standards please visit <http://www.ercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email [convaidsales.us@etac.com](mailto:convaidsales.us@etac.com)**

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166

Prices in USD & subject to change, ex works Torrance, California.

Effective: 6/1/2022

Convaid  
2830 California Street  
Torrance, CA 90503  
1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF0100-ENG-REV.07

