

# Convoid Rodeo Tilt w/out Seating Transit - HCPCS E1234

## Retail Price List / Order Form

Account #: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

Date: \_\_\_\_\_

Contact: \_\_\_\_\_ ATP/RTS: \_\_\_\_\_

Mark For: \_\_\_\_\_

|                 |                 |
|-----------------|-----------------|
| Bill To:        | Ship To:        |
| Address:        | Address:        |
| City/State/Zip: |                 |
| Phone:          | City/State/Zip: |
| Email:          | Phone:          |

**Rodeo Tilt w/out Seating Standard Features:**

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- Adjustable One-Piece Push Handle
- Elevating Legrest / One-Piece, Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires

- Quick Release Rolko Rear Wheels
- Single Foot Activated Wheel Lock Bar
- Seat Depth Growth of 5"
- Made in the USA

**Transit Standard Features:**

- Transit Anchor Brackets
- 2-Point Positioning Belt

**Standard Models:**

|   | <u>Price</u>   |
|---|----------------|
| <input type="checkbox"/> Rodeo Transit w/out Seating RD10TWOS - C904967 | <b>\$2,715</b> |
| <input type="checkbox"/> Rodeo Transit w/out Seating RD12TWOS - C904968 | <b>\$2,715</b> |
| <input type="checkbox"/> Rodeo Transit w/out Seating RD14TWOS - C904969 | <b>\$2,715</b> |
| <input type="checkbox"/> Rodeo Transit w/out Seating RD16TWOS - C904970 | <b>\$2,715</b> |

| <u>Seating Module</u>   | <u>Price</u> | <u>HCPCS</u> | <u>RD10</u> | <u>RD12</u> | <u>RD14</u> | <u>RD16</u> |
|---|--------------|--------------|-------------|-------------|-------------|-------------|
| <input type="checkbox"/> Depth Adjustable Solid Seat Pan (Manufacturer Installation Required) | <b>\$157</b> | <b>E2293</b> | C509549     | C509550     | C509551     | C509552     |

| <u>Frame Option</u>  | <u>Price</u> | <u>HCPCS</u> | <u>RD10</u> | <u>RD12</u> | <u>RD14</u> | <u>RD16</u> |
|--|--------------|--------------|-------------|-------------|-------------|-------------|
| <input type="checkbox"/> Recline Option (Manufacturer Installation Required) | <b>\$481</b> | <b>E1225</b> | C904342     | C904343     | C904344     | C904345     |
| <input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)           | <b>\$238</b> | <b>E0973</b> | C904346     | C904346     | C904346     | C904346     |
| <input type="checkbox"/> Curb Tipper   | <b>\$102</b> |              | C905310     | C905310     | C905310     | C905310     |
| <input type="checkbox"/> Rear Anti-Tip Tubes (Pair)*                         | <b>\$120</b> | <b>E0971</b> | C904649     | C904649     | C904649     | C904649     |

| <u>Tire Options</u>  | <u>Price</u> | <u>HCPCS</u> | <u>RD10</u> | <u>RD12</u> | <u>RD14</u> | <u>RD16</u> |
|--|--------------|--------------|-------------|-------------|-------------|-------------|
| <input type="checkbox"/> 12.5" Rear Solid Knobby Tire                      | <b>\$208</b> |              | C904117     | C904117     | C904117     | C904117     |
| <input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire | <b>\$311</b> |              | C904118     | C904118     | C904118     | C904118     |

| <u>Foot &amp; Leg Positioning</u>   | <u>Price</u> | <u>HCPCS</u> | <u>RD10</u> | <u>RD12</u> | <u>RD14</u> | <u>RD16</u> |
|---|--------------|--------------|-------------|-------------|-------------|-------------|
| <input type="checkbox"/> Individual Angle Adjustable Footplates             | <b>\$186</b> | <b>K0040</b> |             | C904957     | C904958     | C904959     |
| <input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required) | <b>\$339</b> | <b>K0108</b> | C904953     | C904954     | C904955     | C904956     |
| <input type="checkbox"/> Legrest Extension - 2"                             | <b>\$57</b>  |              | C904095     | C904095     | C904096     | C904097     |
| <input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)           | <b>\$163</b> | <b>K0108</b> | C903415     | C903415     | C903415     | C903415     |
| <input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)          | <b>\$163</b> | <b>K0108</b> | C903416     | C903416     | C903416     | C903416     |
| <input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)           | <b>\$163</b> | <b>K0108</b> | C903417     | C903417     | C903417     | C903417     |

**\*2" Legrest Extension interferes with caster swivel when seat is in upright position.  
Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.**

| <u>Canopies</u>  | <u>Price</u> | <u>HCPCS</u> | <u>RD10</u> | <u>RD12</u> | <u>RD14</u> | <u>RD16</u> |
|--|--------------|--------------|-------------|-------------|-------------|-------------|
| <input type="checkbox"/> Headrest Cover (Canopy)                         | <b>\$238</b> |              | C904366     | C902679     | C902659     | C903152     |
| <input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)   | <b>\$260</b> |              | C904728     | C904729     | C904730     | C904731     |
| <input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy) | <b>\$307</b> |              | C904744     | C904745     | C904746     | C904747     |
| Canopy Upholstery (Please Choose One Below)                              |              |              |             |             |             |             |
| <input type="checkbox"/> Panther Black                                   | <b>\$0</b>   |              | C903854     | C903854     | C903854     | C903854     |
| <input type="checkbox"/> Royal Blue                                      | <b>\$0</b>   |              | C903850     | C903850     | C903850     | C903850     |
| <input type="checkbox"/> Sassy Purple                                    | <b>\$0</b>   |              | C903851     | C903851     | C903851     | C903851     |
| <input type="checkbox"/> Princess Pink                                   | <b>\$0</b>   |              | C903759     | C903759     | C903759     | C903759     |
| <input type="checkbox"/> Forest Green                                    | <b>\$0</b>   |              | C903857     | C903857     | C903857     | C903857     |
| <input type="checkbox"/> Cherry Red                                      | <b>\$0</b>   |              | C903855     | C903855     | C903855     | C903855     |
| <input type="checkbox"/> Chocolate Brown                                 | <b>\$0</b>   |              | C904189     | C904189     | C904189     | C904189     |
| <input type="checkbox"/> Aqua Blue                                       | <b>\$0</b>   |              | C904753     | C904753     | C904753     | C904753     |
| <input type="checkbox"/> Apple Green                                     | <b>\$0</b>   |              | C904754     | C904754     | C904754     | C904754     |

| <u>Transit Options</u>  | <u>Price</u> | <u>HCPCS</u> | <u>RD10</u> | <u>RD12</u> | <u>RD14</u> | <u>RD16</u> |
|---|--------------|--------------|-------------|-------------|-------------|-------------|
| <input type="checkbox"/> Q'Straint Transit Lap Belt                           | <b>\$94</b>  | <b>E0978</b> | C903094     | C903094     | C903094     | C903094     |
| <input type="checkbox"/> Sure-Lok Transit Lap Belt                            | <b>\$94</b>  | <b>E0978</b> | C903232     | C903232     | C903232     | C903232     |
| <input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small ^  | <b>\$158</b> | <b>E0950</b> | C903669     | C903669     |             |             |
| <input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium ^ | <b>\$158</b> | <b>E0950</b> |             |             | C903670     | C903670     |

**^ Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.**

| Accessories  | Price | HCPCS |
|--|-------|-------|
| <input type="checkbox"/> Upper Extremity Support Surface (Clear Tray)^               | \$116 | E0950 |
| <input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)           | \$237 | E1028 |
| <input type="checkbox"/> Solid Side Panels   | \$220 | K0108 |
| <input type="checkbox"/> Mesh Side Panels  | \$155 | K0108 |
| <input type="checkbox"/> Medical Necessity Storage Basket                            | \$307 |       |
| <input type="checkbox"/> Under Seat Storage Basket - Tilt-in-Space                   | \$78  |       |
| <input type="checkbox"/> IV Pole (height adjustable and collapsible) <sup>2</sup>    | \$157 | K0105 |
| <input type="checkbox"/> Oxygen Tank Bag   | \$186 | E2208 |
| <input type="checkbox"/> Utility Bag   | \$40  |       |
| <input type="checkbox"/> Attendant Hand Brakes (pair) <sup>1</sup>                   | \$246 |       |
| <input type="checkbox"/> Caster Lock for 2" Wide Front Tires <sup>1</sup>            | \$98  | K0073 |
| <input type="checkbox"/> Rain Cover Size 1 (Canopy Req)                              | \$157 |       |
| <input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req)   | \$157 |       |
| <input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)                     | \$157 |       |
| <input type="checkbox"/> Mosquito Net Size 1 (Canopy req)                            | \$157 |       |
| <input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req) | \$157 |       |
| <input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)                   | \$157 |       |

<sup>1</sup>Flip-Up Armrests and Hardware required.  
<sup>1</sup>Manufacturer Installation Required.  
<sup>2</sup>Anti-tip Tubes Required

| RD10    | RD12    | RD14    | RD16    |
|---------|---------|---------|---------|
| C904367 | C903405 | C903401 | C903402 |
| C903369 | C903369 | C903369 | C903369 |
| C904388 | C902672 | C902652 | C903145 |
| C904389 | C902680 | C902681 | C903146 |
| C904476 | C904476 | C904477 | C904477 |
| C904529 | C904529 | C904530 | C904530 |
| C904566 | C904566 | C904566 | C904566 |
| C902870 | C902870 | C902870 | C902870 |
| C904531 | C904531 | C904531 | C904531 |
|         |         | C904370 | C904371 |
| C904650 | C904650 | C904650 | C904650 |
| C904939 | C904939 | C904939 |         |
|         |         | C904940 | C904940 |
|         |         |         | C905283 |
| C904655 | C904655 | C904655 |         |
|         |         | C904656 | C904656 |
|         |         |         | C905319 |

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo Mobility Base was tested with the Rodeo Seating Module, equipped with a head support and wheelchair-anchored 2-point pelvic lap belt. Do not use the foot positioners on the Convaid Rodeo during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is or can being made.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specifically designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.ercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email [convaidsales.us@etac.com](mailto:convaidsales.us@etac.com)**

Prices in USD & subject to change, ex works Torrance, California.

Effective: 6/1/2022

OF0110-ENG-REV.06

Convaid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

