

Convoid Safari Tilt - HCPCS E1232

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone: Fax:	City/State/Zip:
Email:	Phone:

Safari Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (SFT12/14/16)
- 6" Front, 11" Rear Solid Tires (SFT18)
- Quick Release Rear Wheels (SFT12/14/16)
- Wheel Locks
- 2-Point Positioning Belt
- Heavy-Duty Back Upholstery
- One Piece, Non-Folding/Non-Elevating Footplate
- Support Cushion with Solid Seat Insert
- Made in the USA

Standard Models:

	Price
<input type="checkbox"/> Safari SFT12 - C900493	\$2,466
<input type="checkbox"/> Safari SFT14 - C901104	\$2,466
<input type="checkbox"/> Safari SFT16 - C900526	\$2,466
<input type="checkbox"/> Safari SFT18 - C902432	\$2,466

<u>Upholstery Colors</u>	Price
<input type="checkbox"/> Panther Black	\$0
<input type="checkbox"/> Royal Blue	\$0
<input type="checkbox"/> Sassy Purple	\$0
<input type="checkbox"/> Princess Pink	\$0
<input type="checkbox"/> Forest Green	\$0
<input type="checkbox"/> Cherry Red	\$0
<input type="checkbox"/> Chocolate Brown	\$0
<input type="checkbox"/> Aqua Blue	\$0
<input type="checkbox"/> Apple Green	\$0

<u>Frame Option</u>	Price	HCPCS
<input type="checkbox"/> Anatomic Back	\$349	
<input type="checkbox"/> Recline Option	\$481	E1225
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	\$238	E0973
<input type="checkbox"/> Attendant Hand Brakes (pair)*	\$246	
<input type="checkbox"/> Solid Side Panels (N/A w/ Recline Option)	\$220	
<input type="checkbox"/> Curb Tipper (elevates front of chair to go up curbs)	\$58	

**Available with Knobby Tires only. Not available with Recline Option.*

<u>Transit Required Items</u>	Price	HCPCS
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	\$209	
<input type="checkbox"/> Headrest Extension	\$135	E0966
<input type="checkbox"/> H-Harness with Padded Covers	\$150	E0960
<input type="checkbox"/> 3-Pt Positioning Belt	\$51	E0978
<input type="checkbox"/> Elevating Legrests with Swing-Away Footplates (pair)	\$282	E0990
Foot Positioners - Select Size Below:		
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$163	
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$163	
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$163	

<u>Tire Options</u>	Price	HCPCS
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	\$208	
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$311	
<input type="checkbox"/> 12.5" x 2" Rear Solid Knobby Tire	\$209	

<u>Head Supports</u>	Price	HCPCS
<input type="checkbox"/> Padded Headwings	\$254	E0955
<input type="checkbox"/> Occi Headwings	\$238	E0955
<input type="checkbox"/> Headrest Extension	\$135	E0966

<u>Trunk Positioning</u>	Price	HCPCS
<input type="checkbox"/> H-Harness with Padded Covers (3-Pt Positioning Belt Req)	\$150	E0960
<input type="checkbox"/> SA Adj Lateral Support (single flap) with Scoli Strap	\$142	E0956
<input type="checkbox"/> SA Adj Lateral Support (double flap) with Scoli Strap	\$142	E0956
<input type="checkbox"/> Full Torso Support Vest*	\$275	
<input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")	\$368	

**Not available with Single or Double Flap Lateral Support.*

<u>Pelvic Positioning</u>	Price	HCPCS
<input type="checkbox"/> 3-Pt Positioning Belt	\$51	E0978
<input type="checkbox"/> Position Cushion - Medial Thigh Support/Anti Thrust*	\$362	E2607
<input type="checkbox"/> Align Cushion - Lateral Pelvis/Lateral Thigh Support*	\$362	E2607
<input type="checkbox"/> Incontinence Cushion Liner (Cushion and Solid Seat req)	\$75	E2619
<input type="checkbox"/> Solid Seat Insert (Required with Cushion)	\$124	E0992

**Solid Seat Insert required with any Cushion.*

SFT12	SFT14	SFT16	SFT18
C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855
C904189	C904189	C904189	C904189
C904753	C904753	C904753	C904753
C904754	C904754	C904754	C904754

SFT12	SFT14	SFT16	SFT18
C900697	C901164	C901319	C902476
C900689	C901139	C901321	
C900716	C901143	C901322	C902477
C904064	C904064	C904064	C904064
C904390	C904391	C904392	C904393
C904063	C904063	C904063	C904063

SFT12	SFT14	SFT16	SFT18
C900494	C901141	C900527	C902435
C900712	C901168	C901310	C902471
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C900698	C901117	C901299	C902460
/			
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

SFT12	SFT14	SFT16	SFT18
C904117	C904117	C904117	
C904118	C904118	C904118	
			C903420

SFT12	SFT14	SFT16	SFT18
C900704	C901111	C901309	C902470
C905270	C905271	C905271	C905273
C900712	C901168	C901310	C902471

SFT12	SFT14	SFT16	SFT18
C903988	C903988	C903988	C903988
C900705	C901112	C901306	C902467
C900706	C901144	C901307	C902468
C900707	C901217	C901308	C902469
C904211	C904110	C904110	C904110

SFT12	SFT14	SFT16	SFT18
C904017	C904017	C904017	C904017
C903588	C903589	C903590	C903591
C903592	C903593	C903594	C903595
C902324	C902324	C902324	C902324
C904029	C904029	C904029	C904029

<u>Foot & Leg Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Medial Thigh Support (Abductor)	\$169	E0957
<input type="checkbox"/> Lateral Thigh Support (Adductor)	\$169	E0956
<input type="checkbox"/> Elevating Legrests w/ Swing-Away Footplates (pair)	\$282	E0990
<input type="checkbox"/> Legrest Extension - 3"*^	\$57	
<input type="checkbox"/> Calf Panel	\$64	E0995
<input type="checkbox"/> Angle Adj. Footplates - Std Size (pair)^	\$186	K0040
<input type="checkbox"/> Angle Adj. Footplates - Adult Size (pair)^	\$186	K0040
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)^	\$163	
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)^	\$163	
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)^	\$163	
<input type="checkbox"/> Footplate Securement Strap	\$32	

***3" Legrest Extension interferes with caster swivel when seat is in upright position.
 Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.
 ^Elevating Legrests with Swing-Away Footplates required.**

SFT12	SFT14	SFT16	SFT18
C900708	C901113	C901297	C902458
C900709	C901145	C901298	C902459
C900698	C901117	C901299	C902460
	C901325	C901324	C902478
C904580	C904580	C904580	C904580
C903270	C903270		
		C903271	C903271
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417
C903584	C903584	C903584	C903584

<u>Transit Options</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$94	E0978
<input type="checkbox"/> Sure-Lok Transit Lap Belt	\$94	E0978

SFT12	SFT14	SFT16	SFT18
C903095	C903095	C903095	C903095
C903233	C903233	C903233	C903233

<u>Accessories</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$238	
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$260	
<input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy)	\$307	
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray) ^	\$116	E0950
<input type="checkbox"/> Adj. Removeable Hardware (Required for Clear Tray)	\$237	E1028
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small^	\$158	E0950
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium^	\$158	E0950
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Large^	\$158	E0950
<input type="checkbox"/> Wheelchair Travel Bag	\$215	
<input type="checkbox"/> Utility Bag	\$40	
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$157	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req)	\$157	

SFT12	SFT14	SFT16	SFT18
C900714	C901116	C901313	C902475
C904732	C904733	C904734	C904735
C904748	C904749	C904750	C904751
C903401	C903402	C903403	C903404
C904519	C904519	C904519	C904519
C903669			
	C903670	C903670	
			C903671
C902022	C902022	C902022	C902022
C904531	C904531	C904531	C904531
C904655	C904655		
		C904656	C904656

**^Flip-Up Armrests required.
 ^Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.**

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 6/1/2022

Convoid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0700-ENG-REV.04

