


Page 1: Standard Assembly


SIZE	KIT NAME	ITEM NO.	DESCRIPTION	QTY2	KIT	Price
18	8" Urethane Wheel & Bearing Assy 7/16"	27	8" Urethane Wheel & Bearing Assy 7/16"	1	992647	Call Customer Service
12	Footplate Assembly SMALL - LH	29	CX-10 Footrest Assembly, LH	1	992658	\$ 156.80
12	Footplate Assembly SMALL - RH	30	CX-10 Footrest Assembly, RH	1	992659	\$ 156.80
14	Footplate Assembly STD - LH	29	Footrest Assembly - Standard LHS	1	992660	\$ 202.80
14	Footplate Assembly STD - RH	30	Footrest Assembly - Standard RHS	1	992661	\$ 202.80
16, 18	Footplate Assembly Adult - LH	29	Footrest Assy, Adult, LH	1	992662	\$ 225.00
16, 18	Footplate Assembly Adult - RH	30	Footrest Assy, Adult, RH	1	992663	\$ 225.00
12, 14, 16	11.5" x 2.5" Rear Solid Tire on Grey Rim (Ea)	27	Rear Wheel 290mm x 60mm (Rolko)	1	994746	\$ 129.00




Page 2: Rider Bus Kit

SIZE	KIT NAME	ITEM NO.	DESCRIPTION	QTY2	KIT	Price		
12	EZ12T Transit Conversion Kit	6	Locknut 1/4"-28 Jam Flexlock Black Zinc	2	990461	\$ 447.60		
12		7	Washer 1/4" AN-960-416	10				
12		8	Transport Anchor	2				
12		9	Bolt 1/4"-28 x 3" Socket Button Head Alloy S	2				
12		13	Caution Card WC/19	1				
12		14	Anchor Logo	4				
12		15	Transit Anchor Lock Pin 1-3/8"	2				
12		16	Transit Anchor Lock Pin Sleeve	2				
12		17	Inline Transit Anchor	2				
12		18	BOLT 1/4"-20 X 1-1/2" Socket Button Head	4				
14		6	Locknut 1/4"-28 Jam Flexlock Black Zinc	2				
14		7	Washer 1/4" AN-960-416	12				

14	EZ14T Transit Conversion Kit	8	Transport Anchor	2	990561	\$ 447.60
14		9	Rivet 1/4" X 1-5/8" Stainless Steel	2		
14		13	Caution Card WC/19	1		
14		14	Anchor Logo	4		
14		15	Transit Anchor Lock Pin 1-5/8"	2		
14		16	Transit Anchor Lock Pin Sleeve	2		
14		17	Inline Transit Anchor	2		
14		18	BOLT 1/4"-20 X 1-1/2" Socket Button Head	4		
14		19	Removable Thread Locker (Blue)-Medium S	AR		
16	EZ16T Transit Conversion Kit	6	Locknut 1/4"-28 Jam Flexlock Black Zinc	2	990661	\$ 447.60
16		7	Washer 1/4" AN-960-416	10		
16		8	Transport Anchor	2		
16		9	Rivet 1/4" X 1-5/8" Stainless Steel	2		
16		13	Caution Card WC/19	1		
16		14	Anchor Logo	2		
16		15	Transit Anchor Lock Pin 1-3/4"	2		
16		16	Transit Anchor Lock Pin Sleeve	2		
16		-	Rubber Dampener	2		
16	17	Inline Transit Anchor	2			
16	18	BOLT 1/4"-20 X 1-1/2" Socket Button Head	4			
18	EZ18T Transit Conversion Kit	6	Locknut 1/4"-28 Jam Flexlock Black Zinc	2	990761	\$ 447.60
18		7	Washer 1/4" AN-960-416	10		
18		8	Transport Anchor	2		
18		9	Rivet 1/4" X 1-5/8" Stainless Steel	2		
18		13	Caution Card WC/19	1		
18		14	Anchor Logo	4		
18		7	Transit Anchor Lock Pin 1-3/4"	2		
18		16	Transit Anchor Lock Pin Sleeve	2		
18		17	Inline Transit Anchor	2		
18	18	BOLT 1/4"-20 X 1-1/2" Socket Button Head	4			







Page 3: Rear Frame Assembly, Size  
12-16 Only. For Size 18, please  
scroll to Page 4

SIZE	KIT NAME	ITEM NO.	DESCRIPTION	QTY2	KIT	Price
12	EZ12 Front Assembly - Std	1	EZ-12 Front Frame Assembly	1	990401	\$ 913.20
12	EZ12 Rear Lock Brace Assembly	29	EZ-12 Rear Lock Assembly	1	990410	\$ 258.00
12	EZ12 Seat Tubes Assembly	35	EZ-12 Seat Tube	2	990450	\$ 80.40
12	EZ12 Seat Back Tubes - pair	42	EZ-12 Fabric Tube Assembly	2	990469	\$ 63.00
12	EZ12 Seat Support w/Velcro	38	EZ-12 Seat Support w/Velcro	1	990496	\$ 26.00
14	EZ14 Front Assembly - Std	1	EZ-14 Front frame assy	1	990501	\$ 913.20
14	EZ14 Rear Assembly - Std	2	EZ-14 Rear Frame Assembly	1	990503	\$ 744.00
14		12	EZ-14 Rear Upper Leg Welded	1		
14	EZ14 Rear Lock Brace Assembly	29	EZ-14 Rear Lock Assembly	1	990510	\$ 258.00
14	EZ14 Seat Tube Assembly	35	EZ-14 Seat Tube	2	990550	\$ 80.40
14	EZ14 Fabric Tube Outer Welded - LH	42	EZ-14 Back Tube Welded LH	1	990565	\$ 124.80
14	EZ14 Fabric Tube Outer Welded - RH	49	EZ-14 Back Tube Welded RH	1	990566	\$ 124.80
14	EZ14 One Piece Push Handle	53	EZ-14 Single Pc Handle Assembly	1	990592	\$ 110.00
14		54	Screw #10-24 X 1-1/4" Phil Pan	2		


14		55	Locknut #10-24 Jam Flexlock Black Zinc	2		
14	EZ14 Seat Support w/Velcro	38	EZ-14 Seat Support w/Velcro	1	990596	\$ 33.00
16	EZ16 Front Assembly - Std	1	EZ-16 Front Frame Assy	1	990601	#####
16	EZ16 Rear Assembly - Std	2	EZ-16 Rear Frame Assembly	1	990603	\$ 801.60
16	EZ16 Rear Lock Brace Assembly	29	EZ-16 Rear Lock Assembly	0	990610	\$ 258.00
16	EZ16 Upper Rear Leg Assy - Std LH	12	EZ-16 Rear Leg Upper Welded	1	990651	\$ 42.00
16		39	Locknut 1/4"-28 Flexlock Black Zinc	1		
16	EZ16 Fabric Tube Outer Welded - LH	42	EZ-16 Fabric Tube Outer Welded LH	1	990665	\$ 51.60
16	EZ16 Fabric Tube Outer Welded - RH	49	EZ-16 Fabric Tube Outer Welded RH	1	990666	\$ 51.60
16	EZ16 One Piece Push Handle	53	EZ-16 Single Pc Handle Assembly	1	990692	\$ 110.00
16		54	Screw #10-24 X 1-1/4" Phil Pan	2		
16		55	Locknut #10-24 Jam Flexlock Black Zinc	2		
16	EZ16 Seat Support w/Velcro	38	EZ-16 Seat Support w/Velcro	1	990696	\$ 34.00
12, 16	VV Seat Extension Tube - pair	4	VV-12 Seat Tube Extension Assy	1	994023	\$ 236.40
12, 14, 16	7.5" x 2" Front Solid Tire on Grey Rim (Ea)	28	Front Wheel w/3-Spokes 182.5x51mm Rolk	1	994749	\$ 130.00
12, 14, 16		30	Bolt M8-1.25mm x 90mm Socket Cap Scre	1		
12, 14, 16		31	Locknut M8-1.25 Nylon Insert	1		
12, 14, 16		50	Washer 5/16" X 9/16"	6		
12, 14, 16	Caster Fork Assembly 7.5" X 2" (ROLKO)	25	Balloon Caster Fork Assembly	1	996009	\$ 40.00
12, 14, 17		26	"0" Ring Buna-Nitrile	1		
14	EZ Closure Strap Replacement	11	Locknut 1/4"-28 Jam Flexlock Black Zinc	1	997315	\$ 13.00
14, 16		56	Footplate Strap	1		



Page 4: Rear Frame Assembly, Size  
18 Only

SIZE	KIT NAME	ITEM NO.	DESCRIPTION	QTY2	KIT	Price
18	EZ18 Front Assembly - Std	1	EZ-18 Front Frame Assembly	1	990701	#####
18		33	Handle Grip- 3/4"ID X 7"	2		
18		34	3M Rubber & Gasket Adhesive 5oz (847)	1		
18	EZ18 Rear Assembly - Std	2	EZ-18 Rear Frame Assy.	1	990703	\$ 801.60
18	EZ18 Rear Lock Brace Assembly	29	EZ-18 Rear Lock Brace Assy	1	990710	\$ 258.00
18	EZ18 Seat Tubes Assembly	37	EZ-18 Seat Tube	1	990750	\$ 88.80
18	EZ18 Fabric Tube Outer Welded - LH	47	EZ-18 Back Tube Welded LH	1	990765	\$ 51.60
18	EZ18 Fabric Tube Outer Welded - RH	48	EZ-18 Back Tube Welded RH	1	990766	\$ 51.60
18	EZ18 Seat Support w/Velcro	45	EZ-18 Seat Support w/Velcro	1	990796	\$ 37.00
18	Wheel Assembly with 3/8" Bearing - 8" Transit (Ea)	21	8" Urethane Wheel & Bearing Assy 3/8"	1	992646	\$ 77.00
18		31	Threaded Barrel	1		
18		32	Bolt 1/4"-28 x 7/8" Socket Button Head - blk z	2		
18		N/A	Thread Locker-Medium-Vibra-tite	1		
18	Footplate Assembly Adult - LH	54	Footrest Assy, Adult, LH	1	992662	\$ 225.00
18	Footplate Assembly Adult - RH	55	Footrest Assy, Adult, RH	1	992663	\$ 225.00
18	CM18 Upper Rear Leg Assembly LH	12	EZ-18 Rear Leg Upper Welded	1	994233	\$ 30.00
18	8" Rear Wheel Assy w/3/8" Bearing	21	8" Urethane Wheel & Bearing Assy 3/8"	1	994745	\$ 77.00
18	CX10, CX/CC/CL/CS12, EZ/CV18 Seat Ext Assembly	40	CX-10/12 Seat Extension Tube Assy Short	1	994777	\$ 15.00
18	Caster Fork Assembly 8" x 2"	26	Caster Stem Assembly	1	996309	\$ 40.00
18		27	"0" Ring Buna-Nitrile	1		
18		28	Locknut 1/2"-20 Capped Nylon - blk zinc 098	1		
18	EZ Closure Strap Replacement	11	Locknut 1/4"-28 Jam Flexlock Black Zinc	1	997315	\$ 13.00
18		57	Footplate Strap	1		






Page 5: Front Frame Assembly

SIZE	KIT NAME	ITEM NO.	DESCRIPTION	QTY2	KIT	Price
12	EZ12 Front X Brace Assembly	4	EZ-12 Front "X" Brace Assembly	1	990408	\$ 358.80
14	EZ14 Front X Brace Assembly	4	EZ-14 Front X brace assembly	1	990508	\$ 351.60
16	EZ16 Front X Brace Assembly	4	EZ-16 Front X-Brace Assembly	1	990608	\$ 394.80
18	EZ18 Front X Brace Assembly	4	EZ-18 Front X Brace Assembly	1	990708	\$ 394.80
18	EZ18 Push Handle Tube Assembly - Std	1	EZ-18 Handle Tube Assy, LH	1	990711	\$ 315.60
18		1	EZ-18 Handle Tube Assy, RH	1		
18		2	Caster Stem Bushing	2		
18		3	Flanged Bearing 1/2" shielded Greased	2		
All Sizes		11	Tube Connector 7/8" Heavy Duty	Varies per Kit		
18		-	Bolt 1/4"-28 x 1-1/2" Socket Button Head- b	2		
18		-	Locknut 1/4"-28 Capped Nylon 098	2		
12		EZ/CV12 Front Leg Tube Assy (1pc	1	EZ-12 Front Leg Assy for 1-Pc Handle, LH		
12	1		EZ-12 Front Leg Assy for 1-Pc Handle, RH	1		

12	Handle) - Std	2	Caster Stem Bushing	2		
12		3	Flanged Bearing 1/2" shielded Greased	2		
14	EZ/CV14 Front Leg Tube Assy (1pc Handle) - Std	1	EZ-14 Front Leg Assy for 1-Pc Handle, LH	1	992370	\$ 220.00
14		1	EZ-14 Front Leg Assy for 1-Pc Handle, RH	1		
14		2	Caster Stem Bushing	2		
14		3	Flanged Bearing 1/2" shielded Greased	2		
16	EZ/CV16 Front Leg Tube Assy (1pc Handle) - Std	1	EZ-16 Front Leg Welded Assy for 1-Pc Hand, LH	1	992470	\$ 220.00
16		1	EZ-16 Front Leg Welded Assy for 1-Pc Hand, RH	1		
16		2	Caster Stem Bushing	2		
16		3	Flanged Bearing 1/2" shielded Greased	2		
16		-	Pivot Bracket 7/8"-7/8"	2		
16		-	Rivet 3/16" X 1-1/8" S/T	4		
All Sizes	Detent Pin W/ Webbing Assembly	12	Detent Pin With Webbing Assembly	1	995605	\$ 2.65
All Sizes		11	Tube Connector 7/8" Heavy Duty	1		
All Sizes	Standard Footplate Assembly Kit LH	13	Footplate Adjustment Tube LH	1	996110	\$ 200.00
All Sizes	Standard Footplate Assembly Kit RH	14	Footplate Adjustment Tube RH	1	996111	\$ 200.00
12	Included in 992270 and 990408 Kits	-	Bolt 1/4"-28 x 1-1/2" Socket Button Head- b	2	Included in 992270 and	\$ -
12		-	Locknut 1/4"-28 Capped Nylon 098	2		
All Sizes	Included in 992270, 996110, and 996111 Kits	11	Tube Connector 7/8" Heavy Duty	Varies per Kit	Included in 992270, 996110, and 996111 Kits	\$ -
14	Included in 992370 and 990508 Kits	N/A	Bolt 1/4"-28 x 1-1/2" Socket Button Head- b	2	Included in 992370 and	\$ -
14		N/A	Locknut 1/4"-28 Capped Nylon 098	2		
All Sizes	Included in 992370, 996110, and 996111 Kits	11	Tube Connector 7/8" Heavy Duty	Varies per Kit	Included in 992370, 996110, and 996111 Kits	\$ -
16	Included in 992470 and 990608 Kits	-	Bolt 1/4"-28 x 1-1/2" Socket Button Head- b	2	Included in 992470 and	\$ -
16		-	Locknut 1/4"-28 Capped Nylon 098	2		
All Sizes	Included in 992470, 996110, and 996111 Kits	11	Tube Connector 7/8" Heavy Duty	Varies per Kit	Included in 992470, 996110, and 996111 Kits	\$ -
All Sizes	Included in 996110 and 996111	10	Footrest Extens. Housing Assy.	1	Included in	\$ -

All Sizes	Kits	12	Detent Pin With Webbing Assembly	1	996110 and	P
-----------	------	----	----------------------------------	---	------------	---



Page 6: Rear Frame Assembly

SIZE	KIT NAME	ITEM NO.	DESCRIPTION	QTY2	KIT	Price
12	EZ12 Rear X Brace Assembly	1	EZ-12 Rear "X" Brace Assembly	1	990409	\$ 272.40
14	EZ14 Rear X Brace Assembly	1	EZ-14 Rear X Brace Assembly	1	990509	\$ 272.40
16	EZ16 Rear X Brace Assembly	1	EZ-16 Rear X-Brace Assembly	1	990609	\$ 300.00
18	EZ18 Rear X Brace Assembly	1	EZ-18 Rear X Brace Assembly	1	990709	\$ 300.00
18		9	Washer 1/4" AN-960-416	2		
12, 14, 16		6	R82 Brake Assembly LH	1		
12, 14, 16		7	R82 Brake Assembly RH	1		
12, 14, 16		9	Washer 1/4" AN-960-416	4		
12, 14, 16		19	Locknut 1/4"-28 Flexlock Black Zinc	4		
12, 14, 16		23	Coved Washer - 3/4"	4		

12, 14, 16	R82 Brake Replacement Kit (CX10-12 ME12)	24	Washer 7/16" AN-960-716	12	997310	\$ 150.00
12, 14, 16		27	Rear Axle	2		
12, 14, 16		28	Locknut 7/16"-20 Nylon Insert Black Zinc	2		
12, 14, 16		29	R82 Brake Bracket	2		
12, 14, 16		30	Bolt 1/4"-20 X 1/2" Socket Button Head blak	4		
12, 14, 16		37	Split Lock Washer 1/4" ID	4		