

# R82 Cricket

## USD Retail Price List/ Order Form

Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_


Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### US Configuration

Effective: 01.01.22 - Rev 6

- Cricket
- Cushion
- Arm rests
- Head support
- Tray
- Grip bar
- Folding top
- Transport fittings



- A (Back Height) \_\_\_\_\_
- B (Arm Height) \_\_\_\_\_
- C (Lower Leg Length) \_\_\_\_\_
- D (Seat Depth) \_\_\_\_\_
- E (Chest Width) \_\_\_\_\_
- F (Seat Width) \_\_\_\_\_
- G (Height from Seat to Top of Head) \_\_\_\_\_
- H (Shoulder Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_

### Measurements



#### Size 1 Inches

#### Size 2 Inches

- Effective seat width (A)
- Effective seat depth (B)
- Backrest height (C1)
- Seat surface height at front edge (E)
- Footrest to seat distance (D)
- Leg to seat surface angle
- Backrest width, outside (F1)
- Backrest width, inside (F)
- Overall width (G)
- Length (H)
- Overall height (I)
- Folded width
- Folded height
- Folding length
- Backrest angle
- Minimum turning radius
- Total mass
- Max load
- Max load at transportation
- Required width of right angled corridor
- Static stability (sideways, downhill, uphill)

- 13½"
- 8½-12¼"
- 19¾"
- 19¾-21"
- 7-14¼"
- 360°
- 13½"
- 10¾"
- 20¾"
- 29½"
- 39¼"
- 14"
- 19"
- 30"
- 110°-125°
- 27½"
- 27.3 lbs.
- 77 lbs.
- 59½ lbs.
- 32"
- 10°

- 15½"
- 12¼-17"
- 22½"
- 21¼-22¼"
- 8¾-14½"
- 360°
- 15½"
- 13
- 22¾"
- 29½"
- 41¼"
- 14"
- 19"
- 30"
- 110°-125°
- 30"
- 29½ lbs.
- 110 lbs.
- 94½ lbs.
- 32"
- 10°

### Notes:

**Cricket Standard Features**




- Back Recline
- Adjustable Seat Depth
- Height and Angle Adjustable Foot Support
- Armrest, Fixed (Pr.)

Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Cricket sz 1	87043100	\$1,993	E1236
<input type="checkbox"/> Cricket sz 2	87043200	\$1,993	E1236

Description	Size	Part No.	USD Price	HCPCS
-------------	------	----------	-----------	-------

**Required, must select a color**

- Blue = 1                       Red = 2
- Pink = 5                       Grey = 4



<input type="checkbox"/> Cushion	1	8704311x	\$0	
<input type="checkbox"/> Cushion	2	8704321x	\$0	
<b>◆ Colors Green, Orange &amp; Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item</b>				
<input type="checkbox"/> Seat Side Parts ◆	One Size	87043130	\$54	
<b>◆ Required with Seat Width Reduction Kit</b>				



Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Transit Option		One Size	87043188	\$209	




Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Armrest, Fixed (Pr.) ◆		1	87043120-1	\$0	
<input type="checkbox"/> Armrest, Fixed (Pr.) ◆		2	87043120-2	\$0	
<b>◆ Required with Grip Bar, In Armrests</b>					



Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Seat Width Reduction Kit, Wedges (Pr.) ◆			87043130-1	\$75	K0108
<input type="checkbox"/> Seat Width Reduction Kit, Wedges (Pr.) ◆			87043130-2	\$75	K0108
<b>◆ Requires Seat Side Parts (87043130)</b>					















<input type="checkbox"/> Wedges for Shaping Back (Pr.)	4½" x 6" x 11"	Type 1	87043140-1	\$75	K0108
<input type="checkbox"/> Wedges for Shaping Back (Pr.)	5¾" x 6" x 11"	Type 2	87043140-2	\$75	K0108








<input type="checkbox"/> Chest Belt, Upholstery (Ea.)	A:11¼ x B:7½ x C:9¼-12"	1	85428-1	\$282	E0960
<input type="checkbox"/> Chest Belt, Upholstery (Ea.)	A:14 x B:10 x C:10-12¾"	2	85428-2	\$282	E0960



<input type="checkbox"/> Torso vest (Ea.)	S-P: 7¾-10¾"	XS	85705-XS	\$209	E0960
<input type="checkbox"/> Torso vest (Ea.)	S-P: 10¾-12½"	S	85705-S	\$209	E0960
<input type="checkbox"/> Torso vest (Ea.)	S-P: 12½-16"	M	85705-M	\$209	E0961
<input type="checkbox"/> Torso vest (Ea.)	S-P: 16-18½"	L	85705-L	\$209	E0962

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> 2-Point Pelvic Belt, Black	L: 13-14½"	1	85424-1	\$101	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 6¾-12½"	XS	85725-XS	\$140	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 8¼-13¾"	S	85725-S	\$140	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 9¾-16¼"	M	85725-M	\$140	E0978
	<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs ♦ ♦Required for use of Vests		One Size	862109-35	\$44	
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>USD Price</b>	<b>HCPCS</b>
	<input type="checkbox"/> Pelvic Harness	A:6½ x B:7½"	0	89108-0	\$255	E0978
	<input type="checkbox"/> Pelvic Harness	A:8½ x B:8¾"	1	89108-1	\$255	E0978
	<input type="checkbox"/> Pelvic Harness	A:10½ x B:9¾"	2	89108-2	\$255	E0978
	<input type="checkbox"/> Heel Straps (Pr.)		One Size	87043185	\$56	E0951
	<input type="checkbox"/> Upholstered Ankle Straps	A: 5½-7 x B:8½-10"	1	85426-1	\$193	K0108
	<input type="checkbox"/> Upholstered Ankle Straps	A: 6¼-7¾ x B:10-11¼"	2	85426-2	\$193	K0108
<b>Other types of vests and belts available, please see Vest and Belts Order Form</b>						
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>USD Price</b>	<b>HCPCS</b>
	<input type="checkbox"/> Headrest Cushion, Black (Ea.)		One Size	87043150	\$209	E0955
	<input type="checkbox"/> Sculpted Headrest Cushion, Black (Ea.)		One Size	87043151	\$214	E0955
	<input type="checkbox"/> Backrest Height Extension (Ea.)		1	87043160-1	\$57	K0108
	<input type="checkbox"/> Backrest Height Extension (Ea.)		2	87043160-2	\$57	K0108
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>USD Price</b>	<b>HCPCS</b>
	<input type="checkbox"/> Calf Support Strap (Ea.)		1	87043190-1	\$64	E0995
	<input type="checkbox"/> Calf Support Strap (Ea.)		2	87043190-2	\$64	E0995
	<input type="checkbox"/> Additional Footplate Locking Mechanism (Ea.)		One Size	87043186	\$32	
	<input type="checkbox"/> Grip Bar, In Armrests ♦		1	87043181-1	\$179	
	<input type="checkbox"/> Grip Bar, In Armrests ♦ ♦Requires Armrest		2	87043181-2	\$179	
	<input type="checkbox"/> Grip Bar, In Frame		1	87043182-1	\$179	
	<input type="checkbox"/> Grip Bar, In Frame		2	87043182-2	\$179	
	<input type="checkbox"/> Caster Locks		One Size	87043187	\$98	K0073

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<b>Required, must select a color</b>					
	<b>x = Color</b>					
	<input type="checkbox"/> Blue = 1	<input type="checkbox"/> Red = 2				
	<input type="checkbox"/> Pink = 5	<input type="checkbox"/> Grey = 4				
	<input type="checkbox"/> Folding Top (Canopy) (Ea.) ♦		1	8704317x	\$238	
	<input type="checkbox"/> Folding Top (Canopy) (Ea.) ♦		2	8704327x	\$238	
	<b>♦Required w/ Rain Cover</b>					
	<input type="checkbox"/> Underseat Storage Basket ♦		1	87043184-1	\$73	
	<input type="checkbox"/> Underseat Storage Basket ♦		2	87043184-2	\$73	
	<b>♦22 lbs. weight limit</b>					
	<input type="checkbox"/> Tray ♦		1	87043180-1	\$352	E0950
	<input type="checkbox"/> Tray ♦		2	87043180-2	\$352	E0950
	<b>♦Must select same size as armrest</b>					
	<input type="checkbox"/> Storage Strap		One Size	85100	\$32	
	<input type="checkbox"/> Rain Cape with Hood, Clear	A:11 x B:29½"	4	87465-410	\$147	
	<input type="checkbox"/> Rain Cape with Hood, Clear	A:11½ x B:34½"	5	87465-510	\$147	
	<input type="checkbox"/> Rain Cover, Clear (Ea.) ♦		2	87467-210	\$120	
	<input type="checkbox"/> Rain Cover, Clear (Ea.) ♦		3	87467-310	\$120	
	<b>♦Requires Folding Top (Canopy)</b>					

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

The Cricket complies with ISO7176-19, for safe transportation in motor vehicles. The R82 Cricket was tested with the Cricket back rest height extension and Y-hip belt and R82 transit anchors. The performance and safety of accessories and/or other seating options not tested by R82 is unknown, and no warranty or claim about their performance is or can be made.

When being transported in a vehicle, it is R82's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair; therefore R82 cannot recommend a transit system.

**Prices in USD are subject to change without notice**  
**For Quotes and Orders please send completed Order Form to: sales.us@etac.com or fax: 1-800-760-0537**

OF5100-ENG-REV.06

Etac  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
www.etac.us.com

  
Creating Possibilities