

# R82 Meerkat

## CAD Retail Price List/ Order Form

Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:


### US Configuration

Effective: 01.01.22 - Rev 5

Meerkat (Baseplate and Column)


**Accessories**

- Chest Support, Adjustable
- Pelvic Support, Adjustable
- Knee Supports
- Sandals
- Rocker
- Tray



- A (Total Height) \_\_\_\_\_
- B (Armpit Height) \_\_\_\_\_
- C (Top of Pelvic Height) \_\_\_\_\_
- D (Chest Circumference) \_\_\_\_\_
- E (Pelvic Circumference) \_\_\_\_\_
- F (Knee Circumference) \_\_\_\_\_
- G (Foot Length) \_\_\_\_\_
- H (Foot Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_

### Measurements



### Product Specs

	Size 1	Size 2	Size 3
Base Width (A)	20 <sup>3</sup> / <sub>4</sub> "	30 <sup>3</sup> / <sub>4</sub> "	30 <sup>3</sup> / <sub>4</sub> "
Base Length (B)	25 <sup>1</sup> / <sub>4</sub> "	33 <sup>3</sup> / <sub>4</sub> "	33 <sup>3</sup> / <sub>4</sub> "
Product Height (C)	14 <sup>3</sup> / <sub>4</sub> -29 <sup>1</sup> / <sub>2</sub> "	19 <sup>3</sup> / <sub>4</sub> -39 <sup>1</sup> / <sub>2</sub> "	25 <sup>1</sup> / <sub>2</sub> -51 <sup>1</sup> / <sub>4</sub> "
Product Weight	19 <sup>3</sup> / <sub>4</sub> lbs.	33 <sup>1</sup> / <sub>2</sub> lbs.	34 <sup>3</sup> / <sub>4</sub> lbs.

### User Measurement Guide for Stander

	Size 1	Size 2	Size 3
Approx. Age	1-4 Years	4-10 Years	6-18 Years
Max. User Weight/Load	44 lbs.	143 <sup>1</sup> / <sub>4</sub> lbs.	198 <sup>1</sup> / <sub>2</sub> lbs.
Approx. Total User Height (A)	27 <sup>1</sup> / <sub>2</sub> -41 <sup>1</sup> / <sub>4</sub> "	37 <sup>1</sup> / <sub>2</sub> -53 <sup>1</sup> / <sub>4</sub> "	46-70 <sup>3</sup> / <sub>4</sub> "
Max Armpit Height (B)	15 <sup>3</sup> / <sub>4</sub> -29 <sup>1</sup> / <sub>2</sub> "	28 <sup>1</sup> / <sub>4</sub> -39 <sup>1</sup> / <sub>4</sub> "	33-51 <sup>1</sup> / <sub>4</sub> "
Approx. Top of Pelvis (C)	13 <sup>3</sup> / <sub>4</sub> -21 <sup>3</sup> / <sub>4</sub> "	19 <sup>3</sup> / <sub>4</sub> -31 <sup>1</sup> / <sub>2</sub> "	25 <sup>1</sup> / <sub>2</sub> -41 <sup>1</sup> / <sub>4</sub> "


### User Measurement Guide for Accessories

	Small	Medium	Large	XL
Chest Circumference (D)	17 <sup>3</sup> / <sub>4</sub> -23 <sup>1</sup> / <sub>2</sub> "	23 <sup>1</sup> / <sub>2</sub> -29 <sup>1</sup> / <sub>4</sub> "	27 <sup>1</sup> / <sub>2</sub> -37 <sup>1</sup> / <sub>2</sub> "	37 <sup>1</sup> / <sub>2</sub> - 43 <sup>1</sup> / <sub>4</sub> "
Pelvic Circumference (E)	17 <sup>3</sup> / <sub>4</sub> -23 <sup>1</sup> / <sub>2</sub> "	23 <sup>1</sup> / <sub>2</sub> -29 <sup>1</sup> / <sub>4</sub> "	27 <sup>1</sup> / <sub>2</sub> -37 <sup>1</sup> / <sub>2</sub> "	37 <sup>1</sup> / <sub>2</sub> - 43 <sup>1</sup> / <sub>4</sub> "
Knee Circumference (F)	7 <sup>1</sup> / <sub>2</sub> -10"	10-11 <sup>1</sup> / <sub>2</sub> "	11 <sup>1</sup> / <sub>2</sub> -15 <sup>1</sup> / <sub>4</sub> "	
Knee Width	2 <sup>1</sup> / <sub>4</sub> "	3 <sup>1</sup> / <sub>4</sub> "	4"	
Heel Width (G)	2 <sup>1</sup> / <sub>4</sub> -3"	2 <sup>1</sup> / <sub>2</sub> -3 <sup>1</sup> / <sub>4</sub> "	3 <sup>1</sup> / <sub>4</sub> -3 <sup>3</sup> / <sub>4</sub> "	
Forefoot Width (H)	3 <sup>3</sup> / <sub>4</sub> -4 <sup>1</sup> / <sub>4</sub> "	4-5 <sup>3</sup> / <sub>4</sub> "	4 <sup>3</sup> / <sub>4</sub> -6 <sup>1</sup> / <sub>2</sub> "	

The Meerkat is measured in standard set-up and with the center column in vertical position. Different settings may result in another set of measurements. To accomplish some of the measurements above, a slight adjustment might be required. We supply various types of accessories, which is why the weight is listed without.

**NOTES:**

**Meerkat Standard Features**














- Base Plate with Wheels, Column and Tool
- Rocker + Handle


Description	Size	Part No.	CAD Price
<input type="checkbox"/> Meerkat, Pastel Blue, size 1 Plate, size 1 Column	1	870011-1R	\$1,261
<input type="checkbox"/> Meerkat, Strawberry Pink, size 1 Plate, size 1 Column ●●	1	870011-2R	\$1,261
<input type="checkbox"/> Meerkat, Grey, size 2 Plate, size 2 Column	2	870022-R	\$1,448
<input type="checkbox"/> Meerkat, Grey, size 2 Plate, Grey, size 3 Column	3	870023-R	\$1,521

●●*Custom order; Approx. lead time: 6-8 weeks. Non-returnable item*

Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Pelvic Support	5½" – 6½" (width) 17¾ - 23½" (circumference)	Small	869240	\$527
	<input type="checkbox"/> Pelvic Support	6½" – 8½" (width) 23½ - 29" (circumference)	Medium	869241	\$532
	<input type="checkbox"/> Pelvic Support	8½" – 12"(width) 27½ - 37½" (circumference)	Large	869242	\$550
	<input type="checkbox"/> Pelvic Support	12" – 14¼ (width) 37½ - 43 ¼" (circumference)	XL	869243	\$562
	<input type="checkbox"/> Chest Support	5½" – 6½" (width) 17¾ - 23½" (circumference)	Small	869480	\$527
	<input type="checkbox"/> Chest Support	6½" – 8½" (width) 23½ - 29" (circumference)	Medium	869481	\$532
	<input type="checkbox"/> Chest Support	8½" – 12"(width) 27½ - 37½" (circumference)	Large	869482	\$550
	<input type="checkbox"/> Chest Support	12" – 14¼ (width) 37½ - 43 ¼" (circumference)	XL	869483	\$562

Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Knee Supports	7½ - 9¾"	Small	869840	\$550
	<input type="checkbox"/> Knee Supports	9¾ - 11½"	Medium	869841	\$564
	<input type="checkbox"/> Knee Supports	11½ - 15¼"	Large	869842	\$564
	<input type="checkbox"/> Sandals for size 1 Plate	2¼ - 3" (Heel Width) 3¾ - 4¼" (Forefoot Width)	Small	869490	\$277
	<input type="checkbox"/> Sandals for size 2 Plate	2½ - 3¼" (Heel Width)	Medium	869491	\$290
	<input type="checkbox"/> Sandals for size 2 Plate	4 - 5¾" (Forefoot Width) 3¼ - 3¾" (Heel Width)	Large	869492	\$299
<hr/>					
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>CAD Price</b>
	<input type="checkbox"/> Tray (Only Anterior Use)		One Size	869485	\$416
	<input type="checkbox"/> Displacement kit bracket for mounting on column for pelvic and chest support		One Size	124025	\$94
	<input type="checkbox"/> Daily adjustment kit applicable on knee support <i>*For width adjustment on knee support</i>		One Size	124026	\$48
	<input type="checkbox"/> Daily adjustment kit, height ** <i>**For height adjustment of chest, pelvic or knee support</i>		One Size	124027	\$37
	<input type="checkbox"/> Padding for Meerkat Ratchet Straps (Pr.)		One Size	117818-PADS	\$65

Spare Parts	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Chest Belt		Small	869240-1	\$326
	<input type="checkbox"/> Chest Belt		Medium	869241-1	\$326
	<input type="checkbox"/> Chest Belt		Large	869242-1	\$326
	<input type="checkbox"/> Mounting for Ski Strap Buckle, Chest and Pelvic Support		One Size	117818	\$60
	<input type="checkbox"/> Mounting for Ski Strap Buckle, Knee Support		One Size	117818-1	\$60
	<input type="checkbox"/> Knee Support, Cup		Small	135310	\$136
	<input type="checkbox"/> Knee Support, Cup, Right		Medium	135311	\$136
	<input type="checkbox"/> Knee Support, Cup, Left		Medium	135312	\$136
	<input type="checkbox"/> Knee Support, Cup, Right		Large	135313	\$136
	<input type="checkbox"/> Knee Support, Cup, Left		Large	135314	\$136
	<input type="checkbox"/> Knee Support, Strap and Pad		Small	128254-1	\$150
	<input type="checkbox"/> Knee Support, Strap and Pad		Medium	128255-1	\$150
	<input type="checkbox"/> Knee Support, Strap and Pad		Large	128256-1	\$150
	<input type="checkbox"/> Sandal Strap		Small,Medium	127085	\$56
	<input type="checkbox"/> Sandal Strap		Large	127086	\$56

Spare Parts	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Column Mount Complete		One size	117832	\$127
	<input type="checkbox"/> Rods and Screws for Rocker		1	118598-2	\$163
	<input type="checkbox"/> Rods and Screws for Rocker		2	117047-2	\$163
	<input type="checkbox"/> Wheels (Airplane)		1	3035000	\$59
	<input type="checkbox"/> Wheels (Cow)		1	3035010	\$59
	<input type="checkbox"/> Wheels (Grey)		2	3006825	\$79

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

**Prices in CAD & is subject to change without notice**  
**For Quotes and Orders please send completed Order Form to: [us.info@r82.com](mailto:us.info@r82.com)**

Etac  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF5702-ENG-REV.05

