

# R82 Pony - E8002

## USD Retail Price List/ Order Form


Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	Phone:

### US Configuration

**Effective: 01.01.22 - Rev 4**

- Pony
- Head Support
- Abdominal Pad
- Hip Support
- Arm Support



**A** (Total Height) \_\_\_\_\_  
**B** (Shoulder Height) \_\_\_\_\_  
**C** (Chest Height) \_\_\_\_\_  
**D** (Buttock Height) \_\_\_\_\_  
**E** (Chest Width) \_\_\_\_\_  
**F** (Shoulder Width) \_\_\_\_\_  
 Age \_\_\_\_\_  
 Weight \_\_\_\_\_

### Measurements



	<b>Size 0 Inches</b>	<b>Size 1 Inches</b>	<b>Size 2 Inches</b>
Seat height (L)	9¼ - 14¾"	14¾ - 19½"	18¾ - 25"
Width (G)	23"	27"	30"
Length (H)	26¼"	28"	40½"
Height (I)	20¼ - 27"	23½ - 32"	26½ - 40½"
Angle range (A)	32°	32°	32°
Weight	15½ lbs.	22 lbs.	30 lbs.
Max. load	66 lbs.	88 lbs.	110 lbs.


### Pony Standard Features



- Rear Drag Wheels with Directional Stabilizer
- Leg & Foot Separator
- Height Adjustable Seat with Spring
- Height Adjustable Removable Back Support
- Height Adjustable Chest Support
- 5 Wheels (1 Rear & 2 Front with Locks)
- Height & Depth Adjustable Handle Bar

### NOTES:

Pony	Description	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Pony Walking Aid, Yellow	0	86400	\$2,425	E8002
	<input type="checkbox"/> Pony Walking Aid, Purple	1	86401	\$2,425	E8002
	<input type="checkbox"/> Pony Walking Aid, Lilac	2	86402	\$2,425	E8002

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Hip & Thigh Supports (Pr.)	6" x 3"	0	86408	\$344	
	<input type="checkbox"/> Hip & Thigh Supports (Pr.)	6" x 4"	1	86410	\$344	
	<input type="checkbox"/> Hip & Thigh Supports (Pr.) ••	6" x 4"	2	86412	\$344	

••Custom Order; Approx. lead time: 6-8 weeks

Minimum/Maximum Width	Size 0	Size 1	Size 3
	6" x 8¾"	6" x 8¾"	6" x 10¾"

Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Abdominal Cushion (Ea.) ••		1	86420	\$204	
<input type="checkbox"/> Abdominal Cushion (Ea.) ••		2	86422	\$204	

••Custom Order; Approx. lead time: 6-8 weeks




<input type="checkbox"/> Armrest/Tray (Ea.)		0	86427	\$506	
<input type="checkbox"/> Armrest/Tray (Ea.)		1	86428	\$506	
<input type="checkbox"/> Armrest/Tray (Ea.) ••		2	86429	\$506	

••Custom Order; Approx. lead time: 6-8 weeks



<input type="checkbox"/> Vertical Handle for Mounting in the Middle (Pr.)		One Size	86432	\$120	
---	--	----------	-------	-------	--

Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Head Support No. 1, Black PU		One Size	6171550	\$214	
<input type="checkbox"/> Head Support No. 2, Black PU		One Size	6172550	\$214	
<input type="checkbox"/> Head Support No. 3, Black PU		One Size	6173550	\$214	



<input type="checkbox"/> Anatomic Head Support ♦	A:5½ x B:5½ x C:9½"	Small	5099001	\$266	E0955
<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	N/C	
<input type="checkbox"/> Anatomic Head Support ♦	A:6 x B:5½ x C:10"	Medium	5099002	\$266	E0955
<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	N/C	
<input type="checkbox"/> Anatomic Head Support ♦	A:6½ x B:5½ x C:12½"	Large	5099003	\$266	E0955
<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	N/C	

♦ Turtle Bar System Short or Long Required; 98619-25 or 98619-26

<input type="checkbox"/> Headband System for Anatomic Head Support Only		One Size	99500	\$240	E0955
---	--	----------	-------	-------	-------



Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------








**Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable**

<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$240	E1028
<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$240	E1028

♦ Required with Anatomic Head Support; 5099001, 5099002 or 5099003

<input type="checkbox"/> Connecting Rod for Head Support 7¾"		1	861424-3	\$269	
<input type="checkbox"/> Connecting Rod for Head Support 13"		2	861424-4	\$269	

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Leg Protection		One Size	86450	\$230	
	<input type="checkbox"/> Push Handle		One Size	86416	\$266	

Spare Parts	Description	Measurements	Size	Part No.	USD Price	Qty
	<input type="checkbox"/> Chest Support, Complete		0	86310	\$532	_____
	<input type="checkbox"/> Chest Support, Complete		1	86311	\$532	_____
	<input type="checkbox"/> Chest Support, Complete		2	86312	\$532	_____
	<input type="checkbox"/> Straps for Chest Support		0	5540250	\$38	_____
	<input type="checkbox"/> Straps for Chest Support		1	5540300	\$48	_____
	<input type="checkbox"/> Straps for Chest Support		2	5540350	\$58	_____
	<input type="checkbox"/> Handle		0	86320	\$120	_____
	<input type="checkbox"/> Handle		1	86321	\$120	_____
	<input type="checkbox"/> Handle		2	86322	\$120	_____
	<input type="checkbox"/> Seat		0	86330	\$260	_____
	<input type="checkbox"/> Seat		1	86331	\$260	_____
	<input type="checkbox"/> Seat		2	86332	\$260	_____
	<input type="checkbox"/> Back Support		0	86434	\$214	_____
	<input type="checkbox"/> Back Support		1, 2	86435	\$214	_____
	<input type="checkbox"/> Directional Stabilizer for Rear Wheel		0	86438	\$91	_____
	<input type="checkbox"/> Directional Stabilizer for Rear Wheel		1	86440	\$91	_____
	<input type="checkbox"/> Directional Stabilizer for Rear Wheel		2	86442	\$91	_____
	<input type="checkbox"/> Directional Stabilizer for Front Wheel		0, 1, 2	86448	\$386	_____
	<input type="checkbox"/> Drag Wheel		One Size	86418	\$85	_____

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Prices in USD & is subject to change without notice  
For Quotes and Orders please send completed Order Form to: [us.info@r82.com](mailto:us.info@r82.com)

Etac  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
www.etac.us.com

OF5900-ENG-REV.04  
  
Creating Possibilities