

R82 Wombat Living

CAD Retail Price List/ Order Form


Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

| | |
|-----------------|-----------------|
| Bill To: | Ship To: |
| Address: | Address: |
| City/State/Zip: | City/State/Zip: |
| Phone: | Phone: |
| Email: | |

US Configuration



Effective: 01.01.22 - Rev 7

- Wombat Living
- Head Support
- Side Support
- Armrest
- Hip Support
- Pommel
- Vests and Belts
- Tray
- Helper's Handle



- A** (Back Height) _____
- B** (Arm Height) _____
- C** (Lower Leg Length) _____
- D** (Seat Depth) _____
- E** (Chest Width) _____
- F** (Seat Width) _____
- G** (Height from Seat to Top of Head) _____
- H** (Shoulder Width) _____
- Age _____
- Weight _____

Measurements



| | Size 1 Inches | Size 2 Inches | Size 3 Inches |
|---|--|--|--|
| Effective seat width (A) | 11 ³ / ₄ " | 14 ¹ / ₂ " | 16 ¹ / ₂ " |
| Seat width between hip supports | 4 ¹ / ₄ -11 ³ / ₄ " | 6 ³ / ₄ -14 ¹ / ₂ " | 8 ¹ / ₂ -16 ¹ / ₂ " |
| Seat width, pelvis support w/sides | 6 ³ / ₄ -9 ¹ / ₂ " | 9-11 ³ / ₄ " | 11 ¹ / ₂ -14 ¹ / ₄ " |
| Width, pelvis support without sides | 9 ¹ / ₂ " | 9 ¹ / ₂ " | 11 ¹ / ₂ " |
| Effective seat depth (B) | 7-12 ¹ / ₄ " | 9 ³ / ₄ -15" | 12 ¹ / ₂ -19 ³ / ₄ " |
| Backrest height (C) | 11-13 ¹ / ₂ " | 12 ¹ / ₄ -15" | 15 ¹ / ₄ -19 ³ / ₄ " |
| Backrest height (C) (high back) | 14 ¹ / ₄ -17" | 15 ¹ / ₂ -18" | 18-22 ³ / ₄ " |
| Footrest to seat distance (D) | 6-10" | 6-13 ³ / ₄ " | 6-17 ¹ / ₄ " |
| Seat surface height at front edge (4" Castors) | 15 ¹ / ₄ -26 ³ / ₄ " | 15 ¹ / ₄ -26 ³ / ₄ " | 16 ¹ / ₂ -27 ¹ / ₂ " |
| Back width, flat cushion (F) | 9 ³ / ₄ " | 11 ³ / ₄ " | 13 ³ / ₄ " |
| Back width, contoured cushion, between side supports (F) | 6 ³ / ₄ " | 9" | 10 ¹ / ₂ " |
| Arm rest height above seat cushion long rod and pad | 2 ¹ / ₂ -7" | 2 ¹ / ₂ -7" | 2 ¹ / ₂ -7" |
| Backrest angle (a1) | 90° to 130° | 90° to 130° | 90° to 133° |
| Seat plane angle (a2) | -15° to 20° | -15° to 20° | -10° to 22° |
| Leg to seat surface angle (a3) | -30° to 30° | -30° to 30° | -30° to 30° |
| Overall width (G) | 21 ¹ / ₂ " | 21 ¹ / ₂ " | 24" |
| Length (H) | 23 ¹ / ₂ " | 23 ¹ / ₂ " | 28 ¹ / ₄ " |
| Max load | 99.2 lbs. | 99.2 lbs. | 154.3 bls. |
| Total mass | 34.1 lbs. | 34.1 lbs. | 48 ¹ / ₂ lbs. |

NOTES:

Wombat Living Standard Features

















- Tilt in Space
- Back Recline and Prone Angle
- Adjustable Seat Depth
- Adjustable Back Height
- Height and Angle Adjustable Foot Support
- Without Seat/Back Cushions





| Description | Pricing | |
|--|---------|-----------|
| | | 95XXXX-00 |
| Wombat Living, size 1 | \$4,064 | 1 |
| Wombat Living, size 2 | \$4,064 | 2 |
| Wombat Living, size 3 | \$4,064 | 3 |
| Wombat Living, size 1 with High Back●● | \$4,064 | 5 |
| Wombat Living, size 2 with High Back●● | \$4,064 | 6 |
| Wombat Living, size 3 with High Back●● | \$4,064 | 7 |
| ●● Custom Order; Approx. lead time: 6-8 weeks | | |
| Frame Color: Anthracite (Dark Grey) | \$0 | 2 |
| 4" Front/Rear Wheels | \$0 | 2 |
| 4" Front/Rear Wheels with central brake - upcharge●● | \$371 | 3 |
| ●● Custom Order; Approx. lead time: 6-8 weeks | | |
| w/Gas Spring: Height Pedal on Frame, Tilt Lever Under Seat | \$0 | 1 |
| w/Gas Spring: Height Pedal on Frame, Tilt Lever on Push Brace - upcharge | \$419 | 2 |
| w/power actuator: Powered height setting, tilt lever under seat - upcharge ●● | \$2,247 | 5 |
| w/power actuator: Powered height setting, tilt lever on push brace - upcharge ●● | \$2,247 | 6 |
| ●●Custom Order for frame size 1 & 2; Apporx. lead time: 6-8 weeks | | |

Important!

Not recommended for users with a lot of involuntary movements

| Accessories | Description | Size | Part No. | CAD Price |
|---|---|---|------------------------------------|-----------|
| | Required, must select a color | | | |
| | xx = Color | | | |
| | <input type="checkbox"/> Red = 24 | <input type="checkbox"/> Petroleum = 25 | <input type="checkbox"/> Grey = 28 | |
|  | <input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.) | 1 | 955031-2xx | \$182 |
| | <input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.) | 2 | 955032-2xx | \$182 |
| | <input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.) | 3 | 955033-2xx | \$182 |
|  | <input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.) | 1 | 955031-238 | \$312 |
| | <input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.) | 2 | 955032-238 | \$312 |
| | <input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.) | 3 | 955033-238 | \$312 |
|  | <input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.) | 1 | 955021-2xx | \$351 |
| | <input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.) | 2 | 955022-2xx | \$351 |
| | <input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.) | 3 | 955023-2xx | \$351 |

| Accessories | Description | Size | Part No. | CAD Price | |
|---|--|-----------------------|------------|-------------|-----------|
|  | <input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.) | 1 | 955021-238 | \$351 | |
| | <input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.) | 2 | 955022-238 | \$351 | |
| | <input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.) | 3 | 955023-238 | \$351 | |
|  | <input type="checkbox"/> Seat Cushion, Contoured w/ Incontinence Cover (Ea.) | 1 | 955021-237 | \$377 | |
| | <input type="checkbox"/> Seat Cushion, Contoured w/ Incontinence Cover (Ea.) | 2 | 955022-237 | \$377 | |
| | <input type="checkbox"/> Seat Cushion, Contoured w/ Incontinence Cover (Ea.) | 3 | 955023-237 | \$377 | |
| <hr/> | | | | | |
| | Description | Measurements | Size | Part No. | CAD Price |
|  | <input type="checkbox"/> Lower Back Ext, Black (Ea.) * | H:2½" x W:10½" | 3 | 95140-60 | \$126 |
| | <i>*Only available on Wombat Living sz 3</i> | | | | |
| | <input type="checkbox"/> Lumbar Foam Wedge (Ea.) | | 1 | 955030-1 | \$33 |
|  | <input type="checkbox"/> Lumbar Foam Wedge (Ea.) | | 2 | 955030-2 | \$33 |
| | <input type="checkbox"/> Lumbar Foam Wedge (Ea.) | | 3 | 955030-3 | \$33 |
| | <hr/> | | | | |
| | Description | Measurements | Size | Part No. | CAD Price |
|  | <input type="checkbox"/> Hip Support (Pr.) | H:3¼" x W:7½" | 1,2 | 95110-1 | \$259 |
| | <input type="checkbox"/> Covers for Hip Supports (Required) | | 1,2 | 95110-60 | \$0 |
| | <input type="checkbox"/> Hip Support (Pr.) | H:4" x W:10" | 3 | 95111-1 | \$260 |
| | <input type="checkbox"/> Covers for Hip Supports (Required) | | 3 | 95111-60 | \$0 |
|  | <input type="checkbox"/> Pelvis Support without Sides (Ea.) ♦ | | 1,2 | 95140-61 | \$300 |
| | <input type="checkbox"/> Pelvis Support without Sides (Ea.) ♦ | | 3 | 95140-63 | \$300 |
| | <i>♦ Only available on Wombat Living with High Back</i> | | | | |
|  | <input type="checkbox"/> Pelvis Support with Sides (Ea.) ♦ | | 1 | 95140-64 | \$541 |
| | <input type="checkbox"/> Pelvis Support with Sides (Ea.) ♦ | | 2 | 95140-65 | \$541 |
| | <input type="checkbox"/> Pelvis Support with Sides (Ea.) ♦ | | 3 | 95140-66 | \$541 |
| <i>♦ Only available on Wombat Living with High Back</i> | | | | | |
|  | <input type="checkbox"/> Pommel (Ea.) | H:3½" x W:3½" x D:3½" | 1 | 95130-1 | \$188 |
| | <input type="checkbox"/> Pommel (Ea.) | H:3½" x W:3½" x D:4½" | 2 | 95130-2 | \$188 |
| | <input type="checkbox"/> Pommel (Ea.) | H:3½" x W:3½" x D:6" | 3 | 95130-3 | \$188 |
| | <input type="checkbox"/> Pommel (Ea.) | H:3½" x W:3½" x D:7" | 4 | 95130-4 | \$188 |
|  | <input type="checkbox"/> Knee Abduction Pommel | | One Size | 95230-10 | \$221 |
| | <input type="checkbox"/> Padding for Knee Abduction Pommel | | One Size | 95430-10260 | \$84 |
| <hr/> | | | | | |
| | Description | Measurements | Size | Part No. | CAD Price |
|  | <input type="checkbox"/> Armrest, Short (Pr.) | | 1,2 | 95112-1 | \$456 |
| | <input type="checkbox"/> Armrest, Short (Pr.) | | 3 | 95114 | \$455 |
|  | <input type="checkbox"/> Armrest, Long Rod and Pad (Pr.) | | 1,2 | 95113-1 | \$455 |
| | <input type="checkbox"/> Armrest, Long Rod and Pad (Pr.) | | 3 | 95115 | \$455 |

| Accessories | Description | Size | Part No. | CAD Price | |
|---|---|-----------|----------|-----------|-------|
|  | <input type="checkbox"/> Fixed Side Supports w/Covers (Pr.) | 1,2 | 95121 | \$221 | |
| | <input type="checkbox"/> Swing-Away Lateral Support Bracket, RH (Ea.) | One Size | 83400-R | \$293 | |
| | <input type="checkbox"/> Swing-Away Lateral Support Bracket, LH (Ea.) | One Size | 83400-L | \$293 | |
| | <input type="checkbox"/> Swing-Away Lateral Support Bracket, Offset, RH (Ea.) | One Size | 83600-R | \$293 | |
| | <input type="checkbox"/> Swing-Away Lateral Support Bracket, Offset, LH (Ea.) | One Size | 83600-L | \$293 | |
|  | <input type="checkbox"/> Bar, Straight (Ea.) | 4¾" | One Size | 81932 | \$0 |
| | <input type="checkbox"/> Bar, Straight, Stainless (Ea.) | 3¾" | One Size | 82569 | \$0 |
|  | <input type="checkbox"/> Bar, Offset, 1" (Ea.) | 4¾" | One Size | 81968 | \$0 |
| | <input type="checkbox"/> Bar, Offset, 1¾" (Ea.) | 4¾" | One Size | 82000 | \$0 |
| | <input type="checkbox"/> Bar, Offset, 2½" (Ea.) | 4¾" | One Size | 81969 | \$0 |
|  | <input type="checkbox"/> Curved Cushion (Ea.) | 3½" x 7½" | One Size | 6175406 | \$136 |
| | <input type="checkbox"/> Curved Cushion (Ea.) | 2½" x 6¾" | One Size | 6175405 | \$136 |
| | <input type="checkbox"/> Curved Cushion (Ea.) | 2½" x 6" | One Size | 6175404 | \$136 |
| | <input type="checkbox"/> Standard Cushion (Ea.) | 6" x 4¾" | One Size | 6175403 | \$136 |
| | <input type="checkbox"/> Standard Cushion (Ea.) | 4" x 4" | One Size | 6175402 | \$136 |
| | <input type="checkbox"/> Standard Cushion (Ea.) | 2⅔" x 4½" | One Size | 6175401 | \$136 |

| Description | Measurements | Size | Part No. | CAD Price |
|-------------|--------------|------|----------|-----------|
|-------------|--------------|------|----------|-----------|



■ XX = 60 (Black)



| | | | | |
|---|------------------|----|----------|-------|
| <input type="checkbox"/> Cross Vest (Ea.) | S-P: 7¾–12½" | XS | 85720-XS | \$247 |
| <input type="checkbox"/> Cross Vest (Ea.) | S-P 10¼–13½" | S | 85720-S | \$247 |
| <input type="checkbox"/> Cross Vest (Ea.) | S-P: 13½ - 17¾" | M | 85720-M | \$247 |
| <input type="checkbox"/> Cross Vest (Ea.) | S-P: 17¾–20¾" | L | 85720-L | \$247 |
| <input type="checkbox"/> Cross Vest (Ea.) | S-P: 20¾ - 24 ½" | XL | 85720-XL | \$247 |







| | | | | |
|---|---------------|-----|-----------|-------|
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 7¾–10¾" | XS | 85705-XS | \$260 |
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 10¾–12½" | S | 85705-S | \$260 |
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 12½–16" | M | 85705-M | \$260 |
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 16–18½" | L | 85705-L | \$267 |
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 18½–20½" | XL | 85705-XL | \$260 |
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 20½–23½" | XXL | 85705-XXL | \$260 |



| | | | | |
|--|----------------------|---|---------|-------|
| <input type="checkbox"/> Chest Belt, Upholstered (Ea.) | A:11¼xB:7½xC:9¼-12" | 1 | 85428-1 | \$351 |
| <input type="checkbox"/> Chest Belt, Upholstered (Ea.) | A:14xB:10xC:10-12¾" | 2 | 85428-2 | \$351 |
| <input type="checkbox"/> Chest Belt, Upholstered (Ea.) | A:16½xB:12½xC:10-13" | 3 | 85428-3 | \$351 |



| | | | | |
|--|---------------|----|----------|-------|
| <input type="checkbox"/> 4-Point Padded Hip Belt, 1" | H-H: 6¾–12½" | XS | 85725-XS | \$176 |
| <input type="checkbox"/> 4-Point Padded Hip Belt, 1" | H-H: 8¼–13¾" | S | 85725-S | \$176 |
| <input type="checkbox"/> 4-Point Padded Hip Belt, 1" | H-H: 9¾–16¼" | M | 85725-M | \$175 |
| <input type="checkbox"/> 4-Point Padded Hip Belt, 1" | H-H: 12½–17¼" | L | 85725-L | \$175 |
| <input type="checkbox"/> 4-Point Padded Hip Belt, 1" | H-H: 16¼–20½" | XL | 85725-XL | \$175 |

| Accessories | Description | Measurements | Size | Part No. | CAD Price |
|---|--|--------------|----------|-----------|-----------|
|  | <input type="checkbox"/> 2-Point Pelvic Belt, Black | L: 13-14½" | 1 | 85424-1 | \$126 |
|  | <input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs ** **Required for use of Vests | | One Size | 862109-35 | \$55 |
| | <input type="checkbox"/> Y-Hip Strap, Upholstered, 1" End Girth (Ea.) | | 1 | 85427-1 | \$176 |
| | <input type="checkbox"/> Y-Hip Strap, Upholstered, 1" End Girth (Ea.) | | 2 | 85427-2 | \$176 |
|  | <input type="checkbox"/> Ankle straps | C: 6¼–9" | XS | 85745-XS | \$240 |
| | <input type="checkbox"/> Ankle straps | C: 9–12¼" | M | 85745-M | \$240 |
| | <input type="checkbox"/> Ankle straps | C: 11½–15" | XL | 85745-XL | \$241 |
|  | <input type="checkbox"/> Neoprene Sandals (Pr.) | | One Size | 85696-1 | \$169 |

Other types of vests and belts available, please see Vest and Belts Order Form

| Description | Measurements | Size | Part No. | CAD Price |
|-------------|--------------|------|----------|-----------|
|-------------|--------------|------|----------|-----------|

Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable

| | | | | |
|--|------------------|-------|----------|-------|
| <input type="checkbox"/> Turtle Bar System ♦ | A: 5", B: 4-8½" | Short | 98619-25 | \$299 |
| <input type="checkbox"/> Turtle Bar System ♦ | A: 10", B: 4-8½" | Long | 98619-26 | \$299 |

♦Required with Anatomic Head Support; 5099001, 5099002 or 5099003

| | | | | |
|--|----------------------|--------|---------|-------|
| <input type="checkbox"/> Anatomic Head Support ♦ | A:5½ x B:5½ x C:9½" | Small | 5099001 | \$332 |
| <input type="checkbox"/> Head Support Cover (Required) | | Small | 9923-83 | \$0 |
| <input type="checkbox"/> Anatomic Head Support ♦ | A:6 x B:5½ x C:10" | Medium | 5099002 | \$332 |
| <input type="checkbox"/> Head Support Cover (Required) | | Medium | 9945-83 | \$0 |
| <input type="checkbox"/> Anatomic Head Support ♦ | A:6½ x B:5½ x C:12½" | Large | 5099003 | \$332 |
| <input type="checkbox"/> Head Support Cover (Required) | | Large | 9967-83 | \$0 |

♦Turtle Bar System Short or Long Required; 98619-25 or 98619-26

XX = Turtle Bar Size

- 92 = Short Turtle Bar 93 = Long Turtle Bar

Adjustable Headrest without Dynamic Feature - Requires Bracket (see below)

| | | |
|--|------------|-------|
| <input type="checkbox"/> Adj. Width Headrest, Large Laterals, Black Cover ♦ | 99737xx-85 | \$397 |
| <input type="checkbox"/> Adj. Width Headrest, Large Laterals, w/Speakers, Black Cover ♦♦ | 99717xx-85 | \$527 |

♦♦Custom order; Approx. lead time: 6-8 weeks.
♦ Turtle Bar System Short or Long Required; x98619-25 or x98619-26

Adjustable Headrest with Dynamic Feature - Requires Bracket (see below)

| | | |
|--|------------|-------|
| <input type="checkbox"/> Adj. Width Headrest, Large Laterals, Black Cover ♦ | 99736xx-85 | \$501 |
| <input type="checkbox"/> Adj. Width Headrest, Large Laterals, w/Speakers, Black Cover ♦♦ | 99756xx-85 | \$501 |

♦♦Custom order; Approx. lead time: 6-8 weeks.
♦ Turtle Bar System Short or Long Required; x98619-25 or x98619-26






| | | | |
|---|----------|-------|-------|
| <input type="checkbox"/> Headband System for Anatomic Head Support Only | One Size | 99500 | \$299 |
|---|----------|-------|-------|



| | | | |
|---|----------|-------|------|
| <input type="checkbox"/> Universal Headrest Mount ♦ | One Size | 87590 | \$79 |
|---|----------|-------|------|

♦Required when using other manufactures Head Support

Allows most Whitmyer, Stealth, Otto Bock headrest, I2I to mount to our products (x:panda, Serval, etc.). Please contact Customer Service for confirmation

Other types of headrest available, please see Head Supports Order Form

| Accessories | Description | Measurements | Size | Part No. | CAD Price |
|---|---|--------------|----------|-----------|-----------|
| | <input type="checkbox"/> Directional Stabilizer Kit ●● ●● <i>Custom order; Approx. lead time: 6-8 weeks.</i> | | One Size | 89592-KIT | \$280 |
|  | <input type="checkbox"/> Tray (Ea.) | | 1,2 | 95125 | \$390 |
| | <input type="checkbox"/> Tray (Ea.) | | 3 | 95126 | \$390 |
| | <input type="checkbox"/> Tray, Grey | | 1,2 | 8976350 | \$293 |
|  | <input type="checkbox"/> Push Brace, Seat Mounted | | 1,2 | 95148-1 | \$299 |
| | <input type="checkbox"/> Push Brace, Seat Mounted | | 3 | 95148-3 | \$299 |
|  | <input type="checkbox"/> Angle Adjustable Push Brace (Ea.) | | 3 | 95147 | \$273 |
|  | <input type="checkbox"/> Helper's Handle, Fold Down (Pr.) | | One Size | 95145 | \$104 |
| | <input type="checkbox"/> Helper's Handle, Fixed (Pr.) | | One Size | 95146 | \$104 |
|  | <input type="checkbox"/> Anti-Tip Device | | One Size | 95953 | \$209 |

| Spare Parts | Description | Measurements | Size | Part No. | CAD Price |
|---|---|--------------|------|----------|-----------|
|  | <input type="checkbox"/> Wombat Living Step 1 to 2 Kit, Seat & Back | | 1 | 95560 | \$364 |
| | <input type="checkbox"/> Wombat Living Step 1 to 2 Kit, Seat & Back | | 2 | 95561 | \$364 |
|  | <input type="checkbox"/> Wombat Living Short to High Back Kit | | 1 | 95560-1 | \$205 |
| | <input type="checkbox"/> Wombat Living Short to High Back Kit | | 2 | 95561-1 | \$202 |

| Description | Size | Part No. | CAD Price |
|-------------|------|----------|-----------|
|-------------|------|----------|-----------|

Required, must select a color





xx = Color

Red = 24

Petroleum = 25

Grey = 28

| | | | |
|---|---|------------|-------|
| <input type="checkbox"/> Cover for Standard Back Cushion | 1 | 955611-2xx | \$140 |
| <input type="checkbox"/> Cover for Standard Back Cushion | 2 | 955612-2xx | \$140 |
| <input type="checkbox"/> Cover for Standard Back Cushion | 3 | 955613-2xx | \$140 |
| <input type="checkbox"/> Cover for Contoured Back Cushion, Grey | 1 | 955611-238 | \$177 |
| <input type="checkbox"/> Cover for Contoured Back Cushion, Grey | 2 | 955612-238 | \$177 |
| <input type="checkbox"/> Cover for Contoured Back Cushion, Grey | 3 | 955613-238 | \$177 |
| <input type="checkbox"/> Cover for Seat Cushions | 1 | 955621-2xx | \$140 |
| <input type="checkbox"/> Cover for Seat Cushions | 2 | 955622-2xx | \$140 |
| <input type="checkbox"/> Cover for Seat Cushions | 3 | 955623-2xx | \$140 |

| Spare Parts | Description | Size | Part No. | CAD Price |
|-------------|---|-------------|-----------------|------------------|
| | <input type="checkbox"/> Cover w/Foam for Pelvis Support without Sides | 1,2 | 95122-1 | \$168 |
| | <input type="checkbox"/> Cover w/Foam for Pelvis Support without Sides | 3 | 95122-2 | \$168 |
| | <input type="checkbox"/> Cover w/Foam for Pelvis Support with Sides | 1 | 95123-1 | \$200 |
| | <input type="checkbox"/> Cover w/Foam for Pelvis Support with Sides | 2 | 95123-2 | \$200 |
| | <input type="checkbox"/> Cover w/Foam for Pelvis Support with Sides | 3 | 95123-3 | \$200 |
| | <input type="checkbox"/> Covers for Fixed Side Supports (95121) | 1,2 | 95121-60 | \$110 |
| |  | | | |
| | <input type="checkbox"/> Footplate, Angle & Height Adj. 12" x 7½" | | 95145-1 | \$150 |
| | <input type="checkbox"/> Footplate, Angle & Height Adj. 14" x 9" | | 85145-2 | \$150 |
| |  | | | |
| | <input type="checkbox"/> Footplate, Angle & Height Adj. ♦ 14" x 9" | | 95145-2 | \$150 |
| | <input type="checkbox"/> Footplate, Angle & Height Adj. ♦ 16" x 9½" | | 95145-3 | \$150 |
| | <input type="checkbox"/> Footplate, Angle & Height Adj. ♦ 18¾" x 9½" | | 95145-4 | \$150 |
| | ♦ Includes 4" mounting fitting | | | |
| |  | | | |
| | <input type="checkbox"/> Footplate with Rod 7½" x 12" | 1 | 95144-1 | \$293 |
| | <input type="checkbox"/> Footplate with Rod 9" x 14" | 2 | 95144-2 | \$293 |
| | <input type="checkbox"/> Footplate with Rod 9¼" x 16" | 3 | 95144-3 | \$293 |
| | <input type="checkbox"/> Footplate with Rod 9½" x 19" | 4 | 95144-4 | \$293 |
| | Description | Size | Part No. | CAD Price |
| | <input type="checkbox"/> Batteries for Wombat Living with Power (Set) | One Size | 2222012 | \$397 |
| | ♦ Both batteries <u>must</u> be charged at the same time | | | |
| |  | | | |
| | <input type="checkbox"/> Charger & Plug for USA/Canada (Ea.) | One Size | 2222161 | \$416 |

Prices in CAD & is subject to change without notice

For Quotes and Orders please send completed Order Form to: sales.us@etac.com or Fax: 1-800-760-0537

OF6602-ENG-REV.07

Etac
 13137 Bleinheim Lane
 Matthews, North Carolina 28105
 +1 844-US-MOBILITY (844-876-6245)
www.etac.us.com


 Creating Possibilities