



Fixed Tilt CX Planar


Standard Specifications

MEASUREMENTS:

INCHES

HCPCS:E1236

- * - Plus or minus .5 inches due to fabric
- Δ - Weight of chair is without front rigging
- - Adjustable push handle is not available on 18 inch models
- ◇ - Models with Textilene Fabric, Shoulder Strap Height is 2 inches shorter
- ^ - Seat back height
- (C) - Cordura upholstery
- (T) - Textilene upholstery



	CX 10	CX 12	CX 14	CX 16	CX18
Seat Width*	10	12	14	16	18
Seat Depth	6 – 11	8 – 13	9 – 14	10 – 15	13 – 20
Seat Back Height^	18	18	29.5	27	26
Seat to Footplate*	5 – 28	7 – 22	6.5 – 26	7 – 25	13 – 25
Seat to Floor*	27.5	26	26	26	24.5
Seat-to-Back Angle	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°
Seat Angle	30°	30°	30°	30°	30°
Overall Height	37.5	38	40	42.5	38
Overall Length	53.5	44	49	60	40
Overall Width	21.25	23.25	25.5	27.5	28
Folded Length	47.5	48	47	48.5	48
Folded Width	17	15	16	19	15
Folded Height	16.5	17	17	16	16.5
Headrest Extension	5	4.5	5	8	8
Shoulder Strap Height ◇	12 & 18	11 & 18	13 & 24	15.5 & 22.5	16 & 22
Lateral Support Width Adjustment	3.5 – 10.5	5.5 – 12.5	7.5 – 14.5	9.5 – 16.5	11.5 – 18.5
Weight Capacity Regular/Transit	75/66 lbs	75/66 lbs	100/100 lbs	170/170 lbs	250/170 lbs
Weight of Chair Δ	36 lbs	36 lbs	38 lbs	45 lbs	30 lbs