



# Fixed-Tilt Scout

## Standard Specifications

### MEASUREMENTS:

INCHES

#### HCPCS:E1236

- \* - Plus or minus .5 inches due to fabric
- Δ - Weight of chair is without front rigging
- - Adjustable push handle is not available on 18 inch models

	SC 12	SC 14	SC 16	SC 18
Seat Width* 	12	14	16	18
Seat Depth 	8 – 13	10 – 15	11 – 16	14 – 21
Seat Back Height 	20	22.5	25	24
Seat to Footplate* 	4 – 25	3 – 23	5.5 – 22.5	10 – 24
Seat to Floor* 	22	22	22	23
Seat-to-Back Angle 	85°-90°-95°	85°-90°-95°	85°-90°-95°	85°-90°-95°
Seat Angle 	30°	30°	30°	30°
Overall Height 	37	38	39	39
Overall Length 	51	54	55	56
Overall Width 	24.5	26	28	30.5
Folded Length 	52.5	52	55	57
Folded Width 	16	16	16	16
Folded Height 	16.5	15.5	15	16.5
Headrest Extension 	8	8	8	8
Shoulder Strap Height 	12 – 18	14 – 20	12.5 – 22.5	16 – 22
Push Handle Adjustments ○ 	29 – 46	28 – 45	31 – 48	N/A
Weight Capacity Regular/Transit 	75/66 lbs	100/100 lbs	170/170 lbs	250/170 lbs
Weight of Chair Δ 	31 lbs	32 lbs	35.5 lbs	38 lbs