

# FIXED-TILT CRUISER



- **30° Fixed Tilt**  
Increased tilt angle improves postural control through upper body positioning
- **Customizable**  
Custom-built to your specifications to meet each individual's physical needs
- **Multiple Upholstery Styles & Colors**  
Removable, washable Cordura or Textilene available in various colors
- **5 Years of Growth Capabilities**  
Adjustable seat depth extends chair life for a growing child
- **Self-Tensio™**  
Built in self-tensioning seating system increases sitting tolerance



Positioning options & accessories sold separately

## Transit



Meets or exceeds WC19 safety standards for safe transportation

## Planar



Solid seat & back provides in-depth seating & therapeutic positioning

## Scout



All-terrain option built for maximum comfort and durability

## Push Handle



Two Piece Push Handle (Option)  
or  
One Piece Height Adjustable Push Handle (Standard)

**Lightweight • Compact Folding • Convenient**

 **Convaid**  
*Now You're Going Places*



888-Convaid (266-8243)  
www.convaid.com

### Accessories

- **Extended Headrest Cover with windows** **A**  
(Canopy)
- **Extended Headrest Cover without windows**  
(Canopy)
- **Headrest Cover**  
(Canopy)
- **Mosquito Net** **B**  
(Canopy Required)
- **Rain Cover**  
(Canopy Required)
- **Attendant Hand Brake** **C**
- **Upper Extremity Support Surface Tray** **D**  
(Clear Tray)
- **Utility Bag**
- **Saddle Bags** (Pair)
- **Wheelchair Travel Bag**
- **Under Seat Storage Basket** **E**
- **Rear Anti-Tip Tubes** **F**
- **Curb Tipper** **G**
- **Caster Locks** **H**  
(For 2" Wide Front Wheels)



### Head Supports

-  • **Headrest Extension**  
(Required for Transit)  
- Curved Cordura (Shown)  
- Rectangular Textilene
-  • **Occi Headwings**
-  • **Padded Headwings**

### Trunk Positioning

-  • **H-Harness with Padded Covers**  
(Required for Transit)
-  • **Soft Adjustable Lateral Support w/ Scoliosis Strap**  
(Single or Double Flap)
-  • **Reducer Seat Insert**  
(Shown)
- **Full Torso Support Vest**

### Pelvic Positioning

-  • **Support Cushion**
-  • **Align Cushion**
-  • **Position Cushion**
- **Solid Seat Insert**  
(Required for all Cushions)
- **Incontinence Cushion Liner**  
(Solid Seat Insert & Cushion Required)

### Frame Options

- **One Piece Height Adjustable Push Handle** (Standard)
- **Two Piece Push Handle**
- **Anatomic Back**  
(For Additional Lumbar Support)

### Transit Options

- **5-Point Harness** (CX10/12 Only)
- **Q'Straint Transit Lap Belt**
- **Sure-Lok Transit Lap Belt**

### Upholstery Options



(Zippered Pockets in Cordura Seat & Back Upholstery)

### Cruiser 18



Crash tested model shown

#### Standard Features

- **Two Piece Push Handle**
- **8" Front, 11" Rear Solid Tires**
- **Cordura or Textilene Upholstery**  
(Cordura Shown)
- **Toggle Wheel Locks**

### Upholstery Colors

#### Textilene

- Navy Blue**
- Teal Green**
- Violet Purple**
- Wild Raspberry**
- Cordura**
- Panther Black**
- Cherry Red**
- Sassy Purple**
- Princess Pink**
- Royal Blue**
- Forest Green**
- Chocolate Brown**
- Aqua Blue**
- Tangerine Orange**
- Apple Green**

### Foot & Leg Positioning

-  • **Medial Thigh Support Abductor**
-  • **Lateral Thigh Support Adductor**  
(Shown)
- **Calf Panel**
- **Foot Positioners**  
(Required for Transit)
- **Footplate Securement Strap**
- **Angle Adjustable Footplates**

### Compact Fold



• **Closure Strap**

### Tire Options

**cx** - Cruiser **cx sc** - Cruiser Scout Option

#### Front Tire

- **7.5" x 2" Front Pneumatic** (Rolko) ----- **cx** 10/12/14/16
- **8" x 2" Front Solid Ribbed** ----- **cx sc** 12/14/16/18 **cx** 18
- **8" x 2" Front Pneumatic Ribbed** ----- **cx sc** 12/14/16/18 **cx** 18

#### Rear Tire

- **12.5" x 2.5" Rear Solid Knobby** (Rolko) ----- **cx** 10/12/14/16
- **12.5" x 2.5" Rear Pneumatic Knobby** (Rolko) ----- **cx** 10/12/14/16
- **12.5" x 2.5" Rear Solid Knobby** ----- **cx** 18
- **12.5" x 2.5" Rear Pneumatic Knobby** ----- **cx** 18