



Upright EZ Convertible

Standard Specifications

MEASUREMENTS:

INCHES

HCPCS:E1236

- * - Plus or minus .5 inches due to fabric
- Δ - Weight of chair is without front rigging
- - Adjustable push handle is not available on 18 inch models

	CV 12	CV 14	CV 16	CV 18
Seat Width* 	12	14	16	18
Seat Depth 	9 – 14	11 – 16	12 – 17	14 – 19
Seat Back Height 	21C/19T	24C/22T	26C/24T	27C/25T
Seat to Footplate* 	3 - 20	3 – 21	3 – 22	3 – 22
Seat to Floor* 	17	18.5	21	21
Back Angle 	100°	100°	100°	100°
Seat Angle 	10°	10°	10°	10°
Overall Height 	37.5	37	37	49
Overall Length 	43	43	47.5	48.5
Overall Width 	25.5	26.5	28.5	32
Folded Length 	44	44	45.5	45
Folded Width 	12.5	12	15	14.5
Folded Height 	14.5	14.25	16	15.25
Headrest Extension 	8	8	8	10
Shoulder Strap Height 	11 – 17	14 – 20	15 – 21	16 – 22
Height of Rear Wheels 	20	22	24	24
Weight Capacity Regular/Transit 	75/66 lbs	100/100 lbs	170/170 lbs	250/170 lbs
Weight of Chair Δ 	31 lbs	36 lbs	39 lbs	36.5 lbs