



UPRIGHT EZ RIDER



- **10° Fixed Tilt**
Provides upright posture to enhance functional positioning
- **Customizable**
Custom-built to your specifications to meet each individual's physical needs
- **Multiple Upholstery Styles & Colors**
Removable, washable Cordura or Textilene available in various colors
- **5 Years of Growth Capabilities**
Adjustable seat depth extends chair life for a growing child
- **Self-Tensio®**
Built in self-tensioning seating system increases sitting tolerance



Positioning options & accessories sold separately

Transit



Meets or exceeds WC19 safety standards for safe transportation

Planar



Solid seat & back provides in-depth seating & therapeutic positioning

Convertible



Removable Mag Rear Wheels allows independent mobility

Push Handle



Two Piece Push Handle (Option)
or
One Piece Height Adjustable Push Handle (Standard)

Lightweight • Compact Folding • Convenient

Convaid
Now You're Going Places



844 US MOBILITY (844-876-6245)
www.convaid.com

Accessories

- Extended Headrest Cover with windows **A** (Canopy)
- Extended Headrest Cover without windows (Canopy)
- Headrest Cover (Canopy)
- Mosquito Net **B** (Canopy Required)
- Rain Cover (Canopy Required)
- Upper Extremity Support Surface Tray **C** (Clear Tray) (Armrests Required)
- Utility Bag **D**
- Wheelchair Travel Bag
- Rear Anti-Tip Tubes **E**
- Caster Locks (For 2" Wide Front Wheels)
- Mag Wheel Bag **F** (EZ CV Required)



Frame Options

- One Piece Height Adjustable Push Handle (Standard)
- Two Piece Push Handle
- Height Adjustable Flip-up Armrests
- 4" Extended Two Piece Push Handle

Transit Options

- 5-Point Harness (EZ12 Only)
- Q'Straint Transit Lap Belt
- Sure-Lok Transit Lap Belt
- Upper Extremity Support Surface Tray (Foam Tray) (Armrests Required)

Upholstery Options



(Zippered Pockets in Cordura Seat & Back Upholstery)

EZ Rider 18




Crash tested model shown

Standard Features

- Two Piece Push Handle
- 8" Front, 11" Rear Solid Tires
- Cordura or Textilene Upholstery (Cordura Shown)
- Toggle Wheel Locks

Head Supports



Trunk Positioning



Pelvic Positioning



Upholstery Colors


Textilene

- Navy Blue
- Teal Green
- Violet Purple
- Wild Raspberry

Cordura

- Panther Black
- Cherry Red
- Sassy Purple
- Princess Pink
- Royal Blue
- Forest Green
- Chocolate Brown
- Aqua Blue
- Tangerine Orange
- Apple Green

Foot & Leg Positioning



Compact Fold



Tire Options

EZ -EZ Rider **EZ CV** - EZ Rider Convertible Option

Front Tires

- 7.5" x 2" Front Pneumatic (Rolko)-----EZ 12/14/16
- 8" x 2" Front Solid Ribbed-----EZ 18
- 8" x 2" Front Pneumatic Ribbed -----EZ 18

Rear Tires

- 12.5" x 2.5" Rear Solid Knobby (Rolko)-----EZ 12/14/16
- 12.5" x 2.5" Rear Pneumatic Knobby (Rolko) -EZ 12/14/16
- 12.5" x 2.5" Rear Solid Knobby-----EZ 18
- 12.5" x 2.5" Rear Pneumatic Knobby-----EZ 18
- 20" Mag Rear Solid Tires -----EZ CV 12
- 22" Mag Rear Solid Tires -----EZ CV 14
- 24" Mag Rear Solid Tires -----EZ CV 16/18