

Convoid Carrot 3 Booster Kit - Carrot 3 Child Restraint System

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone: Fax:	City/State/Zip:
Email:	Phone:

Carrot 3 Booster Kit Standard Features:

- 4-Point Positioning Harness
- H-Harness Attaching Bracket
- Booster Seat Base
- Booster Metal Rods (Pair)
- LATCH Connector System

Standard Models:

- | | |
|---|----------------|
| <input type="checkbox"/> Carrot 3 Booster Kit CRT3BLACK - CRT3BK-19 | Price |
| <input type="checkbox"/> Carrot 3 Booster Kit CRT3BLUE - CRT3BK-29 | \$1,150 |
| <input type="checkbox"/> Carrot 3 Booster Kit CRT3RED - CRT3BK-39 | \$1,150 |

Carrot 3 Booster Kit Required Information

Carrot 3 Serial #: _____

Carrot 3 Booster Seat Recommended Item

- Tether strap (Black)

Price
\$89

Black	Blue	Red
CT3TS0-99	CT3TS0-99	CT3TS0-99

Carrot 3 Booster Kit Required Items*

- Backrest Extension - 2"
 Shoulder Protector Wings

Price
\$68
\$527

Black	Blue	Red
CT3BE2-19	CT3BE2-29	CT3BE2-39
CT3SP0-19	CT3SP0-29	CT3SP0-39

* Option to transform an existing Carrot 3 Seat to a Booster Seat

Accessories

- Booster Support Tray (Crash Tested)
 3D Body Balance Pads L/R Set
 Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°)
 Backrest Extension - 4" (Extends back height to 32.8")

Price
\$290
\$150
\$113
\$91

Black	Blue	Red
CRT3BT-19	CRT3BT-29	CRT3BT-39
CT3BBP-99	CT3BBP-99	CT3BBP-99
CT3SA0-19	CT3SA0-29	CT3SA0-39
CT3BE4-19	CT3BE4-29	CT3BE4-39

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-888-CONVAID

2830 California St., Torrance, CA 90503 • 1-888-CONVAID (266-8243), 310-618-0111 • Fax: 310-618-2166


Prices in USD & subject to change, ex works Torrance, California.

www.convoid.com

Effective: 5/2/2016

Rev 1

CLIENT MEASUREMENTS:



Carrot 3 Base Model Specifications:

Seat Width:	16.5"
Seat Depth:	17"
Backrest Height Adjustment:	28" - 33"
Head Support Height:	8.5"
Backrest Height:	11.8"
Backrest Extension Height (Standard):	2"
Backrest Extension Height (Option):	4"
Shoulder Wing Height (Standard):	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"
Product Weight:	21 lbs.
Recline Adjustment:	90° - 180°
Product Overall Width:	20" - 22"

User Measurement Guide:

Minimum User Weight:	79 lbs.
Maximum User Weight:	165 lbs.
Minimum User Height:	55"
Maximum User Height:	69"

2830 California St., Torrance, CA 90503 • 1-888-CONVAID (266-8243), 310-618-0111 • Fax: 310-618-2166