

# Convoid Rodeo Tilt - HCPCS E1232/E1161

## Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

### Rodeo Tilt Standard Features:

- 5 to 45 Degrees of Adjustable Tilt
- Silver Powder Coated Frame
- One Piece Height Adjustable Push Handle
- One Piece Angle Adjustable Footplate
- 2-Point Positioning Belt
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire
- Zippered Pockets in Seat Bottom Upholstery
- HCPCS code E1161 is only valid on size 16 (RD16)
- Quick Release Rear Wheels
- Single Foot Activated Wheel Locks
- Closure Strap
- Seat Depth Growth of 5"
- Elevating Legrest
- Cordura Upholstery
- Tool-less Upholstery Removal
- Made in the USA

### Standard Models:

	Price
<input type="checkbox"/> Rodeo RD10 - 904340	<b>\$2,468</b>
<input type="checkbox"/> Rodeo RD12 - 902689	<b>\$2,468</b>
<input type="checkbox"/> Rodeo RD14 - 902691	<b>\$2,468</b>
<input type="checkbox"/> Rodeo RD16 - 903134	<b>\$2,468</b>

### Upholstery Colors

	Price
<input type="checkbox"/> Royal Blue	\$0
<input type="checkbox"/> Forest Green	\$0
<input type="checkbox"/> Cherry Red	\$0
<input type="checkbox"/> Sassy Purple	\$0
<input type="checkbox"/> Princess Pink	\$0
<input type="checkbox"/> Panther Black	\$0
<input type="checkbox"/> Chocolate Brown	\$0
<input type="checkbox"/> Aqua Blue	\$0
<input type="checkbox"/> Apple Green	\$0
<input type="checkbox"/> Tangerine Orange	\$0

### Transit Required Items

	Price	HCPCS
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	<b>\$200</b>	
<input type="checkbox"/> Headrest Extension	<b>\$130</b>	<b>E0966</b>
<input type="checkbox"/> H-Harness with Padded Covers	<b>\$143</b>	<b>E0960</b>
<input type="checkbox"/> 3-Pt Positioning Belt with Depth Adjustable Crotch Strap	<b>\$49</b>	<b>E0978</b>
Foot Positioners - Select Size Below:		
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	<b>\$156</b>	<b>K0108</b>
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	<b>\$156</b>	<b>K0108</b>
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	<b>\$156</b>	<b>K0108</b>

### Transit Options

	Price	HCPCS
<input type="checkbox"/> Q'Straint Transit Lap Belt	<b>\$90</b>	<b>E0978</b>
<input type="checkbox"/> Sure-Lok Transit Lap Belt	<b>\$90</b>	<b>E0978</b>
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small ^	<b>\$151</b>	<b>E0950</b>
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium ^	<b>\$151</b>	<b>E0950</b>

^ Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

### Pelvic Positioning

	Price	HCPCS
<input type="checkbox"/> 3-Pt Positioning Belt with Depth Adjustable Crotch Strap	<b>\$49</b>	<b>E0978</b>
<input type="checkbox"/> Support - General Use Cushion *	<b>\$241</b>	<b>E2601</b>
<input type="checkbox"/> Position Cushion - Medial Thigh Support/Anti Thrust *	<b>\$466</b>	<b>E2607</b>
<input type="checkbox"/> Align Cushion - Lateral Pelvis/Lateral Thigh Support *	<b>\$466</b>	<b>E2607</b>
<input type="checkbox"/> Incontinence Cushion Liner (Cushion req)	<b>\$72</b>	<b>E2619</b>

\* Solid Seat Insert included with all Cushions.

### Foot & Leg Positioning

	Price	HCPCS
<input type="checkbox"/> Individual Angle Adjustable Footplates	<b>\$179</b>	<b>K0040</b>
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	<b>\$325</b>	<b>K0108</b>
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	<b>\$156</b>	<b>K0108</b>
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	<b>\$156</b>	<b>K0108</b>
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	<b>\$156</b>	<b>K0108</b>
<input type="checkbox"/> Legrest Extension - 2" *	<b>\$55</b>	
<input type="checkbox"/> Medial Thigh Support (Abductor)	<b>\$162</b>	<b>E0957</b>
<input type="checkbox"/> Lateral Thigh Support (Adductor)	<b>\$162</b>	<b>E0956</b>
<input type="checkbox"/> Calf Panel	<b>\$61</b>	<b>E0995</b>

\* 2" Legrest Extension interferes with caster swivel when seat is in upright position.  
 Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

RD10	RD12	RD14	RD16
903850	903850	903850	903850
903857	903857	903857	903857
903855	903855	903855	903855
903851	903851	903851	903851
903759	903759	903759	903759
903854	903854	903854	903854
904189	904189	904189	904189
904753	904753	904753	904753
904754	904754	904754	904754
904752	904752	904752	904752

RD10	RD12	RD14	RD16
904341	903167	903168	903169
903332	902863	900484	900138
903988	903988	903988	903988
904017	904017	904017	904017
903415	903415	903415	903415
903416	903416	903416	903416
903417	903417	903417	903417

RD10	RD12	RD14	RD16
903094	903094	903094	903094
903232	903232	903232	903232
903669	903669		
		903670	903670

RD10	RD12	RD14	RD16
904017	904017	904017	904017
904347	904350	904353	904356
904348	904351	904354	904357
904349	904352	904355	904358
902324	902324	902324	902324

RD10	RD12	RD14	RD16
	904957	904958	904959
904953	904954	904955	904956
903415	903415	903415	903415
903416	903416	903416	903416
903417	903417	903417	903417
904095	904095	904096	904097
904359	902667	902647	903140
904360	902668	902648	903141
904361	904361	904361	904361

- Trunk Positioning**
- H-Harness with Padded Covers (3-Pt Positioning Belt Req) **Price \$143 HCPCS E0960**
  - Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)
    - X Small= 2.75"x 3.75" \* **\$118 E0956**
    - Small = 3"x 5" **\$118 E0956**
    - Medium = 3.75"x 5.5" **\$118 E0956**
    - Large = 4.75"x 6" **\$118 E0956**
  - Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)
    - X Small= 2.75"x 3.75" \* **\$118 E0956**
    - Small = 3"x 5" **\$118 E0956**
    - Medium = 3.75"x 5.5" **\$118 E0956**
    - Large = 4.75"x 6" **\$118 E0956**
  - Right S/A Adj. Hardware for Lateral Trunk Support **\$238 E1028**
  - Left S/A Adj. Hardware for Lateral Trunk Support **\$238 E1028**
  - Full Torso Swing-Away Support Vest **\$238**
  - Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1") **\$353**
- \* RD10 requires X Small laterals.

RD10	RD12	RD14	RD16
903988	903988	903988	903988
904385	904385	904385	904385
	904071	904071	904071
	904072	904072	904072
	904073	904073	904073
904465	904465	904465	904465
	904466	904466	904466
	904467	904467	904467
	904468	904468	904468
904362	904362	904362	904362
904363	904363	904363	904363
904364	903261	903260	903259
904211	904211	904110	904110

- Head Supports**
- Padded Headwings **Price \$244 HCPCS E0955**
  - Occi Headwings **\$228 E0955**
  - Headrest Extension **\$130 E0966**

RD10	RD12	RD14	RD16
904365	902673	902653	903147
904779	904779	904779	904779
903332	902863	900484	900138

- Frame Option**
- Recline Option **Price \$461 HCPCS E1226**
  - Height Adjustable Flip-Up Armrests (pair) **\$228 E0973**
  - Solid Side Panels **\$211**
  - Mesh Side Panels **\$148**

RD10	RD12	RD14	RD16
904342	904343	904344	904345
904346	904346	904346	904346
904388	902672	902652	903145
904389	902680	902681	903146

- Tire Options**
- 12.5" Rear Solid Knobby Tire **Price \$199 HCPCS**
  - 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire **\$299**

RD10	RD12	RD14	RD16
904117	904117	904117	904117
904118	904118	904118	904118

- Accessories**
- Headrest Cover (Canopy) **Price \$228 HCPCS**
  - Upper Extremity Support Surface (Clear Tray) ^ **\$111 E0950**
  - Adj. Removable Hardware (Required for Clear Tray) **\$227 E1028**
  - Medical Necessity Storage Basket **\$295**
  - Under Seat Storage Basket - Tilt-in-Space **\$75**
  - Oxygen Tank Bag **\$179 E2208**
  - Utility Bag **\$38**
  - Attendant Hand Brakes (pair) **\$236**
  - Caster Lock for 2" Wide Front Tires **\$94 K0073**
  - Rear Anti-Tip Tubes (Pair)\* **\$115 E0971**
  - Rain Cover w/Out Windows - Small (Canopy Req) **\$150**
  - Rain Cover w/Out Windows - Large (Canopy Req) **\$150**
  - Mosquito Net (Canopy req) **\$150**
  - Extended Headrest Cover - No Windows (Canopy) **\$250**
  - Extended Headrest Cover - with Windows (Canopy) **\$295**
- ^ Flip-Up Armrests and Hardware required.  
\*Highly recommended for users weighing 50lbs or less.

RD10	RD12	RD14	RD16
904366	902679	902659	903152
904367	903405	903401	903402
903369	903369	903369	903369
904476	904476	904477	904477
904529	904529	904530	904530
902870	902870	902870	902870
904531	904531	904531	904531
		904370	904371
904650	904650	904650	904650
904649	904649	904649	904649
904939	904939	904939	
			904940
904655	904655	904655	904656
904728	904729	904730	904731
904744	904745	904746	904747

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo mobility base was tested with Rodeo seating module equipped with head support, H-Harness with padded Covers, 3-point restraint system, foot positioners, and 4 transit anchor brackets.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specially designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.ercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-888-CONVAID**

2830 California St., Torrance, CA 90503 • 1-888-CONVAID (266-8243), 310-618-0111 • Fax: 310-618-2166

Prices in USD & subject to change, ex works Torrance, California.