

Convoid Rodeo Tilt - HCPCS E1232/E1161

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

Rodeo Tilt Standard Features:

- | | |
|--|--|
| <ul style="list-style-type: none"> • 5 to 45 Degrees of Adjustable Tilt • Silver Powder Coated Frame • One Piece Height Adjustable Push Handle • One Piece Angle Adjustable Footplate • Seat Depth Growth of 5" • 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire • Quick Release Rear Wheels • Single Foot Activated Wheel Locks • Closure Strap • HCPCS code E1161 is only valid on size 16 (RD16) | <ul style="list-style-type: none"> • Made in the USA • Transit Anchor Brackets • H-Harness with Padded Covers • 3-Pt Positioning Belt w/ Depth Adjustable C • Foot Positioners • Headrest Extension • Cordura Upholstery • Tool-less Upholstery Removal • Zippered Pockets in Seat Bottom Upholstery • Elevating Legrest |
|--|--|

Standard Models:

- | | Price |
|---|----------------|
| <input type="checkbox"/> Rodeo Transit RD10T - 904860 | \$3,226 |
| <input type="checkbox"/> Rodeo Transit RD12T - 904861 | \$3,226 |
| <input type="checkbox"/> Rodeo Transit RD14T - 904862 | \$3,226 |
| <input type="checkbox"/> Rodeo Transit RD16T - 904863 | \$3,226 |

Upholstery Colors

- | | Price |
|---|-------|
| <input type="checkbox"/> Royal Blue | \$0 |
| <input type="checkbox"/> Forest Green | \$0 |
| <input type="checkbox"/> Cherry Red | \$0 |
| <input type="checkbox"/> Sassy Purple | \$0 |
| <input type="checkbox"/> Princess Pink | \$0 |
| <input type="checkbox"/> Panther Black | \$0 |
| <input type="checkbox"/> Chocolate Brown | \$0 |
| <input type="checkbox"/> Aqua Blue | \$0 |
| <input type="checkbox"/> Apple Green | \$0 |
| <input type="checkbox"/> Tangerine Orange | \$0 |

Transit Required Items

- | | Price | HCPCS |
|---|-------|-------|
| <input type="checkbox"/> Transit Option (Manufacturer Installation Required) | \$0 | |
| <input type="checkbox"/> Headrest Extension | \$0 | E0966 |
| <input type="checkbox"/> H-Harness with Padded Covers | \$0 | E0960 |
| <input type="checkbox"/> 3-Pt Positioning Belt with Depth Adjustable Crotch Strap | \$0 | E0978 |
| Foot Positioners - Select Size Below: | | |
| <input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair) | \$0 | K0108 |
| <input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair) | \$0 | K0108 |
| <input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair) | \$0 | K0108 |

Transit Options

- | | Price | HCPCS |
|---|-------|-------|
| <input type="checkbox"/> Q'Straint Transit Lap Belt | \$90 | E0978 |
| <input type="checkbox"/> Sure-Lok Transit Lap Belt | \$90 | E0978 |
| <input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small ^ | \$151 | E0950 |
| <input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium ^ | \$151 | E0950 |
- ^ Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

Pelvic Positioning

- | | Price | HCPCS |
|---|-------|-------|
| <input type="checkbox"/> Support - General Use Cushion * | \$241 | E2601 |
| <input type="checkbox"/> Position Cushion - Medial Thigh Support/Anti Thrust * | \$466 | E2607 |
| <input type="checkbox"/> Align Cushion - Lateral Pelvis/Lateral Thigh Support * | \$466 | E2607 |
| <input type="checkbox"/> Incontinence Cushion Liner (Cushion req) | \$72 | E2619 |
- * Solid Seat Insert included with all Cushions.

Foot & Leg Positioning

- | | Price | HCPCS |
|---|-------|-------|
| <input type="checkbox"/> Individual Angle Adjustable Footplates | \$179 | K0040 |
| <input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required) | \$325 | K0108 |
| <input type="checkbox"/> Legrest Extension - 2" * | \$55 | |
| <input type="checkbox"/> Medial Thigh Support (Abductor) | \$162 | E0957 |
| <input type="checkbox"/> Lateral Thigh Support (Adductor) | \$162 | E0956 |
| <input type="checkbox"/> Calf Panel | \$61 | E0995 |

* 2" Legrest Extension interferes with caster swivel when seat is in upright position.
 Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

RD10	RD12	RD14	RD16
903850	903850	903850	903850
903857	903857	903857	903857
903855	903855	903855	903855
903851	903851	903851	903851
903759	903759	903759	903759
903854	903854	903854	903854
904189	904189	904189	904189
904753	904753	904753	904753
904754	904754	904754	904754
904752	904752	904752	904752

RD10	RD12	RD14	RD16
904341B	903167B	903168B	903169B
903332B	902863B	900484B	900138B
903988B	903988B	903988B	903988B
904017B	904017B	904017B	904017B
903415B	903415B	903415B	903415B
903416B	903416B	903416B	903416B
903417B	903417B	903417B	903417B

RD10	RD12	RD14	RD16
903094	903094	903094	903094
903232	903232	903232	903232
903669	903669		
		903670	903670

RD10	RD12	RD14	RD16
904347	904350	904353	904356
904348	904351	904354	904357
904349	904352	904355	904358
902324	902324	902324	902324

RD10	RD12	RD14	RD16
	904957	904958	904959
904953	904954	904955	904956
904095	904095	904096	904097
904359	902667	902647	903140
904360	902668	902648	903141
904361	904361	904361	904361

- Trunk Positioning**
- Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)
- X Small= 2.75"x 3.75" * **\$118**
 - Small = 3"x 5" **\$118**
 - Medium = 3.75"x 5.5" **\$118**
 - Large = 4.75"x 6" **\$118**
- Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)
- X Small= 2.75"x 3.75" * **\$118**
 - Small = 3"x 5" **\$118**
 - Medium = 3.75"x 5.5" **\$118**
 - Large = 4.75"x 6" **\$118**
- Right S/A Adj. Hardware for Lateral Trunk Support **\$238**
 - Left S/A Adj. Hardware for Lateral Trunk Support **\$238**
 - Full Torso Swing-Away Support Vest **\$238**
 - Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1") **\$353**
- * RD10 requires X Small laterals.

Price	HCPCS
\$118	E0956
\$118	E0956
\$118	E0956
\$118	E0956
\$118	E0956
\$118	E0956
\$118	E0956
\$118	E0956
\$238	E1028
\$238	E1028
\$238	E1028
\$353	E1028

RD10	RD12	RD14	RD16
904385	904385	904385	904385
	904071	904071	904071
	904072	904072	904072
	904073	904073	904073
904465	904465	904465	904465
	904466	904466	904466
	904467	904467	904467
	904468	904468	904468
904362	904362	904362	904362
904363	904363	904363	904363
904364	903261	903260	903259
904211	904211	904110	904110

- Head Supports**
- Padded Headwings **\$244**
 - Occi Headwings **\$228**
- Frame Option**
- Recline Option **\$461**
 - Height Adjustable Flip-Up Armrests (pair) **\$228**
 - Solid Side Panels **\$211**
 - Mesh Side Panels **\$148**
- Tire Options**
- 12.5" Rear Solid Knobby Tire **\$199**
 - 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire **\$299**

Price	HCPCS
\$244	E0955
\$228	E0955
\$461	E1226
\$228	E0973
\$211	E0973
\$148	E0973
\$199	E0955
\$299	E0955

RD10	RD12	RD14	RD16
904365	902673	902653	903147
904779	904779	904779	904779
904342	904343	904344	904345
904346	904346	904346	904346
904388	902672	902652	903145
904389	902680	902681	903146
904117	904117	904117	904117
904118	904118	904118	904118

- Accessories**
- Headrest Cover (Canopy) **\$228**
 - Upper Extremity Support Surface (Clear Tray) ^ **\$111**
 - Adj. Removable Hardware (Required for Clear Tray) **\$227**
 - Medical Necessity Storage Basket **\$295**
 - Under Seat Storage Basket - Tilt-in-Space **\$75**
 - Oxygen Tank Bag **\$179**
 - Utility Bag **\$38**
 - Attendant Hand Brakes (pair) **\$236**
 - Caster Lock for 2" Wide Front Tires **\$94**
 - Rear Anti-Tip Tubes (Pair)* **\$115**
 - Rain Cover w/Out Windows - Small (Canopy Req) **\$150**
 - Rain Cover w/Out Windows - Large (Canopy Req) **\$150**
 - Mosquito Net (Canopy req) **\$150**
 - Extended Headrest Cover - No Windows (Canopy) **\$250**
 - Extended Headrest Cover - with Windows (Canopy) **\$295**
- ^ Flip-Up Armrests and Hardware required.
*Highly recommended for users weighing 50lbs or less.

Price	HCPCS
\$228	E0955
\$111	E0950
\$227	E1028
\$295	E1028
\$75	E1028
\$179	E2208
\$38	E2208
\$236	E2208
\$94	K0073
\$115	E0971
\$150	E0971
\$150	E0971
\$150	E0971
\$250	E0971
\$295	E0971

RD10	RD12	RD14	RD16
904366	902679	902659	903152
904367	903405	903401	903402
903369	903369	903369	903369
904476	904476	904477	904477
904529	904529	904530	904530
902870	902870	902870	902870
904531	904531	904531	904531
		904370	904371
904650	904650	904650	904650
904649	904649	904649	904649
904939	904939	904939	
			904940
904655	904655	904655	904656
904728	904729	904730	904731
904744	904745	904746	904747

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo mobility base was tested with Rodeo seating module equipped with head support, H-Harness with padded Covers, 3-point restraint system, foot positioners, and 4 transit anchor brackets.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specially designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.nercwtts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

For expedited shipping rates, please call 1-888-CONVAID
 2830 California St., Torrance, CA 90503 • 1-888-CONVAID (266-8243), 310-618-0111 • Fax: 310-618-2166
 Prices in USD & subject to change, ex works Torrance, California.