

Choose the frame and seat combination to find the max. user load



Product		Panda Futura Active Size 2	Panda Futura Active Size 2½	Panda Futura Active Size 3	Panda Futura Active Size 4	Panda Futura Active Size 5	Panda Futura Standard Size 1	Panda Futura Standard Size 2	Panda Futura Standard Size 2½	Panda Futura Standard Size 3	Panda Futura Standard Size 4	Panda Futura Standard Size 5	x:panda Size 1S	x:panda Size 1M	x:panda Size 1L	x:panda Size 2S	x:panda Size 2M	x:panda Size 2L	x:panda Size 3S	x:panda Size 3M	x:panda Size 3L	x:panda Size 4S	x:panda Size 4M	x:panda Size 4L	Serval seat Size 1	Serval seat Size 2	Nandu seat Size 1	Nandu seat Size 2	Nandu seat Size 3	Nandu seat Size 4	Stingray seat Size 1	Stingray seat Size 2
High-Low Indoor Size 1	Ordinary use	50	50	50	50	-	50	50	50	50	50	-	45*	45*	45*	50*	50*	50*	-	-	-	-	-	-	35	43	50	50	50	50	-	-
	Transportation in vehicles	50	50	50	50	-	50	50	50	50	50	-	45*	45*	45*	50*	50*	50*	-	-	-	-	-	-	35	43	-	-	-	-	-	-
High-Low Outdoor Size 1	Ordinary use	-	-	-	-	-	27	27	27	43	43	-	43*	43*	43*	43*	43*	43*	-	-	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	27	43	43	-	43*	43*	43*	43*	43*	43*	-	-	-	-	-	-	-	-	-	-	-	-	-	
High-Low:x Indoor Size 1	Ordinary use	50*	50*	50*	50*	-	50*	50*	50*	50*	50*	-	45	45	45	50	50	50	-	-	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	50*	50*	50*	50*	43	50*	50*	50*	50*	50*	43	-	-	-	-	-	-	47	47	47	-	-	-	-	-	-	-	-	-	-	
High-Low:x Indoor Size 2	Ordinary use	-	-	-	50*	85	-	-	-	-	50*	83	-	-	-	-	-	-	60	60	60	80	80	80	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27*	27*	27*	43*	43*	-	43	43	43	43	43	43	-	-	-	-	-	-	-	-	-	-	-	-	-	-
High-Low:x Outdoor Size 1	Ordinary use	50*	50*	50*	50*	-	50*	50*	50*	50*	50*	-	45	45	45	50	50	50	-	-	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27*	27*	27*	43*	43*	-	43	43	43	43	43	43	-	-	-	-	-	-	-	-	-	-	-	-	-	
High-Low:x Outdoor Size 2	Ordinary use	50*	50*	50*	50*	43	50*	50*	50*	50*	50*	43	-	-	-	-	-	-	60	60	60	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	43	27*	27*	27*	43*	43*	43	-	-	-	-	-	-	43	43	43	-	-	-	-	-	-	-	-	-	-	
High-Low:v Indoor Size 1	Ordinary use	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	40	
	Transportation in vehicles	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	40	40
Nandu Frame Size 1	Ordinary use	50*	50*	-	-	-	50*	50*	50*	-	-	-	45*	45*	45*	50*	50*	50*	-	-	-	-	-	-	-	-	50	50	-	-	-	
	Transportation in vehicles	-	-	50*	50*	-	-	-	-	50*	50*	-	-	-	-	50*	50*	50*	50*	50*	50*	-	-	-	-	-	-	-	50	50	-	-
Combi Frame:x Size 1 20" 12½"	Ordinary use	50	50	50	-	-	50	50	50	50	-	-	45	45	45	50	50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	43	43	-	-	45	45	45	50	50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
Combi Frame:x Size 1 22" 24"	Ordinary use	50	50	50	-	-	50	50	50	50	-	-	45	45	-	50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	43	43	-	-	45	45	-	50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
Combi Frame:x Size 2 20" 22" 24" 12½"	Ordinary use	50	50	50	50	58	50	50	50	50	50	-	45	45	45	50	50	50	60	-	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	43	43	43	-	45	45	45	50	50	50	57	-	-	-	-	-	-	-	-	-	-	-	-	
Combi Frame:x Size 3 20" 12½"	Ordinary use	50	50	50	50	58	50	50	50	50	60	-	45	45	45	50	50	50	60	60	60	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	43	43	43	-	45	45	45	50	50	50	57	57	57	-	-	-	-	-	-	-	-	-	-	
Combi Frame:x Size 3 22" 24"	Ordinary use	50	50	50	50	58	50	50	50	50	60	-	45	45	45	50	50	50	60	60	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	43	43	43	-	45	45	45	50	50	50	57	57	-	-	-	-	-	-	-	-	-	-	-	
Combi Frame:x Size 4 20" 22" 24" 12½"	Ordinary use	50	50	50	50	83	50	50	50	50	60	85	45	45	45	50	50	50	60	60	60	80	80	80	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	43	43	43	85	45	45	45	50	50	50	57	57	57	57	57	57	-	-	-	-	-	-	-	
Multi Frame Size 1 - 12½"	Ordinary use	50	50	50	-	-	50	50	50	50	-	-	45	45	45	50	50	50	60	60	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	27	50	-	-	45	45	45	50	50	50	57	57	-	-	-	-	-	-	-	-	-	-	-	
Multi Frame Size 1 - 22"	Ordinary use	50	50	50	-	-	50	50	50	50	-	-	45	45	-	50	50	50	-	-	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	27	50	-	-	45	45	-	50	50	50	-	-	-	-	-	-	-	-	-	-	-	-	-	
Multi Frame Size 2 - 12½"	Ordinary use	50	50	50	50	-	50	50	50	50	60	-	45	45	45	50	50	50	60	60	60	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	27	50	57	-	45	45	45	50	50	50	57	57	57	-	-	-	-	-	-	-	-	-	-	
Multi Frame Size 2 - 22"	Ordinary use	50	50	50	50	-	50	50	50	50	60	-	45	45	-	50	50	-	60	-	-	-	-	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	27	50	57	-	45	45	-	50	50	-	57	-	-	-	-	-	-	-	-	-	-	-	-	
Multi Frame Size 3 - 12½"	Ordinary use	50	50	50	50	83	50	50	50	50	60	85	45	45	45	50	50	50	60	60	60	80	80	-	-	-	-	-	-	-	-	
	Transportation in vehicles	-	-	-	-	-	27	27	27	50	57	85	45	45	45	50	50	50	57	57	57	57	57	-	-	-	-	-	-	-	-	
Multi Frame Size 3 - 22"	Ordinary use	50	50	50	50	83	50	50	50	50	60	85	45	45	45	-	-	-	60	60	-	-	-	-	-	-	-	-	-	-		
	Transportation in vehicles	-	-	-	-	-	27	27	27	50	57	85	45	45	45	-	-	-	57	57	-	-	-	-	-	-	-	-	-	-	-	
Multi Frame:x Size 3 - 12½"	Ordinary use	-	-	-	-	-	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	80	80	80	-	-	-	-	-	-		
	Transportation in vehicles	-	-	-	-	-	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	57	57	57	-	-	-	-	-	-		
Multi Frame:x Size 3 - 22"	Ordinary use	-	-	-	-	-	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
	Transportation in vehicles	-	-	-	-	-	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
Serval Frame Size 1	Ordinary use	35	35	-	-	-	35	35	35	-	-	-	35	35	-	35	-	-	-	-	-	-	-	-	35	-	-	-	-	-		
	Transportation in vehicles	-	-	-	-	-	27	27	27	-	-	-	35	35	-	35	-	-	-	-	-	-	-	-	35	-	-	-	-	-		
Serval Frame Size 2	Ordinary use	43	43	-	-	-	43	43	43	-	-	-	45	45	45	50	50	-	-	-	-	-	-	-	50	-	-	-	-	-		
	Transportation in vehicles	-	-	-	-	-	27	27	27	-	-	-	43	43	43	43	43	-	-	-	-	-	-	-	43	-	-	-	-	-		
Stingray Frame	Ordinary use	-	-	-	-	-	40	40	40	-	-	-	40	40	40	40	40	-	-	-	-	-	-	-	-	-	-	-	-	40	40	
	Transportation in vehicles	-	-	-	-	-	40	40	40	-	-	-	40	40	40	40	40	-	-	-	-	-	-	-	-	-	-	-	-	40	40	
Strong Base	Ordinary use	-	-	-	-	-	-	-	-	-	-	85	-	-	-	-	-	-	-	-	-	80	80	80	-	-	-	-	-	-		

* Special interface needed