

Convoid Metro - HCPCS E1037

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

Metro Standard Features:

- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (ME12/14/16)
- Quick release Rolko rear wheels (ME12/14/16)
- 8" Front, 11" Rear Solid Tire (ME18)
- Wheel Locks
- Adjustable Swing-Away Footplates
- 30° Fixed Tilt
- Silver Powder Coated Frame
- 2-Pt Positioning Belt
- Made in the USA

Standard Models:

	Price
<input type="checkbox"/> Metro ME12 - 901402	\$1,050
<input type="checkbox"/> Metro ME14 - 901400	\$1,050
<input type="checkbox"/> Metro ME16 - 901335	\$1,050
<input type="checkbox"/> Metro ME18 - 903127	\$1,050

<u>Upholstery Colors</u>	<u>Price</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> Cherry Red	\$0	903855	903855	903855	903855
<input type="checkbox"/> Royal Blue	\$0	903850	903850	903850	903850

<u>Transit Required Items</u>	<u>Price</u>	<u>HCPCS</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	\$200		901403	901401	901336	903128
<input type="checkbox"/> Headrest Extension	\$130	E0966	903263	903264	903265	903266
<input type="checkbox"/> 3-Pt Positioning Belt	\$49	E0978	904017	904017	904017	904017
Foot Positioners - Select Size Below:			903415	903415	903415	903415
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$156		903416	903416	903416	903416
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$156		903417	903417	903417	903417
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$156					

<u>Transit Options</u>	<u>Price</u>	<u>HCPCS</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$90	E0978	903094	903094	903094	903094
<input type="checkbox"/> Sure-Lok Transit Lap Belt	\$90	E0978	903232	903232	903232	903232

<u>Tire Options</u>	<u>Price</u>	<u>HCPCS</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	\$199		904117	904117	904117	
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$299		904118	904118	904118	
<input type="checkbox"/> 8"x2" Front Ribbed, 12.5" Rear Solid Knobby Tire	\$299					903084
<input type="checkbox"/> 8"x2" Front Ribbed, 12.5" Rear Pneumatic Knobby Tire	\$299					903083

<u>Accessories</u>	<u>Price</u>	<u>HCPCS</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> 3-Pt Positioning Belt	\$49	E0978	904017	904017	904017	904017
<input type="checkbox"/> Wheelchair Travel Bag	\$206		902022	902022	902022	902022
<input type="checkbox"/> Under Seat Storage Basket - Fixed Tilt	\$68		904046	904047	904047	904048

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

The Metro Transit option includes four manufacturer installed wheelchair transportation anchors. The Metro transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are ME12: 66 lbs./30kg; ME14: 100 lbs./45.5kg/ME16 and ME18: 170 lbs./77 kg.

The Convoid Metro was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is being, or can be, made.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166

Prices in USD & subject to change, ex works Torrance, California.

Effective: 2/21/2018