

# Convoid Carrot 3 Booster Seat - Child Restraint System

## Canadian Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

**Carrot 3 Booster Seat Standard Features:**

- |  |  |
|--|--|
| <ul style="list-style-type: none"> <li>• Carrot 3 Backrest</li> <li>• Recline Adjustment (90° - 180°)</li> <li>• Head Support: 8.5" x 10.5" (H x L)</li> <li>• 4 Point Positioning Harness</li> <li>• H-Harness Attaching Bracket</li> </ul> | <ul style="list-style-type: none"> <li>• Booster Seat Base</li> <li>• Lower Anchors Connector System</li> <li>• 2" Backrest Extension</li> <li>• Shoulder Protector Wings</li> </ul> |
|--|--|

**Standard Models:**

	<u>Price</u>
<input type="checkbox"/> Carrot 3 Booster Seat CRT3BLACK - CRT3BS-10	<b>\$3,601</b>
<input type="checkbox"/> Carrot 3 Booster Seat CRT3BLUE - CRT3BS-20	<b>\$3,601</b>
<input type="checkbox"/> Carrot 3 Booster Seat CRT3RED - CRT3BS-30	<b>\$3,601</b>

All prices in CAD

**Carrot Booster Seat Recommended Item**

	<u>Price</u>
<input type="checkbox"/> Tether strap (Black)	<b>\$117</b>
<input type="checkbox"/> Top Securement Strap *	<b>\$70</b>
<input type="checkbox"/> Lower Securement Strap *	<b>\$85</b>

Black	Blue	Red
CT3TS0-10	CT3TS0-10	CT3TS0-10
CT3HSF-10	CT3HSF-10	CT3HSF-10
CT3LFS-10	CT3LFS-10	CT3LFS-10

*\*For vehicles not equipped with anchor points.*

**Accessories**

	<u>Price</u>
<input type="checkbox"/> Booster Support Tray (Crash Tested)	<b>\$380</b>
<input type="checkbox"/> 3D Body Balance Pads L/R Set	<b>\$197</b>
<input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°)	<b>\$148</b>
<input type="checkbox"/> Backrest Extension - 4" (Extends back height to 32.8")	<b>\$117</b>

Black	Blue	Red
CRT3BT-10	CRT3BT-20	CRT3BT-30
CT3BBP-00	CT3BBP-00	CT3BBP-00
CT3SA0-10	CT3SA0-20	CT3SA0-30
CT3BE4-10	CT3BE4-20	CT3BE4-30

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245**

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166


Prices in CAD & subject to change, ex works Torrance, California.

www.convoid.com

Effective: 10/3/2018

Rev 6

**CLIENT MEASUREMENTS:**



**Carrot 3 Base Model Specifications:**

Seat Width:	16.5"
Seat Depth:	15" - 16.5" *
Backrest Height Adjustment:	28" - 33"
Head Support Height:	8.5"
Backrest Height:	29" - 34"
Backrest Extension Height (Standard):	2"
Backrest Extension Height (Option):	4"
Shoulder Wing Height (Standard):	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"
Product Weight:	21 lbs.
Recline Adjustment:	90° - 180°
Product Overall Width:	20" - 22"

**User Measurement Guide:**

Minimum User Weight:	79 lbs.
Maximum User Weight:	165 lbs.
Minimum User Height:	57"
Maximum User Height:	70"

\*Factory Setting

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166