

# R82 Pony - E8002

## CAD Retail Price List/ Order Form




Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

|                 |                 |
|-----------------|-----------------|
| Bill To:        | Ship To:        |
| Address:        | Address:        |
| City/State/Zip: |                 |
| Phone:          | City/State/Zip: |
| Email           | Phone:          |

**US Configuration** **Effective: 10.1.18 - Rev 1**

- Pony
- Head Support
- Abdominal Pad
- Hip Support
- Arm Support



- A** (Total Height) \_\_\_\_\_
- B** (Shoulder Height) \_\_\_\_\_
- C** (Chest Height) \_\_\_\_\_
- D** (Buttock Height) \_\_\_\_\_
- E** (Chest Width) \_\_\_\_\_
- F** (Shoulder Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_

**Measurements**



|                 | <b>Size 0</b><br><b>Inches</b> | <b>Size 1</b><br><b>Inches</b> | <b>Size 2</b><br><b>Inches</b> |
|-----------------|--------------------------------|--------------------------------|--------------------------------|
| Seat height (L) | 9¼ - 14¾"                      | 14¾ - 19½"                     | 18¾ - 25"                      |
| Width (G)       | 23"                            | 27"                            | 30"                            |
| Length (H)      | 26¼"                           | 28"                            | 40½"                           |
| Height (I)      | 20¼ - 27"                      | 23½ - 32"                      | 26½ - 40½"                     |
| Angle range (A) | 32°                            | 32°                            | 32°                            |
| Weight          | 15½ lbs.                       | 22 lbs.                        | 30 lbs.                        |
| Max. load       | 66 lbs.                        | 88 lbs.                        | 110 lbs.                       |

**Pony Standard Features**



- Rear Drag Wheels with Directional Stabilizer
- Leg & Foot Separator
- Height Adjustable Seat with Spring
- Height Adjustable Removable Back Support
- Height Adjustable Chest Support
- 5 Wheels (1 Rear & 2 Front with Locks)
- Height & Depth Adjustable Handle Bar

**NOTES:**

| Pony | Description                                       | Size | Part No. | CAD Price |
|------|---|------|----------|-----------|
|      | <input type="checkbox"/> Pony Walking Aid, Yellow | 0    | 86400    | \$2,899   |
|      | <input type="checkbox"/> Pony Walking Aid, Purple | 1    | 86401    | \$2,899   |
|      | <input type="checkbox"/> Pony Walking Aid, Lilac  | 2    | 86402    | \$2,899   |

| Accessories | Description | Measurements | Size | Part No. | CAD Price |
|-------------|-------------|--------------|------|----------|-----------|
|-------------|-------------|--------------|------|----------|-----------|



- Hip & Thigh Supports (Pr.) 6" x 3" 0 86408 \$412
  - Hip & Thigh Supports (Pr.) 6" x 4" 1 86410 \$412
  - Hip & Thigh Supports (Pr.) •• 6" x 4" 2 86412 \$412
- Custom Order; Approx. lead time: 6-8 weeks

| Minimum/Maximum Width | Size 0   | Size 1   | Size 3    |
|-----------------------|----------|----------|-----------|
|                       | 6" x 8¾" | 6" x 8¾" | 6" x 10¾" |

| Description | Measurements | Size | Part No. | CAD Price |
|-------------|--------------|------|----------|-----------|
|-------------|--------------|------|----------|-----------|



- Abdominal Cushion (Ea.) •• 1 86420 \$243
  - Abdominal Cushion (Ea.) •• 2 86422 \$243
- Custom Order; Approx. lead time: 6-8 weeks




- Armrest/Tray (Ea.) 0 86427 \$605
  - Armrest/Tray (Ea.) 1 86428 \$605
  - Armrest/Tray (Ea.) •• 2 86429 \$605
- Custom Order; Approx. lead time: 6-8 weeks



- Vertical Handle for Mounting in the Middle (Pr.) One Size 86432 \$144

| Description | Measurements | Size | Part No. | CAD Price |
|-------------|--------------|------|----------|-----------|
|-------------|--------------|------|----------|-----------|


- Head Support No. 1, Black PU One Size 6171550 \$256
- Head Support No. 2, Black PU One Size 6172550 \$256
- Head Support No. 3, Black PU One Size 6173550 \$256



- Anatomic Head Support ♦ A:5½ x B:5½ x C:9½" Small 5099001 \$318
  - Head Support Cover (Required) Small 9923-83 N/C
  - Anatomic Head Support ♦ A:6 x B:5½ x C:10" Medium 5099002 \$318
  - Head Support Cover (Required) Medium 9945-83 N/C
  - Anatomic Head Support ♦ A:6½ x B:5½ x C:12½" Large 5099003 \$318
  - Head Support Cover (Required) Large 9967-83 N/C
- ♦ Turtle Bar System Short or Long Required; 98619-25 or 98619-26

- Headband System for Anatomic Head Support Only One Size 99500 \$287

| Description | Measurements | Size | Part No. | CAD Price |
|-------------|--------------|------|----------|-----------|
|-------------|--------------|------|----------|-----------|



- Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable**
- Turtle Bar System ♦ A: 5", B: 4-8½" Short 98619-25 \$287
  - Turtle Bar System ♦ A: 10", B: 4-8½" Long 98619-26 \$287
- ♦ Required with Anatomic Head Support; 5099001, 5099002 or 5099003

- Connecting Rod for Head Support 7¾" 1 861424-3 \$322
- Connecting Rod for Head Support 13" 2 861424-4 \$322

| Accessories | Description | Measurements | Size | Part No. | CAD Price |
|-------------|-------------|--------------|------|----------|-----------|
|-------------|-------------|--------------|------|----------|-----------|










- Leg Protection One Size 86450 \$274

| Description | Measurements | Size | Part No. | CAD Price |
|-------------|--------------|------|----------|-----------|
|-------------|--------------|------|----------|-----------|



- Push Handle One Size 86416 \$318

| Spare Parts   | Description   | Measurements        | Size        | Part No.        | CAD Price        | Qty        |
|---|---|---------------------|-------------|-----------------|------------------|------------|
|    | <input type="checkbox"/> Chest Support, Complete                |                     | 0           | 86310           | \$636            | _____      |
|   | <input type="checkbox"/> Chest Support, Complete                |                     | 1           | 86311           | \$636            | _____      |
|   | <input type="checkbox"/> Chest Support, Complete                |                     | 2           | 86312           | \$636            | _____      |
|   | <input type="checkbox"/> Straps for Chest Support               |                     | 0           | 5540250         | \$45             | _____      |
|   | <input type="checkbox"/> Straps for Chest Support               |                     | 1           | 5540300         | \$57             | _____      |
|   | <input type="checkbox"/> Straps for Chest Support               |                     | 2           | 5540350         | \$70             | _____      |
| <hr/>   |   |                     |             |                 |                  |            |
| <b>Description</b>  |   | <b>Measurements</b> | <b>Size</b> | <b>Part No.</b> | <b>CAD Price</b> | <b>Qty</b> |
|    | <input type="checkbox"/> Handle                                 |                     | 0           | 86320           | \$144            | _____      |
|   | <input type="checkbox"/> Handle                                 |                     | 1           | 86321           | \$144            | _____      |
|   | <input type="checkbox"/> Handle                                 |                     | 2           | 86322           | \$144            | _____      |
| <hr/>   |   |                     |             |                 |                  |            |
| <b>Description</b>  |   | <b>Measurements</b> | <b>Size</b> | <b>Part No.</b> | <b>CAD Price</b> | <b>Qty</b> |
|    | <input type="checkbox"/> Seat                                   |                     | 0           | 86330           | \$312            | _____      |
|   | <input type="checkbox"/> Seat                                   |                     | 1           | 86331           | \$312            | _____      |
|   | <input type="checkbox"/> Seat                                   |                     | 2           | 86332           | \$312            | _____      |
|    | <input type="checkbox"/> Back Support                           |                     | 0           | 86434           | \$256            | _____      |
|   | <input type="checkbox"/> Back Support                           |                     | 1, 2        | 86435           | \$256            | _____      |
| <hr/>   |   |                     |             |                 |                  |            |
| <b>Description</b>  |   | <b>Measurements</b> | <b>Size</b> | <b>Part No.</b> | <b>CAD Price</b> | <b>Qty</b> |
|    | <input type="checkbox"/> Directional Stabilizer for Rear Wheel  |                     | 0           | 86438           | \$109            | _____      |
|   | <input type="checkbox"/> Directional Stabilizer for Rear Wheel  |                     | 1           | 86440           | \$109            | _____      |
|   | <input type="checkbox"/> Directional Stabilizer for Rear Wheel  |                     | 2           | 86442           | \$109            | _____      |
|  | <input type="checkbox"/> Directional Stabilizer for Front Wheel |                     | 0, 1, 2     | 86448           | \$461            | _____      |
|  | <input type="checkbox"/> Drag Wheel                             |                     | One Size    | 86418           | \$102            | _____      |

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

**For shipping rates, please call 1-844-876-6245**

12801 E Independence Blvd, Matthews, NC 28105 • 1-844-876-6245 • Web: www.r82.com

Prices in CAD & is subject to change without notice

**For Quotes and Orders please send completed Order Form to: us.info@r82.com**