

R82 Stingray

CAD Retail Price List/ Order Form



Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

US Configuration


Effective: 2/01/2019 - Rev 7

- Frame (Stingray, Trekker, High-Low:v)
- Stingray Seat
- Cushions
- Head Support
- Side Support
- Hip Support
- Pommel
- Vests and Belts
- Tray
- Transport Fittings



- A (Back Height) _____
- B (Arm Height) _____
- C (Lower Leg Length) _____
- D (Seat Depth) _____
- E (Chest Width) _____
- F (Seat Width) _____
- G (Height: Seat to Top of Head) _____
- H (Shoulder Width) _____
- Age _____
- Weight _____

Measurements



	Size 1 Inches	Size 2 Inches
Effective Seat Width (A)	7-11 ³ / ₄ "	9-13 ³ / ₄ "
Seat Width Between Knee Supports	11"	12 ³ / ₄ "
Seat Width Between Hip Supports (Accessories)	7", 8 ¹ / ₂ ", 10 ¹ / ₄ "	9", 10 ¹ / ₂ ", 12"
Effective Seat Depth (B)	7 - 11 ³ / ₄ "	9 ³ / ₄ - 15 ¹ / ₂ "
Back Height at Shoulders (C)	10 ¹ / ₂ - 16"	14 - 19 ³ / ₄ "
Backrest Height (C1)	16 - 21 ¹ / ₂ "	19 ¹ / ₂ - 25 ¹ / ₄ "
Footrest to Seat Distance (D)	5 ¹ / ₂ - 17 ¹ / ₂ "	5 ¹ / ₂ - 17 ¹ / ₂ "
Seat Surface Height at Front Edge (E)	20 ¹ / ₄ "	20 ¹ / ₄ "
Back Width, Inside (F)	11 ³ / ₄ "	14"
Back Width, Outside (F1)	17 ¹ / ₂ "	20 ¹ / ₄ "
Overall Width (G)	23"	23"
Length with Turnable Front Wheels (H1)	28 ³ / ₄ "	28 ³ / ₄ "
Length with Large Front Wheels	32 ³ / ₄ "	32 ³ / ₄ "
Overall Height, Push Brace 15° (I)	39", 40", 41 ¹ / ₄ "	39", 40", 41 ¹ / ₄ "
Foot Plate, Width x Depth	12" x 7 ³ / ₄ "	12" x 7 ³ / ₄ "
Folded Length	35"	35"
Folded Height	19"	19"
Total Mass, Frame	29 ¹ / ₂ lbs.	29 ¹ / ₂ lbs.
Total Mass, Seat	12 lbs.	14 ¹ / ₂ lbs.
Max Load	88 lbs.	88 lbs.
Back Recline / Seat Plane Angle	Back: 170° / Seat: 45°	Back: 170° / Seat: 45°
Approx. Minimum User Height	29 ¹ / ₂ "	35 ¹ / ₂ "
Approx. Maximum User Height	50"	57 ¹ / ₂ "


NOTES:

Stingray Standard Features

- Tilt-in-Space - 45°
- Back Recline - 170°
- Turnable Seat - 180°
- One-Size Frame
- Accepts Two Seat Sizes for Maximum Growth Potential
- Adjustable Seat Depth
- Adjustable Back Height
- Suspension



Description	Part No.	CAD Price
<input type="checkbox"/> Solid Tires, Individual Wheel Locks	8724212	\$3,049
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Hand Brakes ••	8724211	\$3,517
<input type="checkbox"/> Solid Tires, Central Wheel Locks	8724222	\$3,174
<input type="checkbox"/> Solid Tires, Central Wheel Locks & Hand Brakes	8724221	\$3,641
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Quick Release Rear Axels ••Custom order; Approx. lead time: 6-8 weeks. Non-returnable item	8722212	\$3,361
<input type="checkbox"/> Transit Option* *Required if using O2 Tank Holder	One Size 8720902	\$249
<input type="checkbox"/> Medical Essential Carrier	One Size 8720600	\$121
<input type="checkbox"/> Vent Tray ♦* 13" x 17" ♦Requires Central Wheel Locks *44 lbs. weight capacity	One Size 8776290	\$574
<input type="checkbox"/> O2 Tank Holder * •• *Transit Option (8720902) is required ••Custom order; Approx. lead time: 6-8 weeks. Non-returnable item	8776295	\$350
<input type="checkbox"/> IV Pole (height adj. & collapsible) (TR) Anti-Tipper Required	One Size C904566	\$250
<input type="checkbox"/> IV Pole Holder** **Required if using IV Pole	8776291	\$175



Trekker Frame Standard Features

- One Piece Height Adjustable Push Handle (Inside Mount)
- Tilt-in-Space (-10° to +45°)
- 180° Reversible Seating Configuration
- Single Foot Activated Wheel Lock Bar
- 7½" x 2" Front and 11½" x 2½" Rear Solid Tires
- Quick Release Rear Wheels
- Made in the USA
- Transit Option - Transit Ready Frame Only





Description	Part No.	CAD Price
<input type="checkbox"/> Trekker Frame - Transit Ready (TR)^^	C905362	\$3,236
<input type="checkbox"/> Stingray Seat Module - Interface Adaptor Kit for Trekker Base ♦ ♦Required with Trekker Frame ^ ^Stingray Seat Size 1 only compatible with Trekker Frame	C905363	N/C


Description	Part No.	CAD Price
<input type="checkbox"/> Using other manufacturer's Head Support ♦Requires Head Support & Pelvic Positioning Belt	C995681	N/C

Description	Part No.	CAD Price
<input type="checkbox"/> 12½ x 2" Rear Solid Knobby Tire (TR)••	C904117	\$248
<input type="checkbox"/> 7½" x 2" Front & 11½" x 2½" Rear Pneumatic Knobby Tires (TR)•• ••Custom order; Approx. lead time: 6-8 weeks	C904118	\$373

Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> IV Pole (height adj. & collapsible) (TR) Anti-Tipper Required		One Size	C904566	\$250
	<input type="checkbox"/> O2 Tank Holder (TR) *°		One Size	C904565	\$218
	<input type="checkbox"/> LTV Ventilator Holder (TR) °		One Size	C904671	\$561
	<input type="checkbox"/> Ventilator/Suction/Hard Tray (TR)		One Size	C904595	\$561
	<input type="checkbox"/> Extended Ventilator/Suction/Hard Tray (TR)		One Size	C905326	\$655
	<input type="checkbox"/> Curb Tipper (TR)		One Size	C905310	\$122
	<input type="checkbox"/> Rear Anti-Tipper (TR)		One Size	C904649	\$143
	<input type="checkbox"/> Under Seat Storage Basket (TR)		One Size	C904605	\$150
	<input type="checkbox"/> Medical Necessity Storage Basket (TR)		One Size	C904606	\$368
	* Must be removed during transit				
	° Manufacturer Installation Required				

High-Low Frame (Trekker)	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> High-Low:x sz 1, Manual Lift^		1,2	8911603-02	\$1,970
	<input type="checkbox"/> Push Brace - Fixed Angle		One Size	8910545	\$200
	<input type="checkbox"/> Push Brace - Angle Adjustable		1,2	8910555	\$480
	^ Interface Adaptor Required				

High-Low Frame	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> High-Low:v Indoor Frame		1	8911871-20	\$2,207
	<input type="checkbox"/> High-Low:v Indoor Frame		1, 2	8911872-20	\$2,735
	<input type="checkbox"/> Push Brace - Fixed Angle		One Size	8910545	\$200
	<input type="checkbox"/> Push Brace, Angle Adjustable		One Size	8910555	\$480

Stingray Seat	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Seat and Back size 1^^		1	8711120	
	<input type="checkbox"/> Seat size 1		1	x8711000	\$511
	<input type="checkbox"/> Back size 1		1	x8712000	\$630
	<input type="checkbox"/> Seat and Back size 2		2	8712220	
	<input type="checkbox"/> Seat size 2		2	x8721000	\$511
	<input type="checkbox"/> Back size 2		2	x8722000	\$630
	<input checked="" type="checkbox"/> Fully Reclining Back (Required)		One Size	x8702000	\$505
	<i>Adjustable seat depth and height adjustable back</i>				
	^ Stingray Seat Size 1 only compatible with Trekker Frame				

Description	Size	Part No.	CAD Price
-------------	------	----------	-----------

Required, must select a color

X = Color

- Blue = 1 Red = 2
 Pink = 5 Grey = 4

- Seat and Back Cushion (Pr.) 1 875051x N/C
 Seat and Back Cushion (Pr.) 2 875052x N/C












***Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item**

Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------

- Footrest Rod, Elevating ♦ 8¼" Short x8711100 \$137
 Footrest Rod, Elevating ♦ 11" Long x8712200 \$137
♦ **Required with Angle Adjustable Footplate (x8710010)**

- Angle Adjustable Footplate ♦ One Size x8710010 \$119
♦ **Requires Footrest Rod Short or Long (x8711100 or x8712200)**



Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.) ♦		1	8740250	\$281
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) (Required) ♦		1	x8740250	\$281
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.) ♦		2	8740260	\$281
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) (Required) ♦		2	x8740260	\$281
	♦ Hip Support and Multi-Adjustable Hardware must be ordered together				
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	4" x 4"	1	8720301	\$243
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	5" x 5¾"	2	8720302	\$243
<hr/>					
	<input type="checkbox"/> Covers for Knobs on Side Supports (Pr.) ♦♦		One Size	8720360	\$38
	♦♦ Custom order; Approx. lead time: 6-8 weeks. Non-returnable item				
	<input type="checkbox"/> Medial Thigh Support (Ea.) ♦	3½" x 1½" x 2¾"	1	8720501	\$168
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		1	x8720501	\$206
	<input type="checkbox"/> Medial Thigh Support (Ea.) **	3½" x 3½" x 3½"	2	8720502	\$168
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		2	x8720502	\$206
	<input type="checkbox"/> Medial Thigh Support (Ea.) **	3½" x 3½" x 4¾"	3	8720503	\$168
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		3	x8720503	\$206
	♦ Medial Thigh Support and Adjustable, Removable Hardware must be ordered together				
	<input type="checkbox"/> Anti-Thrust Wedge, Black (Pr.)	4" x 6"	1, 2, 3	8910400-3	\$57
<hr/>					
	<input type="checkbox"/> Contoured Headrest, Black		1	8740810	\$249
	<input type="checkbox"/> Contoured Headrest, Black		2	8740820	\$249
	<input type="checkbox"/> Sculpted Head Support ①	9½" x 4½"	Small	8720831	\$256
	<input type="checkbox"/> Sculpted Head Support ①	11" x 5¾"	Large	8720832	\$256
	① Preferred option for Stingray seating				
	<input type="checkbox"/> Anatomic Head Support ♦	A:5½ x B:5½ x C:9½"	Small	5099001	\$318
	<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	N/C
	<input type="checkbox"/> Anatomic Head Support ♦	A:6 x B:5½ x C:10"	Medium	5099002	\$318
	<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	N/C
	<input type="checkbox"/> Anatomic Head Support ♦	A:6½ x B:5½ x C:12½"	Large	5099003	\$318
	<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	N/C
	♦ Turtle Bar System Short or Long Required; 98619-25 or 98619-26				
	Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable				
	<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$287
	<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$287
	♦ Required with Anatomic Head Support; 5099001, 5099002 or 5099003				
	<input type="checkbox"/> Headrest Fittings ♦		One Size	8720850	\$64
	♦ Required with Turtle Bar System or Universal Headrest Mount				

Accessories	Description	Measurements	Size	Part No.	CAD Price
-------------	-------------	--------------	------	----------	-----------



2-Part Headrest
(Mounts directly onto backrest) One Size 8720870 \$337

Universal Headrest Mount One Size 87590 \$76


◆ Required when using other manufactures Head Support

Allow most Whitmyer, Stealth, Otto Bock headrest, I2I to mount to our products (x:panda, Serval, etc.). Please contact Customer Service for confirmation


Other types of headrest available, please see Head Supports Order Form

Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------


■ XX = 60 (Black)



Cross Vest (Ea.) A: 4¼ x B:½ x C:11" 0 89625-0XX \$237
 Cross Vest (Ea.) A: 5 x B:½ x C:12" 1 89625-1XX \$237
 Cross Vest (Ea.) A: 5½ x B:½ x C:13" 2 89625-2XX \$237
 Cross Vest (Ea.) A: 6¼ x B:2½ x C:15¼" 3 89625-3XX \$237
 Cross Vest (Ea.) A: 7½ x B:2½ x C:18½" 4 89625-4XX \$237




Vest (Ea.) 1¼" x 2½" x 5½" x 3½" x 0 89111-0XX \$249
 Vest (Ea.) 1¼" x 2½" x 5¾" x 4¼" x ¾ 89111-075X \$249
 Vest (Ea.) 2¼" x 2½" x 7" x 5" x 13½ 1 89111-1XX \$249
 Vest (Ea.) 2½" x 3½" x 9¼" x 5½" x 2 89111-2XX \$249
 Vest (Ea.) 2¾" x 4¾" x 11 x 5¾" x 1 3 89111-3XX \$249
 Vest (Ea.) 2¾" x 6¾" x 13¾" x 7½" x 4 89111-4XX \$249



Neoprene Vest (Ea.) A:7 x B:5½ x C:4 x D:4½' 1 89113-1 \$249
 Neoprene Vest (Ea.) A:10½ x B:11 x C:6½ x D 2 89113-2 \$249
 Neoprene Vest (Ea.) A:11¼ x B:9 x C:7 x D:6½ 3 89113-3 \$249




4-Point Padded Hip Belt L:10¼-11½" 1 85427-1 \$168
 4-Point Padded Hip Belt L: 12½-13½" 2 85427-2 \$168



2-Point Pelvic Belt, Black L: 13-14½" 1 85424-1 \$121



Fix-Locks, Black, Bag of 6 Pcs ◆ One Size 862109-35 \$52
 ◆ Required for use of Vests



Foot Straps, Plastic Lock (Pr.) Max: 22", A: 13½-20½" One Size 85096 \$70
 Foot Straps, Velcro (Pr.) One Size 85102 \$64

Other types of vests and belts available, please see Vest and Belts Order Form

Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Folding Top (Canopy), Black		1	8740110	\$412
	<input type="checkbox"/> Folding Top (Canopy), Black		2	8740120	\$412
	<input type="checkbox"/> Foam Back Plate ●●	5" x 12"	1	8776325	\$29
	<input type="checkbox"/> Foam Back Plate ●●	5½" x 13¾"	2	8776326	\$29
	Future Growth - Fills out the gap between seat and the back				
	●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item				
	<input type="checkbox"/> Tray with Edge, Grey ♦	17.3" x 13.3"	One Size	8776350	\$125
	<input type="checkbox"/> Adjustable Removable Tray Hardware ♦		One Size	x8776350	\$206
	♦ Tray and Adjustable Removable Hardware must be ordered together				
	<input type="checkbox"/> Rain Cover ♦		1	8720141	\$70
	<input type="checkbox"/> Rain Cover ♦		2	8720142	\$70
	♦ Requires Folding Top (Canopy)				
	<input type="checkbox"/> Handle Grip	14"	1	8720401	\$249
	<input type="checkbox"/> Handle Grip	16"	2	8720402	\$249
	<input type="checkbox"/> Backpack, Gray ●●		1	8720711	\$70
	<input type="checkbox"/> Backpack, Gray ●●		2	8720721	\$70
	●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item				

Spare Parts	Description	Measurements	Size	Part No.	CAD Price
	Required, must select a color				
	X = Color				
	<input type="checkbox"/> Blue = 1	<input type="checkbox"/> Red = 2			
	<input type="checkbox"/> Pink = 5	<input type="checkbox"/> Grey = 4			
	<input type="checkbox"/> Seat and Back Cushion (Pr.) ●●		1	875051x	\$433
	<input type="checkbox"/> Seat and Back Cushion (Pr.) ●●		2	875052x	\$433
	●●Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item				
	<input type="checkbox"/> Adductor Covers, Black (Pr.)		1	8750410	\$70
	<input type="checkbox"/> Adductor Covers, Black (Pr.)		2	8750420	\$70
	<input type="checkbox"/> Hardware Kit for Hip Supports, Fixed		One Size	872023-KIT	\$44

Spare Parts	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Release Strap for Back Angling		1	87102250	\$57
	<input type="checkbox"/> Release Strap for Back Angling		2	87102255	\$57
	<input type="checkbox"/> Footrest Rod, Elevating	8¼"	Short	8775317	\$137
	<input type="checkbox"/> Footrest Rod, Elevating	11"	Long	8775318	\$137
	<input type="checkbox"/> Angle Adjustable Footplate		One Size	8775010	\$324

When being transported in a vehicle, it is R82's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair; therefore R82 cannot recommend a transit system.

The R82 Stingray mobility base was tested with a footrest and R82 transit anchors. The performance and safety of accessories and/or other seating options not tested by R82 is unknown, and no warranty or claim about their performance is or can be made.

For shipping rates, please call 1-844-876-6245
12801 E Independence Blvd, Matthews, NC 28105 • 1-844-876-6245 • Web: www.r82.com
Prices in CAD & is subject to change without notice
For Quotes and Orders please send completed Order Form to: us.info@r82.com