


FOOT & LEG POSITIONING

Footplate Securement Straps

996619

Item	Description	Qty
1	Footplate Securement Strap	1



If you have any questions about attaching the Footplate Securement Strap, please call our Sales Department at 888 CONVAID (266-8243) or +1 310 618-0111.

 **Convaid**
Now You're Going Places

INSTRUCTIONS

Footplate Securement Straps

1. The securement strap holds individual footplates together to prevent them from flipping up and down.
2. To secure the footplates, buckle the left and right strap together.