

FOOT & LEG POSITIONING

Foot Straps Positioners - Pair

996616


Item	Description	Qty
1	Foot Positioners - Small	1

996617

Item	Description	Qty
1	Foot Positioners - Medium	1

996618

Item	Description	Qty
1	Foot Positioners - Large	1




If you have any questions about attaching the Foot Positioners, please call our Sales Department at 888 CONVAID (266-8243) or +1 310 618-0111.

INSTRUCTIONS

Foot Straps Positioners - Pair



1. Foot Positioners may be criss-crossed over the foot to secure the whole foot or can be converted into simple ankle straps.
2. To attach foot positioner, thread strap through footplate then bolt strap to underside of footplate



If you have any questions about attaching the One Piece Handle, please call our Sales Department at 888 CONVAID (266-8243) or +1 310 618-0111.

 **Convaid**
Now You're Going Places