


HEAD SUPPORT

Headrest Extension - Curved Cordura

995927

Item	Description	Qty
1	Headrest Extension - Curved Cordura	1



If you have any questions about attaching the Headrest Extension, please call our Sales Department at 888 CONVAID (266-8243) or +1 310 618-0111.

INSTRUCTIONS

Headrest Extension

1. Insert Headrest Extension on Tubes into the Headrest Holding Brackets and push down.

If you have any questions about attaching the Headrest Extension, please call our Sales Department at 888 CONVAID (266-8243) or +1 310 618-0111.

 **Convaid**
Now You're Going Places