


# HEAD SUPPORT

## Occi Headwings

# 995932

Item	Description	Qty
1	Occi Headwings	1



If you have any questions about attaching the Occi Headwings, please call our Sales Department at 888 CONVAID (266-8243) or +1 310 618-0111.

## **INSTRUCTIONS**

---

### **Occi Headwings**

1. Occi Headwings supports and positions the head at the occipital area. This headrest provides added comfort, allows side to side head movement without obstruction, and does not block the individual's ears and line of vision.
2. The Occi headwing can be attached at any height. Secure headwing by tying laces through the grommet holes at the back of the seat fabric or headrest extension.