


PELVIC POSITIONING

Adjustable Crotch Strap for 3 Point Belt

994716

Item	Description	Qty
1	Crotch Pad Assembly	1



If you have any questions about attaching the Crotch Strap, please call our Sales Department at 888 CONVAID (266-8243) or +1 310 618-0111.

 **Convaaid**
Now You're Going Places

INSTRUCTIONS

3 Point Belt with Depth Adjustable Crotch Strap

1. The crotch strap can be adjusted by threading the strap through the desired slot.

If you have any questions about attaching the Crotch Strap,
please call our Sales Department at 888 CONVAID (266-8243) or +1 310 618-0111.

 **Convaaid**
Now You're Going Places