

TRUNK POSITIONING


Full Torso Support Vest

995925

Item	Description	Qty
1	Full Torso Support Vest for Cordura Upholstery	1

995926

Item	Description	Qty
1	Full Torso Support Vest for Textilene Upholstery	1



If you have any questions about attaching the Full Torso Support Vest, please call our Sales Department at 888 CONVAID (266-8243) or +1 310 618-0111.

INSTRUCTIONS

Full Torso Support Vest

- An Adjustable support vest keeps the user in place comfortably and securely. It helps to maintain midline seating position and attached to the chair with straps that wrap around the seat back and connect with Velcro[®], and shoulder straps that screw into the seat back. Adjust the side straps for proper fit.