

# R82 Cricket

## CAD Retail Price List/ Order Form



Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### US Configuration


Effective: 05.08.19 - Rev 4

- Cricket
- Cushion
- Arm rests
- Head support
- Tray
- Grip bar
- Folding top
- Transport fittings



- A** (Back Height) \_\_\_\_\_
- B** (Arm Height) \_\_\_\_\_
- C** (Lower Leg Length) \_\_\_\_\_
- D** (Seat Depth) \_\_\_\_\_
- E** (Chest Width) \_\_\_\_\_
- F** (Seat Width) \_\_\_\_\_
- G** (Height from Seat to Top of Head) \_\_\_\_\_
- H** (Shoulder Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_

### Measurements



**Size 1**  
Inches



**Size 2**  
Inches

Effective seat width (A)	13½"	15½"
Effective seat depth (B)	8½-12¼"	12¼-17"
Backrest height (C1)	19¾"	22½"
Seat surface height at front edge (E)	19¾-21"	21¼-22¼"
Footrest to seat distance (D)	7-14¼"	8¾-14½"
Leg to seat surface angle	360°	360°
Backrest width, outside (F1)	13½"	15½"
Backrest width, inside (F)	10¾"	13"
Overall width (G)	20¾"	22¾"
Length (H)	29½"	29½"
Overall height (I)	39¼"	41¼"
Folded width	14"	14"
Folded height	19"	19"
Folding length	30"	30"
Backrest angle	110°-125°	110°-125°
Minimum turning radius	27½"	30"
Total mass	27.3 lbs.	29½ lbs.
Max load	77 lbs.	110 lbs.
Max load at transportation	59½ lbs.	94½ lbs.
Required width of right angled corridor	32"	32"
Static stability (sideways, downhill, uphill)	10°	10°

### Notes:

### Cricket Standard Features




- Back Recline
- Adjustable Seat Depth
- Height and Angle Adjustable Foot Support
- Armrest, Fixed (Pr.)
- Angle Adjustable Pushbrace
- Foot Operated Wheel Lock
- Castors - 7" Front 11" Rear Solid

Description	Part No.	CAD Price	HCPCS
<input type="checkbox"/> Cricket sz 1	87043100	\$2,382	
<input type="checkbox"/> Cricket sz 2	87043200	\$2,383	

Description	Size	Part No.	CAD Price	HCPCS
-------------	------	----------	-----------	-------

**Required, must select a color**

- Blue = 1       Red = 2  
 Pink = 5       Grey = 4



<input type="checkbox"/> Cushion	1	8704311x	\$0	
<input type="checkbox"/> Cushion	2	8704321x	\$0	

◆ Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item

<input type="checkbox"/> Seat Side Parts ◆	One Size	87043130	\$65	
--	----------	----------	------	--

◆ Required with Seat Width Reduction Kit

Description	Measurements	Size	Part No.	CAD Price	HCPCS
-------------	--------------	------	----------	-----------	-------

<input type="checkbox"/> Transit Option		One Size	87043188	\$249	
---	--	----------	----------	-------	--



Description	Measurements	Size	Part No.	CAD Price	HCPCS
-------------	--------------	------	----------	-----------	-------

<input type="checkbox"/> Armrest, Fixed (Pr.) ◆		1	87043120-1	\$0	
<input type="checkbox"/> Armrest, Fixed (Pr.) ◆		2	87043120-2	\$0	


◆ Required with Grip Bar, In Armrests



Description	Measurements	Size	Part No.	CAD Price	HCPCS
-------------	--------------	------	----------	-----------	-------

<input type="checkbox"/> Seat Width Reduction Kit, Wedges (Pr.) ◆			87043130-1	\$90	K0108
<input type="checkbox"/> Seat Width Reduction Kit, Wedges (Pr.) ◆			87043130-2	\$90	K0108


◆ Requires Seat Side Parts (87043130)




<input type="checkbox"/> Wedges for Shaping Back (Pr.)	4½" x 6" x 11"	Type 1	87043140-1	\$90	K0108
<input type="checkbox"/> Wedges for Shaping Back (Pr.)	5¾" x 6" x 11"	Type 2	87043140-2	\$90	K0108



<input type="checkbox"/> Chest Belt, Upholstery (Ea.)	A:11¼ x B:7½ x C:9¼-12"	1	85428-1	\$337	E0960
<input type="checkbox"/> Chest Belt, Upholstery (Ea.)	A:14 x B:10 x C:10-12¾"	2	85428-2	\$337	E0960




<input type="checkbox"/> Neoprene Vest (Ea.)	A:7 x B:5½ x C:4 x D:4½"	1	89113-1	\$249	E0960
<input type="checkbox"/> Neoprene Vest (Ea.)	A:10½ x B:11 x C:6½ x D:5⁹"	2	89113-2	\$249	E0960




Accessories	Description	Measurements	Size	Part No.	CAD Price	HCPCS
-------------	-------------	--------------	------	----------	-----------	-------

<input type="checkbox"/> 2-Point Pelvic Belt, Black	L: 13-14½"	1	85424-1	\$121	E0978
---	------------	---	---------	-------	-------




<input type="checkbox"/> 4-Point Padded Hip Belt	L:10¼-11½"	1	85427-1	\$168	E0978
<input type="checkbox"/> 4-Point Padded Hip Belt	L: 12½-13½"	2	85427-2	\$168	E0978



<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs ◆		One Size	862109-35	\$52	
---	--	----------	-----------	------	--

◆ Required for use of Vests






Description	Measurements	Size	Part No.	CAD Price	HCPCS
<input type="checkbox"/> Pelvic Harness	A:6½ x B:7½"	0	89108-0	\$306	E0978
<input type="checkbox"/> Pelvic Harness	A:8½ x B:8¾"	1	89108-1	\$306	E0978
<input type="checkbox"/> Pelvic Harness	A:10½ x B:9¾"	2	89108-2	\$306	E0978



<input type="checkbox"/> Heel Straps (Pr.)		One Size	87043185	\$67	E0951
--	--	----------	----------	------	-------



<input type="checkbox"/> Upholstered Ankle Straps	A: 5½-7 x B:8½-10"	1	85426-1	\$231	K0108
<input type="checkbox"/> Upholstered Ankle Straps	A: 6¼-7¾ x B:10-11¾"	2	85426-2	\$231	K0108

Other types of vests and belts available, please see Vest and Belts Order Form

Description	Measurements	Size	Part No.	CAD Price	HCPCS
-------------	--------------	------	----------	-----------	-------



<input type="checkbox"/> Headrest Cushion, Black (Ea.)		One Size	87043150	\$249	E0955
--	--	----------	----------	-------	-------




<input type="checkbox"/> Sculpted Headrest Cushion, Black (Ea.)		One Size	87043151	\$256	E0955
---	--	----------	----------	-------	-------




<input type="checkbox"/> Backrest Height Extension (Ea.)		1	87043160-1	\$69	K0108
<input type="checkbox"/> Backrest Height Extension (Ea.)		2	87043160-2	\$69	K0108

Description	Measurements	Size	Part No.	CAD Price	HCPCS
-------------	--------------	------	----------	-----------	-------

<input type="checkbox"/> Calf Support Strap (Ea.)		1	87043190-1	\$76	E0995
<input type="checkbox"/> Calf Support Strap (Ea.)		2	87043190-2	\$76	E0995



<input type="checkbox"/> Additional Footplate Locking Mechanism (Ea.)		One Size	87043186	\$38	
---	--	----------	----------	------	--




<input type="checkbox"/> Grip Bar, In Armrests ♦		1	87043181-1	\$215	
<input type="checkbox"/> Grip Bar, In Armrests ♦		2	87043181-2	\$215	

♦ Requires Armrest



<input type="checkbox"/> Grip Bar, In Frame		1	87043182-1	\$215	
<input type="checkbox"/> Grip Bar, In Frame		2	87043182-2	\$215	



<input type="checkbox"/> Caster Locks		One Size	87043187	\$117	K0073
---------------------------------------	--	----------	----------	-------	-------

Accessories	Description	Measurements	Size	Part No.	CAD Price	HCPCS
-------------	-------------	--------------	------	----------	-----------	-------



Required, must select a color

x = Color

- Blue = 1       Red = 2  
 Pink = 5       Grey = 4

<input type="checkbox"/> Folding Top (Canopy) (Ea.) ♦		1	8704317x	\$284	
<input type="checkbox"/> Folding Top (Canopy) (Ea.) ♦		2	8704327x	\$284	

♦ Required w/ Rain Cover



<input type="checkbox"/> Underseat Storage Basket ♦		1	87043184-1	\$87	
<input type="checkbox"/> Underseat Storage Basket ♦		2	87043184-2	\$87	

♦ 22 lbs. weight limit



<input type="checkbox"/> Tray ♦		1	87043180-1	\$422	E0950
<input type="checkbox"/> Tray ♦		2	87043180-2	\$422	E0950

♦ Must select same size as armrest





Storage Strap

One Size 85100

\$38



Rain Cape with Hood, Clear

A:11 x B:29½"

4

87465-410

\$175

Rain Cape with Hood, Clear

A:11½ x B:34½"

5

87465-510

\$175

Rain Cover, Clear (Ea.) ♦

2

87467-210

\$144

Rain Cover, Clear (Ea.) ♦

3

87467-310

\$144

♦Requires Folding Top (Canopy)

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

The Cricket complies with ISO7176-19, for safe transportation in motor vehicles. The R82 Cricket was tested with the Cricket back rest height extension and Y-hip belt and R82 transit anchors. The performance and safety of accessories and/or other seating options not tested by R82 is unknown, and no warranty or claim about their performance is or can be made.

When being transported in a vehicle, it is R82's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair; therefore R82 cannot recommend a transit system.

**For shipping rates, please call 1-844-876-6245**

12801 E Independence Blvd, Matthews, NC 28105 • 1-844-876-6245 • Web: www.r82.com

Prices in CAD & is subject to change without notice

**For Quotes and Orders please send completed Order Form to: us.info@r82.com**