



# Our Family of Transportation Products



Main Features	Quokka™	Wallaroo®	Carrot 3	Pilot Booster	Carrot 3 Booster Seat
Seat Width	11"	13 ¾"	11 ½"	12"	16 ½"
Seat Depth	11"	10 ½" – 15"	9 ¾" – 22 ¼"	14" – 19"	15" – 16 ½"
Shoulder Harness Range	7 ½" – 13 ½"	11" – 17"	11 ¾" – 22 ¾"	22" – 28 ½"	28 ¾" – 32 ¾"
Top of Head Max	21"	26 ¾"	31 ¼"	28 ½"	33"
User Height Range	Up to 41"	Up to 56"	37" – 60"	38" – 63"	55" – 69"
User Weight Range	5 – 39.6 lb	22 – 106 lb	30 – 108 lb	30 – 120 lb	79 – 165 lb
Product Weight	15.8 lb	27 lb	19.1 lb	15.7 lb	21 lb
Price	\$700	\$995	\$2,143	\$605	\$2,750

\*Pricing & specifications are subject to change



# Our Family of Transportation Products



Main Features	Quokka™	Wallaroo®	Carrot 3*	Pilot Booster*	Carrot 3 Booster Seat*
Rear Facing	✓	X	X	X	X
Forward Facing	✓	✓	✓	✓	✓
Seat Depth Extender	X	✓	✓	✓	X
Back Height Growth	✓	✓	✓	✓	✓
Spica Compatible	✓	✓	X	X	X
5 Point Harness	✓	✓	✓	X	X
Tray	X	✓	✓	X	✓

\*Requires the use of the vehicle seat belt