

Convoid Metro

Canadian Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

Metro Standard Features:

- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (ME12/14/16)
- Quick release Rolko rear wheels (ME12/14/16)
- 8" Front, 11" Rear Solid Tire (ME18)
- Wheel Locks
- Adjustable Swing-Away Footplates
- 30° Fixed Tilt
- Silver Powder Coated Frame
- 2-Pt Positioning Belt
- Made in the USA

Standard Models:

	<u>Price</u>
<input type="checkbox"/> Metro ME12 - 901402	\$1,383
<input type="checkbox"/> Metro ME14 - 901400	\$1,383
<input type="checkbox"/> Metro ME16 - 901335	\$1,383
<input type="checkbox"/> Metro ME18 - 903127	\$1,383

All prices in CAD

<u>Upholstery Colors</u>	<u>Price</u>
<input type="checkbox"/> Cherry Red	\$0
<input type="checkbox"/> Royal Blue	\$0

ME12	ME14	ME16	ME18
903855	903855	903855	903855
903850	903850	903850	903850

<u>Transit Required Items</u>	<u>Price</u>
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	\$263
<input type="checkbox"/> Headrest Extension	\$172
<input type="checkbox"/> 3-Pt Positioning Belt	\$65
Foot Positioners - Select Size Below:	
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$206
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$206
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$206

ME12	ME14	ME16	ME18
901403	901401	901336	903128
903263	903264	903265	903266
904017	904017	904017	904017
/			
903415	903415	903415	903415
903416	903416	903416	903416
903417	903417	903417	903417

<u>Transit Options</u>	<u>Price</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$118
<input type="checkbox"/> Sure-Lok Transit Lap Belt	\$118

ME12	ME14	ME16	ME18
903094	903094	903094	903094
903232	903232	903232	903232

<u>Tire Options</u>	<u>Price</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	\$260
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$392
<input type="checkbox"/> 8"x2" Front Ribbed, 12.5" Rear Solid Knobby Tire	\$392
<input type="checkbox"/> 8"x2" Front Ribbed, 12.5" Rear Pneumatic Knobby Tire	\$392

ME12	ME14	ME16	ME18
904117	904117	904117	
904118	904118	904118	
			903084
			903083

<u>Accessories</u>	<u>Price</u>
<input type="checkbox"/> 3-Pt Positioning Belt	\$65
<input type="checkbox"/> Wheelchair Travel Bag	\$270
<input type="checkbox"/> Under Seat Storage Basket - Fixed Tilt	\$89

ME12	ME14	ME16	ME18
904017	904017	904017	904017
902022	902022	902022	902022
904046	904047	904047	904048

The Metro Transit option includes four manufacturer installed wheelchair transportation anchors. The Metro transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are ME12: 66 lbs./30kg; ME14: 100 lbs./45.5kg/ME16 and ME18: 170 lbs./77 kg.

The Convoid Metro was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is being, or can be, made.

Standard Shipping Rates: Ground \$69 • For expedited shipping rates, please call 1-844-876-6245

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166

Prices in CAD & subject to change, ex works Torrance, California.

www.convoid.com

Effective: 6/22/2020

OF2002-ENG-REV.01