

# Convaid Rodeo Tilt w/out Seating - HCPCS E1234

## Retail Price List / Order Form

Account #: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

Date: \_\_\_\_\_

Contact: \_\_\_\_\_

Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

### Rodeo Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- Adjustable One-Piece Push Handle
- Elevating Legrest / One-Piece, Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires
- Quick Release Rolko Rear Wheels
- Single Foot Activated Wheel Lock Bar
- Seat Depth Growth of 5"
- Made in the USA

#### Standard Models:

	Price
<input type="checkbox"/> Rodeo w/out Seating RD10WOS - 904963	<b>\$2,318</b>
<input type="checkbox"/> Rodeo w/out Seating RD12WOS - 904964	<b>\$2,318</b>
<input type="checkbox"/> Rodeo w/out Seating RD14WOS - 904965	<b>\$2,318</b>
<input type="checkbox"/> Rodeo w/out Seating RD16WOS - 904966	<b>\$2,318</b>

- Seating Module** Price HCPCS
- Depth Adjustable Solid Seat Pan (Manufacturer Installation Required) **\$150** **E2293**

RD10	RD12	RD14	RD16
904971	904972	904973	904974

- Transit Required Items** Price HCPCS
- Transit Option (Manufacturer Installation Required) **\$200**
- Using other manufacturer's Head Support **\$0**
- 2-Point Positioning Belt **\$45** **E0978**

RD10	RD12	RD14	RD16
904341	903167	903168	903169
995681	995681	995681	995681
903931	903931	904559	904559

- Transit Options** Price HCPCS
- Q'Straint Transit Lap Belt **\$90** **E0978**
- Sure-Lok Transit Lap Belt **\$90** **E0978**
- Upper Extremity Support Surface (Foam Tray) Small ^ **\$151** **E0950**
- Upper Extremity Support Surface (Foam Tray) Medium ^ **\$151** **E0950**
- ^ Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

RD10	RD12	RD14	RD16
903094	903094	903094	903094
903232	903232	903232	903232
903669	903669		
		903670	903670

- Foot & Leg Positioning** Price HCPCS
- Individual Angle Adjustable Footplates **\$179** **K0040**
- Full Padded Footbox (One-Piece Footplate Required) **\$325** **K0108**
- Legrest Extension - 2" \* **\$55**
- Foot Positioners - Small 7" - 10" (pair) **\$156** **K0108**
- Foot Positioners - Medium 8" - 11" (pair) **\$156** **K0108**
- Foot Positioners - Large 9" - 12" (pair) **\$156** **K0108**

RD10	RD12	RD14	RD16
	904957	904958	904959
904953	904954	904955	904956
904095	904095	904096	904097
903415	903415	903415	903415
903416	903416	903416	903416
903417	903417	903417	903417

\* 2" Legrest Extension interferes with caster swivel when seat is in upright position.  
Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

- Frame Option** Price HCPCS
- Recline Option (Manufacturer Installation Required) **\$461** **E1225**
- Height Adjustable Flip-Up Armrests (pair) **\$228** **E0973**
- Solid Side Panels **\$211** **K0108**
- Mesh Side Panels **\$148** **K0108**

RD10	RD12	RD14	RD16
904342	904343	904344	904345
904346	904346	904346	904346
904388	902672	902652	903145
904389	902680	902681	903146

- Tire Options** Price HCPCS
- 12.5" Rear Solid Knobby Tire **\$199**
- 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire **\$299**

RD10	RD12	RD14	RD16
904117	904117	904117	904117
904118	904118	904118	904118

- Accessories** Price HCPCS
- Headrest Cover (Canopy) **\$228**
- Extended Headrest Cover - No Windows (Canopy) **\$250**
- Extended Headrest Cover - with Windows (Canopy) **\$295**
- Upper Extremity Support Surface (Clear Tray) ^ **\$111** **E0950**
- Adj. Removable Hardware (Required for Clear Tray) **\$227** **E1028**
- Medical Necessity Storage Basket **\$295**
- Under Seat Storage Basket - Tilt-in-Space **\$75**
- IV Pole (height adjustable and collapsible)<sup>2</sup> **\$150** **K0105**
- Oxygen Tank Bag **\$179** **E2208**
- Utility Bag **\$38**
- Attendant Hand Brakes (pair) <sup>1</sup> **\$236**
- Curb Tipper **\$98**
- Caster Lock for 2" Wide Front Tires <sup>1</sup> **\$94** **K0073**
- Rear Anti-Tip Tubes (Pair)\* **\$115** **E0971**
- Rain Cover Size 1 (Canopy Req) **\$150**
- Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req) **\$150**
- Rain Cover Size 3 (Extended Canopy Req) **\$150**
- Mosquito Net Size 1 (Canopy req) **\$150**
- Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req) **\$150**
- Mosquito Net Size 3 (Extended Canopy req) **\$150**

RD10	RD12	RD14	RD16
904366	902679	902659	903152
904728	904729	904730	904731
904744	904745	904746	904747
904367	903405	903401	903402
903369	903369	903369	903369
904476	904476	904477	904477
904529	904529	904530	904530
904566	904566	904566	904566
902870	902870	902870	902870
904531	904531	904531	904531
		904370	904371
905310	905310	905310	905310
904650	904650	904650	904650
904649	904649	904649	904649
904939	904939	904939	
		904940	904940
			905283
904655	904655	904655	
		904656	904656
			905319

^ Flip-Up Armrests and Hardware required.  
\*Highly recommended for users weighing 50lbs or less.  
<sup>1</sup> Manufacturer Installation Required.  
<sup>2</sup> Anti-tip Tubes Required

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convoid Rodeo Mobility Base was tested with the Rodeo Seating Module, equipped with a head support and wheelchair-anchored 2-point pelvic lap belt. Do not use the foot positioners on the Convoid Rodeo during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is or can being made.

It is Convoid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convoid cannot recommend a transit system.

The WC19 standard is specifically designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.rercwts.pitt.edu/WC19.html>.

*HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.*

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245**

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166

Prices in USD & subject to change, ex works Torrance, California.

Effective: 8/24/2020

[www.convoid.com](http://www.convoid.com)

OF0109-ENG-REV.02